

SmartPeg Reference Guide

Valid from September 2017 | 25016-65

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
A.B. DENTAL DEVICES	Screw Type, Conical Groovy				3.75		100440	Type 32	
	Screw Type, Conical Groovy				4.5		100431	Type 27	
ADDENTA	Bifasico				4.0		100374	Type 14	
ADIN DENTAL IMPLANT SYSTEM	One				3.3		100390	Type 17	
	One				4.2		100390	Type 17	
	Swell				3.3		100480	Type 49	
	Swell				3.75		100480	Type 49	
	Swell				4.2		100480	Type 49	
	Swell				5.0		100480	Type 49	
	Swell				6.0		100480	Type 49	
	Touareg CloseFit				NP		100545	Type 62	
	Touareg CloseFit				RP		100543	Type 60	
	Touareg CloseFit				WP		100480	Type 49	
	Touareg-S				4.2		100480	Type 49	
Touareg-S				5.0		100480	Type 49		
AGS MEDIKAL	Implance Dental Implant System, Bone Level						100378	Type 06	
	Implance Dental Implant System, Mini Type						100455	Type 38	
	Implance Dental Implant System, Tissue Level						100378	Type 06	
ALPHA BIO	ATI						100425	Type 26	
	ATID				3.75		100425	Type 26	
	Neo CS				3.75		100583	Type 78	New
	NICE 3.2 Conical connection				3.2		100554	Type 63	
	SPI						100425	Type 26	
ALPHA-DENT	Active 3D				3.25		100480	Type 49	
	Active 3D				3.75		100480	Type 49	
	Active 3D				4.2		100480	Type 49	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
ALPHA-DENT	Active 3D				5.0	100480	Type 49	
	Active Conus				3.25	100480	Type 49	
	Active Conus				3.75	100480	Type 49	
	Active Conus				4.2	100480	Type 49	
	Active Conus				5.0	100480	Type 49	
	Active Plus				4.2	100480	Type 49	
	Active Plus				5.0	100480	Type 49	
	Classic 3D				3.3	100480	Type 49	
	Classic 3D				3.75	100480	Type 49	
	Classic 3D				4.2	100480	Type 49	
	Classic 3D				5.0	100480	Type 49	
	Classic Conus				3.3	100480	Type 49	
	Classic Conus				3.75	100480	Type 49	
	Classic Conus				4.2	100480	Type 49	
	Classic Conus				5.0	100480	Type 49	
ALPHA-GATE	Bioactive				3.75	100446	Type 35	
ALPHATECH					3.4	100477	Type 46	
					3.8	100477	Type 46	
					4.3	100477	Type 46	
					5.0	100477	Type 46	
					6.0	100477	Type 46	
					7.0	100477	Type 46	
ANKYLOS	C/X		A			100388	Type 16	
	C/X		B			100388	Type 16	
	C/X		C			100388	Type 16	
ANTHOGYR	Anthofit External Hexagon HE					100353	Type 01	
	Anthofit Internal Octagon OI					100398	Type 18	
	Axiom PX				3.4	100455	Type 38	
	Axiom PX				4.0	100455	Type 38	
	Axiom PX				4.6	100455	Type 38	
	Axiom PX				5.2	100455	Type 38	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
ANTHOGYR	Axiom REG				3.4	100455	Type 38	
	Axiom REG				4.0	100455	Type 38	
	Axiom REG				4.6	100455	Type 38	
	Axiom REG				5.2	100455	Type 38	
	Ossfit					100350	Type 04	
ARDS IMPLANTS	Classic				All	100431	Type 27	
	Premium				All	100431	Type 27	
	Smart				All	100431	Type 27	
AVINENT IMPLANT SYSTEM	Coral External				3.5	100436	Type 30	
	Coral External				4.1	100353	Type 01	
	Coral External				5.1	100357	Type 03	
	Internal Hexagon Coral				3.5	100440	Type 32	
	Internal Hexagon Coral				4.1	100440	Type 32	
	Internal Hexagon Ocean				3.5	100440	Type 32	
	Internal Hexagon Ocean				4.1	100440	Type 32	
	Ocean External				3.5	100436	Type 30	
	Ocean External				4.1	100353	Type 01	
	AXIS	Ceramic				4.0	100455	Type 38
B&B DENTAL	Duravit SLIM Implant				3.0	100404	Type 22	New
	Duravit SLIM Implant				3.4	100404	Type 22	New
	Duravit 3P Implant				3.5	100425	Type 26	New
	Duravit 3P Implant				4.0	100425	Type 26	New
	Duravit 3P Implant				4.5	100425	Type 26	New
	Duravit 3P Implant				5.0	100425	Type 26	New
	DURAVIT EV Implant				4.5	100425	Type 26	New
	DURAVIT EV Implant				4.0	100425	Type 26	New
	DURAVIT EV Implant				4.5	100425	Type 26	New
	DURAVIT EV Implant				5.0	100425	Type 26	New
	DURAVIT WIDE Implant				5.0	100425	Type 26	New
	DURAVIT WIDE Implant				6.0	100425	Type 26	New
	DURAVIT WIDE Implant				6.5	100425	Type 26	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
B&B DENTAL	DURAVIT WIDE Implant			7.0		100425	Type 26	New
BEGO IMPLANT SYSTEMS	Semados Mini			2.7		100355	Type 02	
	Semados Mini			2.9		100355	Type 02	
	Semados Mini			3.1		100355	Type 02	
	Semados RI-line			3.75		100425	Type 26	
	Semados RI-line			4.1		100425	Type 26	
	Semados RI-line			4.5		100425	Type 26	
	Semados RI-line			5.5		100425	Type 26	
	Semados RS/RSX			3.0		100481	Type 50	
	Semados RS/RSX			3.75		100425	Type 26	
	Semados RS/RSX			4.1		100425	Type 26	
	Semados RS/RSX			4.5		100425	Type 26	
	Semados RS/RSX			5.5		100425	Type 26	
	Semados S-line			3.25		100425	Type 26	
	Semados S-line			3.75		100425	Type 26	
	Semados S-line			4.1		100425	Type 26	
	Semados S-line			4.5		100425	Type 26	
	Semados S-line			5.5		100425	Type 26	
BIO3	Implant Advanced ATI			3.3		100480	Type 49	
	Implant Advanced ATI			3.8		100480	Type 49	
	Implant Advanced ATI			4.2		100480	Type 49	
	Implant Advanced ATI			5.0		100480	Type 49	
	Implant Progressive PTI			3.3		100480	Type 49	
	Implant Progressive PTI			3.8		100480	Type 49	
	Implant Progressive PTI			4.2		100480	Type 49	
	Implant Progressive PTI			5.0		100480	Type 49	
BIODENTA	Tissue level		NP			100402	Type 21	
	Tissue level		RP			100478	Type 47	
	Tissue level		WP			100478	Type 47	
BIOHORIZONS	External			3.5		100436	Type 30	
	External			4.0		100353	Type 01	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
BIOHORIZONS	External				5.0		100353	Type 01	
	External				6.0		100353	Type 01	
	Internal		Yellow				100440	Type 32	
	Internal		Green				100431	Type 27	
	Internal		Blue				100442	Type 33	
	Internal Laser Lok				3.0		100425	Type 26	
	Single Stage		Yellow				100440	Type 32	
	Single Stage		Green				100431	Type 27	
	Single Stage		Blue				100442	Type 33	
	Tapered Internal		Yellow				100440	Type 32	
	Tapered Internal		Green				100431	Type 27	
	Tapered Internal		Blue				100442	Type 33	
	BIOINFINITY	Junior Platform				3.2		100455	Type 38
		Standard Platform				4.8		100425	Type 26
BIOMET 3i	OSSEOTITE Certain		XP		4/5		100386	Type 15	
	OSSEOTITE Certain		XP		5/6		100386	Type 15	
	OSSEOTITE Tapered Certain		Purple				100472	Type 45	
	OSSEOTITE Tapered Certain		Blue				100386	Type 15	
	NanoTite Certain MicroMiniplant				3.25		100472	Type 45	
	NanoTite External Hex		Purple				100353	Type 01	
	NanoTite External Hex		Blue				100353	Type 01	
	NanoTite External Hex		Yellow				100353	Type 01	
	NanoTite External Hex		Green				100353	Type 01	
	NanoTite External Hex				XP 3/4		100353	Type 01	
	NanoTite External Hex				XP 4/5		100353	Type 01	
	NanoTite External Hex				XP 5/6		100353	Type 01	
	NanoTite External Hex MicroMiniplant		Purple				100353	Type 01	
	NanoTite Tapered Certain		Purple				100472	Type 45	
	NanoTite Tapered Certain		Blue				100386	Type 15	
	NanoTite Tapered Certain		Yellow				100386	Type 15	
	NanoTite Tapered Certain		Green				100386	Type 15	
	NanoTite Tapered Certain Prevail		Purple				100472	Type 45	
	NanoTite Tapered Certain Prevail		Blue				100386	Type 15	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
BIOMET 3i	NanoTite Tapered Certain Prevail		Yellow			100386	Type 15	
	OSSEOTITE Certain Microminiplant			3.25		100472	Type 45	
	OSSEOTITE External Hex		Purple			100353	Type 01	
	OSSEOTITE External Hex		Blue			100353	Type 01	
	OSSEOTITE External Hex		Yellow			100353	Type 01	
	OSSEOTITE External Hex		Green			100353	Type 01	
	OSSEOTITE External Hex			XP 3/4		100353	Type 01	
	OSSEOTITE External Hex			XP 4/5		100353	Type 01	
	OSSEOTITE External Hex			XP 5/6		100353	Type 01	
	OSSEOTITE External Hex MicroMiniplant		Purple			100353	Type 01	
	OSSEOTITE Tapered Certain		Yellow			100386	Type 15	
	OSSEOTITE Tapered Certain		Green			100386	Type 15	
	OSSEOTITE Tapered Certain Prevail		Purple			100472	Type 45	
	OSSEOTITE Tapered Certain Prevail		Blue			100386	Type 15	
	OSSEOTITE Tapered Certain Prevail		Yellow			100386	Type 15	
	T3 Non-Platform Switched		Purple			100472	Type 45	
	T3 Non-Platform Switched		Blue			100386	Type 15	
	T3 Non-Platform Switched		Yellow			100386	Type 15	
	T3 Non-Platform Switched		Green			100386	Type 15	
	T3 Platform Switched		Purple			100472	Type 45	
	T3 Platform Switched		Blue			100386	Type 15	
	T3 Platform Switched		Yellow			100386	Type 15	
BIOTEC SRL	EN				All	100436	Type 30	
	ER				All	100353	Type 01	
	IM				All	100431	Type 27	
	IR				All	100440	Type 32	
	KR				3.0	100455	Type 38	
	KW				3.6	100440	Type 32	
	SR				All	100378	Type 06	
TWO				All	100440	Type 32		
BIOTECH DENTAL	Kontakt K30				3.0	100581	Type 75	New
	Kontakt				3.6	100582	Type 76	New
	Kontakt				4.2	100582	Type 76	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
BIOTECH DENTAL	Kontakt			4.8		100582	Type 76	New
	Kontakt			5.4		100582	Type 76	New
	Weego.FR			3.6		100583	Type 78	New
	Weego.FR			4.2		100583	Type 78	New
	Weego.FR			4.8		100583	Type 78	New
BONE SYSTEM	2P			5.0		100447	Type 36	
BRAINBASE			C			100478	Type 47	
	Mytis Arrow Implants		B			100350	Type 04	
	Mytis Arrow Implants		BW			100350	Type 04	
	Mytis Arrow Implants		E			100350	Type 04	
	Mytis Arrow Implants		EW			100350	Type 04	
BREDDENT	blueSKY					100480	Type 49	
BTI BIOTECHNOLOGY INSTITUTE	Externa		Universal			100353	Type 01	
	Externa		Universal plus			100353	Type 01	
	Interna		Universal			100431	Type 27	
	Interna		Universal plus			100431	Type 27	
	Interna		Wide			100431	Type 27	
	Tiny				2.5	100368	Type 09	
	Tiny				3.75	100368	Type 09	
BTLOCK	CV1				3.75	100480	Type 49	
	CV1				4.5	100480	Type 49	
	CV1				5.5	100480	Type 49	
	CV1				6.5	100480	Type 49	
	CV2				3.75	100480	Type 49	
	CV2				4.5	100480	Type 49	
	CV2				5.5	100480	Type 49	
	CV2				6.5	100480	Type 49	
	CV3				3.75	100480	Type 49	
	CV3				4.5	100480	Type 49	
	CV3				5.5	100480	Type 49	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
BTLOCK	CV3			6.5		100480	Type 49	
CAMLOG	Camlog Cylinder-line			3.3		100432	Type 28	
	Camlog Cylinder-line			3.8		100411	Type 23	
	Camlog Cylinder-line			4.3		100411	Type 23	
	Camlog Cylinder-line			5.0		100412	Type 24	
	Camlog Cylinder-line			6.0		100412	Type 24	
	Camlog Root-line Promote			3.8		100411	Type 23	
	Camlog Root-line Promote			4.3		100411	Type 23	
	Camlog Root-line Promote			5.0		100412	Type 24	
	Camlog Root-line Promote			6.0		100412	Type 24	
	Camlog Screw-line			3.3		100432	Type 28	
	Camlog Screw-line			3.8		100411	Type 23	
	Camlog Screw-line			4.3		100411	Type 23	
	Camlog Screw-line			5.0		100412	Type 24	
	Camlog Screw-line			6.0		100412	Type 24	
	Conelog Screw-line			3.3		100557	Type 65	
	Conelog Screw-line			3.8		100541	Type 58	
	Conelog Screw-line			4.3		100541	Type 58	
	Conelog Screw-line			5.0		100380	Type 07	
	iSy (direct to implant)			All		100497	Type 57	
	iSy (to implant base)			All		100446	Type 35	
CHANG GUNG	ComMed			3.5		100481	Type 50	
	ComMed			4.0		100481	Type 50	
	ComMed			4.5		100380	Type 07	
	ComMed			5.0		100380	Type 07	
	ComMed			5.5		100380	Type 07	
	ComMed			6.0		100380	Type 07	
CHAORUM	External			3.3		100353	Type 01	
	Internal			4.0		100478	Type 47	
	Internal			4.8		100478	Type 47	
COCHLEAR	Baha & Vistafix Flange Fixture Hexagonal Connection			3.0		100368	Type 09	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
COCHLEAR	Baha & Vistafix Flange Fixture Hexagonal Connection				4.0		100368	Type 09	
	Baha BI 300				3.0		100436	Type 30	
	Baha BI 300				4.0		100436	Type 30	
	Baha BIA 300				6.0		100487	Type 55	
	Baha BIA 300				9.0		100487	Type 55	
	Baha BIA 400				6.0		100487	Type 55	
	Baha BIA 400				8.0		100487	Type 55	
	Baha BIA 400				10.0		100487	Type 55	
	Baha BIA 400				12.0		100487	Type 55	
	Baha BIA 400				14.0		100487	Type 55	
	Vistafix VXI 300				3.0		100436	Type 30	
	Vistafix VXI 300				4.0		100436	Type 30	
COMPRESS					4.0		100355	Type 02	
CORTEX	Classix				3.3		100425	Type 26	
	Classix				3.8		100425	Type 26	
	Classix				4.2		100425	Type 26	
	Classix				5.0		100425	Type 26	
	Classix				6.0		100425	Type 26	
	Dynamix				3.3		100425	Type 26	
	Dynamix				3.8		100425	Type 26	
	Dynamix				4.2		100425	Type 26	
	Dynamix				5.0		100425	Type 26	
	Dynamix				6.0		100425	Type 26	
	Saturn				3.3		100425	Type 26	
	Saturn				3.8		100425	Type 26	
	Saturn				4.2		100425	Type 26	
	Saturn				5.0		100425	Type 26	
	Saturn				6.0		100425	Type 26	
COWELLMEDI	Atlas External				3.3		100355	Type 02	
	Atlas External				4.1		100353	Type 01	
	Atlas External				5.0		100353	Type 01	
	Atlas Submerged				4.0		100402	Type 21	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
C-TECH	EL				3.5		100425	Type 26
	EL				4.3		100425	Type 26
	EL				5.1		100425	Type 26
DENTAL TECH			WP				100478	Type 47
DENTAURUM	Tiologic		S				100472	Type 45
	Tiologic		M				100472	Type 45
	Tiologic		L				100472	Type 45
DENTEGRIS	Straight & Tapered				3.3		100425	Type 26
	Straight & Tapered				3.75		100440	Type 32
	Straight & Tapered				4.5		100440	Type 32
	Straight & Tapered				4.75		100440	Type 32
	Straight & Tapered				5.5		100478	Type 47
DENTIS DENTAL	Internal				3.7		100350	Type 04
	Internal				4.3		100350	Type 04
	Internal				4.8		100350	Type 04
	OneQ S-clean Narrow (SLA surface/Submerged)				3.0		100404	Type 22
	OneQ S-clean Narrow (SLA surface/Submerged)				3.3		100404	Type 22
	OneQ S-clean Regular (SLA surface/Submerged)				3.7		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.2		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				4.7		100372	Type 11
	OneQ S-clean Regular (SLA surface/Submerged)				5.2		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				6.0		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				7.0		100372	Type 11
	OneQ S-clean Wide (SLA surface/Submerged)				8.0		100372	Type 11
	Save Submerged				5.5		100380	Type 07
	Save Submerged				6.0		100380	Type 07
	S-clean Mini (Submerged)				3.7		100372	Type 11
	S-clean Regular (Submerged)				4.1		100372	Type 11
S-clean Regular (Submerged)				4.3		100372	Type 11	
S-clean Wide (Submerged)				4.8		100372	Type 11	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
DENTIUM	Genoss				3.0		100554	Type 63
	Genoss				3.6		100402	Type 21
	Genoss				4.3		100402	Type 21
	Genoss				5.0		100402	Type 21
	Implantium		Yellow 3.6		3.4		100378	Type 06
	Implantium		Green 4.0		3.8		100378	Type 06
	Implantium		Blue 4.5		4.3		100378	Type 06
	Simple Line				3.8		100380	Type 07
	Simple Line				4.3		100380	Type 07
	Simple Line				4.8		100380	Type 07
	Super Line				4.0		100380	Type 07
	Super Line				7.0		100380	Type 07
DENTSPLY IMPLANTS	ANKYLOS C/X Implants		A				100388	Type 16
	ANKYLOS C/X Implants		B				100388	Type 16
	ANKYLOS C/X Implants		C				100388	Type 16
	ASTRA TECH Implant system MicroThread				3.5		100376	Type 05
	ASTRA TECH Implant system MicroThread				4.0		100378	Type 06
	ASTRA TECH Implant system MicroThread				4.5		100380	Type 07
	ASTRA TECH Implant system MicroThread				5.0		100380	Type 07
	ASTRA TECH Implant system OsseoSpeed Aqua				3.5 S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed Aqua				4.0 S		100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25212	Green		3.0 S	8	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25213	Green		3.0 S	9	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25214	Green		3.0 S	11	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25215	Green		3.0 S	13	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25216	Green		3.0 S	15	100481	Type 50
	ASTRA TECH Implant system OsseoSpeed EV	25221	Purple		3.6 S	6	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25222	Purple		3.6 S	8	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25223	Purple		3.6 S	9	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25224	Purple		3.6 S	11	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25225	Purple		3.6 S	13	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25226	Purple		3.6 S	15	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25227	Purple		3.6 S	17	100455	Type 38
	ASTRA TECH Implant system OsseoSpeed EV	25231	Yellow		4.2 S	6	100480	Type 49

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
DENTSPLY IMPLANTS	ASTRA TECH Implant system OsseoSpeed EV	25232	Yellow	4.2 S	8	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25233	Yellow	4.2 S	9	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25234	Yellow	4.2 S	11	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25235	Yellow	4.2 S	13	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25236	Yellow	4.2 S	15	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25237	Yellow	4.2 S	17	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25241	Blue	4.8 S	6	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25242	Blue	4.8 S	8	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25243	Blue	4.8 S	9	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25244	Blue	4.8 S	11	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25245	Blue	4.8 S	13	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25246	Blue	4.8 S	15	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25247	Blue	4.8 S	17	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25251	Brown	5.4 S	6	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25252	Brown	5.4 S	8	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25253	Brown	5.4 S	9	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25254	Brown	5.4 S	11	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25255	Brown	5.4 S	13	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25256	Brown	5.4 S	15	100462	Type 40	
	ASTRA TECH Implant system OsseoSpeed EV	25262	Yellow	4.2 C	8	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25263	Yellow	4.2 C	9	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25264	Yellow	4.2 C	11	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25265	Yellow	4.2 C	13	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25266	Yellow	4.2 C	15	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25267	Yellow	4.2 C	17	100480	Type 49	
	ASTRA TECH Implant system OsseoSpeed EV	25272	Blue	4.8 C	8	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25273	Blue	4.8 C	9	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25274	Blue	4.8 C	11	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25275	Blue	4.8 C	13	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25276	Blue	4.8 C	15	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed EV	25277	Blue	4.8 C	17	100478	Type 47	
	ASTRA TECH Implant system OsseoSpeed Lilac			4.5		100380	Type 07	
	ASTRA TECH Implant system OsseoSpeed Lilac			5.0		100380	Type 07	
	ASTRA TECH Implant system OsseoSpeed Lilac			5.0S		100380	Type 07	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25452	Yellow	4.2 PC	8	100370	Type 10	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
DENTSPLY IMPLANTS	ASTRA TECH Implant system OsseoSpeed Profile EV	25453	Yellow	4.2 PC	9	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25454	Yellow	4.2 PC	11	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25455	Yellow	4.2 PC	13	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25456	Yellow	4.2 PC	15	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25457	Yellow	4.2 PC	17	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25458	Blue	4.8 PC	8	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25459	Blue	4.8 PC	9	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25460	Blue	4.8 PC	11	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25461	Blue	4.8 PC	13	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25462	Blue	4.8 PC	15	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25463	Blue	4.8 PC	17	100440	Type 32	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25440	Yellow	4.2 PS	8	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25441	Yellow	4.2 PS	9	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25442	Yellow	4.2 PS	11	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25443	Yellow	4.2 PS	13	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25444	Yellow	4.2 PS	15	100370	Type 10	
	ASTRA TECH Implant system OsseoSpeed Profile EV	25445	Yellow	4.2 PS	17	100370	Type 10	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25446	Blue	4.8 PS	8	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25447	Blue	4.8 PS	9	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25448	Blue	4.8 PS	11	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25449	Blue	4.8 PS	13	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25450	Blue	4.8 PS	15	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed Profile EV	25451	Blue	4.8 PS	17	100440	Type 32	New
	ASTRA TECH Implant system OsseoSpeed TX Aqua			3.5S		100455	Type 38	New
	ASTRA TECH Implant system OsseoSpeed TX Aqua			4.0S		100455	Type 38	
	ASTRA TECH Implant system OsseoSpeed TX Lilac			4.5		100380	Type 07	New
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0		100380	Type 07	New
	ASTRA TECH Implant system OsseoSpeed TX Lilac			5.0 S		100380	Type 07	New
	ASTRA TECH Implant system OsseoSpeed TX Profile			4.5		100378	Type 06	New
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0		100378	Type 06	New
	ASTRA TECH Implant system OsseoSpeed TX Profile			5.0 S		100378	Type 06	New
	ASTRA TECH Implant system OsseoSpeed TX Yellow			3.0 S		100481	Type 50	New
	ASTRA TECH Implant system OsseoSpeed Yellow			3.0 S		100481	Type 50	
	XiVe		TG			100366	Type 08	
	XiVe			3.0		100404	Type 22	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
DENTSPLY IMPLANTS	XiVe				3.4		100402	Type 21
	XiVe				3.8		100472	Type 45
	XiVe				4.5		100472	Type 45
	XiVe				5.5		100472	Type 45
DIO	External FSN				3.3		100355	Type 02
	External FSN				3.75		100353	Type 01
	External FTN				5.0		100357	Type 03
	SM Submerged				3.8		100480	Type 49
	SM Submerged				4.5		100478	Type 47
	SM Submerged				5.3		100478	Type 47
DTI					3.5		100497	Type 57
					4.2		100380	Type 07
DYNA DENTAL	Dyna Helix DC				3.2		100431	Type 27
	Dyna Helix DC				3.6		100431	Type 27
	Dyna Helix DC				4.2		100431	Type 27
	Dyna Helix DC				5.0		100431	Type 27
	Dyna Helix ST				3.6		100431	Type 27
	Dyna Helix ST				4.2		100431	Type 27
	Dyna Helix ST				5.0		100431	Type 27
ECKERMANN	All Spiral						100390	Type 17
	Hexagono Ext Hex						100353	Type 01
	Tripto						100480	Type 49
EMFILS	Colosso IH ECI Evolution				3.3		100372	Type 11
	Colosso IH ECI Evolution				4.0		100372	Type 11
	Colosso IH ECI Evolution				5.0		100372	Type 11
	ICM-4010				4.0		100484	Type 53
	PAT5-4010				4.0		100484	Type 53
ENTIFIC	Flange Fixture						100353	Type 01

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
EUROTEKNIKA	Aesthetica+2		3.6		4.2	100544	Type 61	
	Aesthetica+2		3.6		4.8	100544	Type 61	
	Aesthetica+2		4.1		4.2	100544	Type 61	
	Aesthetica+2		4.1		4.8	100544	Type 61	
	Aesthetica+2		4.8		4.8	100544	Type 61	
	Aesthetica+2		4.8		6.5	100544	Type 61	
	Naturall				3.5	100455	Type 38	
	Naturall				4.0	100378	Type 06	
	Naturall				4.5	100378	Type 06	
	Naturall				5.0	100378	Type 06	
	Uneva+				3.6	100353	Type 01	
	Uneva+				4.1	100353	Type 01	
	Uneva+				4.8	100353	Type 01	
	Uneva+				6.0	100353	Type 01	
FAIR IMPLANT	FairTwo				3.5	100484	Type 53	
	FairTwo				4.2	100484	Type 53	
	FairTwo				5.0	100378	Type 06	
	FairTwo				6.0	100378	Type 06	
GC IMPLANT	Aadva		Narrow			100455	Type 38	
	Aadva		Regular			100541	Type 58	
	Aadva		Wide			100541	Type 58	
	External Plus				3.0	100402	Type 21	
	External Plus				3.4	100402	Type 21	
	External Plus				3.8	100353	Type 01	
	External Plus				5.0	100353	Type 01	
	GENESIO					100402	Type 21	
	Internal Plus				3.0	100402	Type 21	
	Internal Plus				3.4	100402	Type 21	
	SETIO					100353	Type 01	
GHIMAS	PFON EVO SERIES				3.6	100455	Type 38	
	PFON EVO SERIES				4.0	100455	Type 38	
	PFON EVO SERIES				5.0	100380	Type 07	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
HIOSSSEN	ETIII		Mini		3.5		100376	Type 05	
	ETIII				4.0		100478	Type 47	
	ETIII				4.5		100478	Type 47	
	ETIII				5.0		100478	Type 47	
	ETIII		UW		6.0		100378	Type 06	
	ETIII		UW		7.0		100378	Type 06	
	HGII		Mini				100376	Type 05	
	HGII		Regular				100378	Type 06	
	HGII		Wide				100378	Type 06	
	HGII		Ultra Wide				100378	Type 06	
	HGIII		Mini				100376	Type 05	
	HSII		Regular				100350	Type 04	
	HSII		Wide				100350	Type 04	
	IBS INNO BIO SURG	Magicore				5.5		100455	Type 38
IDI EVOLUTION	Vite				3.75		100398	Type 18	
	Vite				4.0		100398	Type 18	
IMMEDIATELOAD	IM40 Immediate Conico				4.0		100480	Type 49	
IMPLANT DIRECT	Legacy1		Blue				100455	Type 38	New
	Legacy1		Green				100440	Type 32	New
	Legacy1		Purple				100431	Type 27	New
	Legacy1		Yellow				100442	Type 33	New
	Legacy2		Blue				100455	Type 38	New
	Legacy2		Green				100440	Type 32	New
	Legacy2		Purple				100431	Type 27	New
	Legacy2		Yellow				100442	Type 33	New
	Legacy3		Blue				100455	Type 38	
	Legacy3		Green				100440	Type 32	New
	Legacy3		Purple				100431	Type 27	New
	Legacy3		Yellow				100442	Type 33	New
	Legacy4		Blue				100455	Type 38	New
	Legacy4		Green				100440	Type 32	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
IMPLANT DIRECT	Legacy4		Purple			100431	Type 27	New
	Legacy4		Yellow			100442	Type 33	New
	ScrewPlant				3.7	100425	Type 26	
	ScrewPlant				4.7	100431	Type 27	
	ScrewPlant				5.7	100431	Type 27	
IMTEC	Hexed-head				3.75	100353	Type 01	
	Hexed-head				4.75	100353	Type 01	
INTRA-LOCK SYSTEM	Internal		SQ			100378	Type 06	
	Internal		Standard			100478	Type 47	
	Internal		Wide			100478	Type 47	
IRES	IMAX NHSE				3.3	100355	Type 02	
	IMAX NHSE				3.7	100353	Type 01	
	IMAX NHSE				4.1	100353	Type 01	
	IMAX NHSE				5.2	100353	Type 01	
	IMAX NHSI				3.7	100431	Type 27	
	IMAX NHSI				4.1	100431	Type 27	
	IMAX NHSI				5.2	100431	Type 27	
	Shapeone S1B				3.1	100431	Type 27	
	Shapeone S1B				4.1	100431	Type 27	
	Shapeone S1B				4.7	100431	Type 27	
	Shapeone S1EH				3.1	100353	Type 01	
	Shapeone S1EH				4.1	100353	Type 01	
	Shapeone S1EH				4.7	100353	Type 01	
	Shapeone S1T				3.7	100350	Type 04	
	Shapeone S1T				4.1	100350	Type 04	
	Shapeone S1T				4.7	100350	Type 04	
	Shapeone S1TN				3.7	100431	Type 27	
	Shapeone S1TN				4.1	100431	Type 27	
	Shapeone S1TN				4.7	100431	Type 27	
	Volution SVB				3.7	100431	Type 27	
Volution SVB				4.1	100431	Type 27		
Volution SVB				4.7	100431	Type 27		

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
ITERUM	Conical			3.75		100425	Type 26	
JDENTALCARE	JDEvolution			3.2		100455	Type 38	
	JDEvolution			3.7		100440	Type 32	
	JDEvolution			4.3		100440	Type 32	
	JDEvolution			5.0		100440	Type 32	
JMM	POI 2 piece			3.2		100398	Type 18	
	POI 2 piece			3.7		100398	Type 18	
	POI 2 piece			4.2		100398	Type 18	
	POIEX			3.7		100386	Type 15	
	POIEX			4.2		100480	Type 49	
	POIEX			4.7		100478	Type 47	
	POIEX			5.2		100478	Type 47	
KENTEC	SB1			3.3		100378	Type 06	
	SB2			4.2		100378	Type 06	
KEYSTONE DENTAL	Genesis			3.8		100431	Type 27	
	Genesis			4.5		100442	Type 33	
	Genesis			5.5		100442	Type 33	
	Genesis			6.5		100442	Type 33	
	Prima Connex		SD	3.5		100431	Type 27	
	Prima Connex		RD	4.1		100442	Type 33	
	Prima Connex		WD	5.0		100442	Type 33	
	Renova RDI			4.5		100431	Type 27	
	Renova SDI			3.75		100425	Type 26	
	Restore			3.3		100436	Type 30	
	Restore			3.75		100353	Type 01	
	Restore			4.0		100353	Type 01	
	Restore			5.0		100357	Type 03	
	Restore			6.0		100357	Type 03	
KLOCKNER IMPLANT SYSTEM	MV VEGA Implant			3.0		100497	Type 57	
	NV VEGA Implant			3.5		100497	Type 57	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
KLOCKNER IMPLANT SYSTEM	RV VEGA Implant				4.0		100425	Type 26	
	RV VEGA Implant				4.5		100425	Type 26	
	Essential Cone				3.5		100480	Type 49	
	Essential Cone				4.0		100480	Type 49	
	Essential Cone				4.5		100480	Type 49	
	Essential Cone				4.8		100480	Type 49	
	KL Implant NP				3.3		100355	Type 02	
	KL Implant NP				3.7		100353	Type 01	
	KL Implant NP				4.2		100353	Type 01	
	KL Implant WP				4.7		100357	Type 03	
	S3M						100436	Type 30	
	SK2						100436	Type 30	
KYOCERA	FINESIA BL		NP				100545	Type 62	
	FINESIA BL		RP				100431	Type 27	
	FINESIA BL		WP				100350	Type 04	
	FINESIA TL		RP				100350	Type 04	
	FINESIA TL		WP				100350	Type 04	
LASAK	BioniQ, QN (Q-Lock Narrow)						100404	Type 22	
	BioniQ, QR (Q-Lock Regular)						100497	Type 57	
	Impladent				2.9		100355	Type 02	
	Impladent				3.7		100447	Type 36	
	Impladent				5.0		100448	Type 37	
LEADER	Implus EH				3.3		100353	Type 01	
	Implus EH				3.75		100353	Type 01	
	Implus EH				5.0		100353	Type 01	
	Implus IH				3.3		100425	Type 26	
	Tixos IH				3.75		100425	Type 26	
	Tixos IH				4.5		100425	Type 26	
	Tixos IH				5.5		100425	Type 26	
	Tixos ITMC				3.75		100497	Type 57	
	Tixos ITMC				4.5		100497	Type 57	
	Tixos ITMC				5.5		100497	Type 57	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
MEDENTIKA	M-Implant®				3.5		100455	Type 38
	M-Implant®				4.0		100455	Type 38
	M-Implant®				4.5		100455	Type 38
	M-Implant®				5.0		100455	Type 38
MEDICAL INSTINCT	Bone Trust				3.4		100466	Type 44
	Bone Trust				4.0		100466	Type 44
	Bone Trust				5.0		100466	Type 44
	Bone Trust Plus Cone				3.4		100466	Type 44
	Bone Trust Plus Cone				4.0		100466	Type 44
MEDIGMA TECHNOLOGIES	Fix 1						100390	Type 17
	Fix-A-Dent						100480	Type 49
	FixTite						100480	Type 49
	FixTite-S						100480	Type 49
MEGAGEN	AnyOne				3.5		100380	Type 07
	AnyOne				4.0		100380	Type 07
	AnyOne External				4.5		100350	Type 04
	AnyRidge						100431	Type 27
	EZ Plus Internal		Small				100455	Type 38
	EZ Plus Internal		Regular				100380	Type 07
	EZ Plus Internal		Wide				100380	Type 07
	EZ Plus Internal				4.0		100380	Type 07
	EZ Plus Internal				4.5		100380	Type 07
	EZ Plus Internal				5.0		100380	Type 07
	Mini Internal				3.0		100481	Type 50
	Resque External				7.0		100357	Type 03
MICRODENT	Duraplant 2.2				4.1		100370	Type 10
	Microdent system						100353	Type 01
	MK				5.6		100357	Type 03
	MS				3.5		100436	Type 30
	MST				3.5		100436	Type 30
	Trylogic				3.8		100431	Type 27

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
MICRODENT	Trylogic			4.3		100431	Type 27	
	Universal			4.1		100353	Type 01	
	Universal			5.1		100353	Type 01	
MIS	C1 Conical Connection		NP			100455	Type 38	New
	C1 Conical Connection		SP			100480	Type 49	
	C1 Conical Connection		WP			100584	Type 79	New
	V3 Conical Connection		NP			100455	Type 38	New
	V3 Conical Connection		SP			100480	Type 49	New
	Internal Hexagon (seven, biocom)			3.75		100440	Type 32	
	Internal Octagon (mistral)			4.75		100350	Type 04	
	Seven		Narrow	3.3		100554	Type 63	
MOZO-GRAU	MG-Inhex		Standard			100455	Type 38	
	MG-Inhex		Maxi			100380	Type 07	
	MG-Osseous		Mini			100353	Type 01	
	MG-Osseous		Standard			100353	Type 01	
	MG-Osseous		Maxi			100357	Type 03	
ND IMPLANTS	Niko			3.6		100431	Type 27	New
	Niko			4.6		100380	Type 07	New
NEMRIS	Aesthura Classic			3.75		100446	Type 35	
	Aesthura Classic			4.75		100446	Type 35	
NEOBIOTECH	IS Fixture (RBM Surface)			4.0		100376	Type 05	
	IS Fixture (RBM Surface)			4.5		100376	Type 05	
	IS Fixture (RBM Surface)			4.8		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			5.5		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			6.0		100376	Type 05	
	IS S-Wide Fixture (RBM Surface)			7.0		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			3.5		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			4.0		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			4.5		100376	Type 05	
	IS-II active Fixture (S.L.A. Surface)			5.0		100376	Type 05	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
NEOBIOTECH	IS-II active S-Wide Fixture (S.L.A. Surface)			5.5		100376	Type 05	
	IS-II active S-Wide Fixture (S.L.A. Surface)			6.0		100376	Type 05	
	IS-II active S-Wide Fixture (S.L.A. Surface)			7.0		100376	Type 05	
	IS-II Fixture (RBM Surface)			3.5		100376	Type 05	
	IS-II Fixture (RBM Surface)			4.0		100376	Type 05	
	IS-II Fixture (RBM Surface)			4.5		100376	Type 05	
	IS-II Fixture (RBM Surface)			5.0		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			5.5		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			6.0		100376	Type 05	
	IS-II S-Wide Fixture (RBM Surface)			7.0		100376	Type 05	
	IS-III active Fixture			3.5		100376	Type 05	
	IS-III active Fixture			4.0		100376	Type 05	
	IS-III active Fixture			4.5		100376	Type 05	
	IS-III active Fixture			5.0		100376	Type 05	
	IS-III active S-Wide Fixture			5.5		100376	Type 05	
	IS-III active S-Wide Fixture			6.0		100376	Type 05	
	IS-III active S-Wide Fixture			7.0		100376	Type 05	
	Neoplant External			4.0		100353	Type 01	
	Neoplant External			5.0		100353	Type 01	
NEODENT	Alvim CM			3.5		100388	Type 16	
	Alvim CM			4.3		100388	Type 16	
	Alvim CM			5.0		100388	Type 16	
	Alvim II Plus			3.75		100402	Type 21	
	Alvim II Plus			4.3		100402	Type 21	
	Alvim II Plus			5.0		100402	Type 21	
	Alvim Ti			4.3		100353	Type 01	
	Drive CM			4.3		100479	Type 48	
	Drive CM			5.0		100479	Type 48	
	Titamax CM			3.5		100388	Type 16	
	Titamax CM			3.75		100388	Type 16	
	Titamax CM			4.0		100388	Type 16	
	Titamax CM			5.0		100388	Type 16	
	Titamax GT			3.7		100368	Type 09	
	Titamax GT			4.0		100368	Type 09	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
NEODENT	Titamax GT				4.5		100368	Type 09	
	Titamax II Plus				3.75		100402	Type 21	
	Titamax II Plus				5.0		100402	Type 21	
	Titamax Ti				3.75		100353	Type 01	
	Titamax Ti				4.0		100353	Type 01	
	Titamax Ti				4.5		100357	Type 03	
	Titamax Ti				5.0		100357	Type 03	
NEOSS	Neoss Implant				3.25		100455	Type 38	
	ProActive Straight				3.5		100398	Type 18	
	ProActive Straight				4.0		100398	Type 18	
	ProActive Straight				4.5		100398	Type 18	
	ProActive Straight				5.0		100398	Type 18	
	ProActive Straight				5.5		100398	Type 18	
	ProActive Straight				6.0		100398	Type 18	
	ProActive Tapered				3.5		100398	Type 18	
	ProActive Tapered				4.0		100398	Type 18	
	ProActive Tapered				4.5		100398	Type 18	
	ProActive Tapered				5.0		100398	Type 18	
	ProActive Tapered				5.5		100398	Type 18	
	ProActive Tapered				6.0		100398	Type 18	
NIPPON PISTON RING	1 Stage Screw RH				3.3		100398	Type 18	
	1 Stage Screw RH				4.0		100398	Type 18	
	1 Stage Screw RL				3.3		100398	Type 18	
	1 Stage Screw RL				4.0		100398	Type 18	
	1 Stage Screw WH				5.0		100398	Type 18	
	1 Stage Screw WL				5.0		100398	Type 18	
	2 Stage Cylinder				3.3		100398	Type 18	
	2 Stage Cylinder				4.0		100398	Type 18	
	2 Stage Cylinder				5.0		100398	Type 18	
	2 Stage Screw				3.3		100398	Type 18	
	2 Stage Screw				4.0		100398	Type 18	
	2 Stage Screw				5.0		100398	Type 18	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
NOBEL BIOCARE	Brånemark system		NP			100355	Type 02	
	Brånemark system		RP			100353	Type 01	
	Brånemark system		WP			100357	Type 03	
	Brånemark system Novum					100353	Type 01	
	Brånemark system Shorty		RP			100353	Type 01	
	Brånemark system Shorty		WP			100357	Type 03	
	NobelActive		NP	3.5		100543	Type 60	
	NobelActive		RP	4.3		100544	Type 61	
	NobelActive		RP	5.0		100544	Type 61	
	NobelActive		WP	5.5		100350	Type 04	
	NobelActive			3.0		100545	Type 62	
	NobelDirect		NP			100366	Type 08	
	NobelDirect		RP			100368	Type 09	
	NobelDirect		WP			100368	Type 09	
	NobelParallel Conical Connection		NP	3.75		100543	Type 60	
	NobelParallel Conical Connection		RP	4.3		100544	Type 61	
	NobelParallel Conical Connection		RP	5.0		100544	Type 61	
	NobelParallel Conical Connection		WP	5.5		100350	Type 04	
	NobelReplace Conical Connection		NP	3.5		100543	Type 60	
	NobelReplace Conical Connection		RP	4.3		100544	Type 61	
	NobelReplace Conical Connection		RP	5.0		100544	Type 61	
	NobelReplace Conical Connection PMC		NP			100543	Type 60	
	NobelReplace Conical Connection PMC		RP			100544	Type 61	
	NobelReplace Platform Shift		NP			100359	Type 12	
	NobelReplace Platform Shift		RP			100361	Type 13	
	NobelReplace Platform Shift		WP			100361	Type 13	
	NobelReplace Straight Groovy		NP			100359	Type 12	
	NobelReplace Straight Groovy		RP			100361	Type 13	
	NobelReplace Straight Groovy		WP			100361	Type 13	
	NobelReplace Tapered Groovy		NP			100359	Type 12	
	NobelReplace Tapered Groovy		RP			100361	Type 13	
	NobelReplace Tapered Groovy		WP			100361	Type 13	
	NobelReplace Tapered Groovy			6.0		100462	Type 40	
	NobelSpeedy Groovy		NP			100355	Type 02	
	NobelSpeedy Groovy		RP			100353	Type 01	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
NOBEL BIOCARE	NobelSpeedy Groovy		WP			100357	Type 03	
	NobelSpeedy Replace		NP			100359	Type 12	
	NobelSpeedy Replace		RP			100361	Type 13	
	NobelSpeedy Replace		WP			100361	Type 13	
	NobelSpeedy Replace			6.0		100462	Type 40	
	NobelSpeedy Shorty		RP			100353	Type 01	
	NobelSpeedy Shorty		WP			100357	Type 03	
	Replace Select Straight		NP			100359	Type 12	
	Replace Select Straight		RP			100361	Type 13	
	Replace Select Straight		WP			100361	Type 13	
	Replace Select Tapered		NP			100359	Type 12	
	Replace Select Tapered		RP			100361	Type 13	
	Replace Select Tapered		WP			100361	Type 13	
	Replace Select Tapered			6.0		100462	Type 40	
	Replace Select Tapered PMC		NP			100359	Type 12	
	Replace Select Tapered PMC		RP			100361	Type 13	
	Replace Select Tapered PMC		WP			100361	Type 13	
	Replace Select Tapered PMC			6.0		100462	Type 40	
	Replace Select TC		NP			100359	Type 12	
	Replace Select TC		RP			100361	Type 13	
NUCLEOSS	Nucleoss T3			3.4		100376	Type 05	
	Nucleoss T3			3.8		100376	Type 05	
	Nucleoss T3			4.2		100376	Type 05	
	Nucleoss T3			5.0		100376	Type 05	
	Nucleoss T4			3.4		100376	Type 05	
	Nucleoss T4			3.8		100376	Type 05	
	Nucleoss T4			4.2		100376	Type 05	
	Nucleoss T5			5.0		100380	Type 07	
OCO BIOMEDICAL	Engage		3.0		3.2	100455	Type 38	
	Engage		3.5		4.0	100425	Type 26	
	Engage		4.5		5.0	100425	Type 26	
	Engage		4.5		6.0	100431	Type 27	
	ERI				4.0	100380	Type 07	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
OCO BIOMEDICAL	ERI				5.0		100380	Type 07	
	Macro				6.0		100478	Type 47	
	TSI				3.25		100446	Type 35	
	TSI				4.0		100380	Type 07	
	TSI				5.0		100380	Type 07	
OSSEOLINK					3.5		100402	Type 21	
					4.0		100402	Type 21	
					5.0		100402	Type 21	
OSSTEM	GSII		Mini				100376	Type 05	
	GSII		Regular				100378	Type 06	
	GSII		Wide				100378	Type 06	
	GSII		Ultra Wide				100378	Type 06	
	GSIII		Mini				100376	Type 05	
	GSIII		Regular				100378	Type 06	
	SSI		Regular				100350	Type 04	
	SSI		Wide				100350	Type 04	
	SSII		Regular				100350	Type 04	
	SSII		Wide				100350	Type 04	
	SSIII		Regular				100350	Type 04	
	SSIII		Wide				100350	Type 04	
	TSIII CA				3.5		100376	Type 05	
	TSIII SA				3.0		100372	Type 11	
	TSIII SA				4.0		100478	Type 47	
	TSIII SA				4.5		100478	Type 47	
	TSIII SA				5.0		100478	Type 47	
	TSIII SA Fixture – NoMount				4.5		100478	Type 47	
	TSIII SA Ultra-Wide				6.0		100350	Type 04	
	TSIII SA Ultra-Wide				7.0		100350	Type 04	
	USII			Mini				100355	Type 02
	USII			Regular				100353	Type 01
	USII			Wide				100357	Type 03
	USII			Ultra Wide				100357	Type 03
USIII			Mini				100355	Type 02	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
OSSTEM	USIII		Regular			100353	Type 01	
	USIII		Wide			100357	Type 03	
	USIII		Ultra Wide			100357	Type 03	
	USIV Conical		Mini			100366	Type 08	
OSTEOPANT	Two-piece					100370	Type 10	
OSTEO-TI	Combi Dental Implant			3.25		100425	Type 26	
	Combi Dental Implant			3.75		100425	Type 26	
	Combi Dental Implant			4.5		100425	Type 26	
	Combi Dental Implant			5.5		100425	Type 26	
OTICON MEDICAL	Ponto Implant			3.0		100368	Type 09	
	Ponto Implant			4.0		100368	Type 09	
	Ponto Implant		Wide			100368	Type 09	
PALTOP	PALTOP		NP			100455	Type 38	
	PALTOP		SP			100480	Type 49	
PHIBO	Aurea EVO		NP			100554	Type 63	
	Aurea EVO		RP			100425	Type 26	
	Aurea EVO		WP			100425	Type 26	
	Aurea		NP			100554	Type 63	
	Aurea		RP			100431	Type 27	
	Aurea		WP			100440	Type 32	
	BNT				S3	100353	Type 01	
	BNT				S4	100353	Type 01	
	BNT				S5	100353	Type 01	
	TSA				S3	100455	Type 38	
	TSA				S4	100431	Type 27	
	TSA				S5	100431	Type 27	
	TSA Advance				S3	100455	Type 38	
	TSA Advance				S4	100431	Type 27	
	TSA Advance				S5	100431	Type 27	
	TSH				S3	100353	Type 01	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
PHIBO	TSH			54		100353	Type 01	
	TSH			55		100353	Type 01	
P-I BRÄNEMARK PHILOSOPHY	Internal Grip Amplified			3.5		100554	Type 63	
	Internal Grip Amplified			4.3		100482	Type 51	
	Internal Grip Amplified			5.1		100482	Type 51	
	Internal Grip External Hex			3.5		100455	Type 38	
	Internal Grip External Hex			3.75		100398	Type 18	
	Internal Grip External Hex			4.1		100398	Type 18	
	Morse Taper			3.3		100402	Type 21	
	Morse Taper			3.75		100402	Type 21	
	Morse Taper			4.8		100402	Type 21	
PLAN 1 HEALTH	GTB			3.3		100497	Type 57	
	GTB			3.6		100497	Type 57	
	GTB			4.3		100497	Type 57	
	Hexcel EHC, External Hex					100353	Type 01	
	Hexcel IHC, Internal Hex					100425	Type 26	
	Hexcel IHT, Internal Hex					100425	Type 26	
PLATON	EL Implant Narrow			3.7		100543	Type 60	
	EL Implant Regular			5.1		100478	Type 47	
	Type I					100350	Type 04	
	Type II					100350	Type 04	
	Type IV Pro					100478	Type 47	
PRODENT ITALIA	Perform SM			3.6		100480	Type 49	
	Perform SM			4.0		100480	Type 49	
	Perform SM			4.5		100380	Type 07	
	Perform SM			5.0		100380	Type 07	
PW PLUS	Dif			3.8		100480	Type 49	New
	Dif			4.2		100380	Type 07	New
REVIOS				3.8		100425	Type 26	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
REVIOS				4.3		100425	Type 26	
				5.0		100425	Type 26	
REX IMPLANTS	Rex TL			1.8		100350	Type 04	
	Rex TL			2.9		100480	Type 49	
S.I.N IMPLANT SYSTEM	External Hex			3.4		100353	Type 01	
	External Hex			4.1		100353	Type 01	
	External Hex			5.0		100353	Type 01	
	Internal Hex			3.8		100386	Type 15	
	Internal Hex			4.5		100386	Type 15	
	Internal Hex			5.5		100386	Type 15	
	Strong SW HI			4.5		100479	Type 48	
	Strong SW Morse			4.5		100484	Type 53	
SDS	SDS 2.0			3.8		100589	Type 69	
SERVO-DENTAL	Trias			3.8		100378	Type 06	
SHINHUNG	Luna		NP	3.5		100376	Type 05	
	Luna		RP	4.0		100478	Type 47	
	Luna		RP	4.5		100478	Type 47	
	Luna		RP	5.0		100478	Type 47	
	Sola		NP	3.5		100355	Type 02	
	Sola		RP	4.0		100353	Type 01	
	Sola		RP	4.5		100353	Type 01	
	Sola		WP	5.0		100357	Type 03	
	Stella		RP	4.0		100544	Type 61	
	Stella		RP	4.5		100544	Type 61	
	Stella		WP	4.5		100544	Type 61	
	Stella		WP	5.0		100544	Type 61	
SHOFU	BFIMP30-OSTL			3.0		100545	Type 62	
	BFIMP35-OSTL			3.5		100404	Type 22	
	BFIMP40-OSTL			4.0		100543	Type 60	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
SIC IMPLANT SYSTEM	Ace				3.4		100466	Type 44	
	Ace				4.0		100466	Type 44	
	Ace				5.0		100466	Type 44	
SOUTHERN IMPLANTS	External Hex						100353	Type 01	
	Extra-oral				3.75		100353	Type 01	
	IT Connection						100350	Type 04	
	MAX				7.0		100353	Type 01	
	MAX				8.0		100353	Type 01	
	MAX				9.0		100353	Type 01	
	MAX-IT				7.0		100380	Type 07	
	MAX-IT				8.0		100380	Type 07	
	MAX-IT				9.0		100380	Type 07	
	Tri-MAX				7.0		100412	Type 24	
	Tri-MAX				8.0		100462	Type 40	
	Tri-MAX				9.0		100462	Type 40	
	Tri-Nex				3.5		100359	Type 12	
	Tri-Nex				4.3		100361	Type 13	
	Tri-Nex				5.0		100361	Type 13	
STERNGOLD DENTAL	External Hexagon		RP		3.75		100353	Type 01	
	External Hexagon		RP		4.0		100353	Type 01	
	External Hexagon		RP		5.0		100353	Type 01	
	Stern IC Internal Octagon		RN		4.1		100350	Type 04	
STRAUMANN	Bone Level Implants Roxolid SLActive		NC		3.3		100484	Type 53	
	Bone Level Implants Roxolid SLActive		RC		4.1		100485	Type 54	
	Bone Level Implants Roxolid SLActive		RC		4.8		100485	Type 54	
	Bone Level Implants Titanium SLA		NC		3.3		100484	Type 53	
	Bone Level Implants Titanium SLA		RC		4.1		100485	Type 54	
	Bone Level Implants Titanium SLA		RC		4.8		100485	Type 54	
	Bone Level Tapered Implants Roxolid SLA	021.3508	NC		3.3	8	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3510	NC		3.3	10	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3512	NC		3.3	12	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3514	NC		3.3	14	100484	Type 53	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
STRAUMANN	Bone Level Tapered Implants Roxolid SLA	021.3516	NC	3.3	16	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.3518	NC	3.3	18	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLA	021.5508	RC	4.1	8	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5510	RC	4.1	10	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5512	RC	4.1	12	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5514	RC	4.1	14	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5516	RC	4.1	16	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.5518	RC	4.1	18	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7508	RC	4.8	8	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7510	RC	4.8	10	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7512	RC	4.8	12	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7514	RC	4.8	14	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7516	RC	4.8	16	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLA	021.7518	RC	4.8	18	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.3308	NC	3.3	8	100484	Type 53	
	Bone Level Tapered Implants Roxolid SLActive	021.3310	NC	3.3	10	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3312	NC	3.3	12	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3314	NC	3.3	14	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3316	NC	3.3	16	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.3318	NC	3.3	18	100484	Type 53	New
	Bone Level Tapered Implants Roxolid SLActive	021.5308	RC	4.1	8	100485	Type 54	
	Bone Level Tapered Implants Roxolid SLActive	021.5310	RC	4.1	10	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5312	RC	4.1	12	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5314	RC	4.1	14	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5316	RC	4.1	16	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.5318	RC	4.1	18	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7308	RC	4.8	8	100485	Type 54	
	Bone Level Tapered Implants Roxolid SLActive	021.7310	RC	4.8	10	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7312	RC	4.8	12	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7314	RC	4.8	14	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7316	RC	4.8	16	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive	021.7318	RC	4.8	18	100485	Type 54	New
	Bone Level Tapered Implants Roxolid SLActive		SC	2.9		100404	Type 22	
	Bone Level Tapered Implants Titanium SLA	021.3408	NC	3.3	8	100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3410	NC	3.3	10	100484	Type 53	New

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
STRAUMANN	Bone Level Tapered Implants Titanium SLA	021.3412	NC	3.3	12	100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3414	NC	3.3	14	100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3416	NC	3.3	16	100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.3418	NC	3.3	18	100484	Type 53	New
	Bone Level Tapered Implants Titanium SLA	021.5408	RC	4.1	8	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5410	RC	4.1	10	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5412	RC	4.1	12	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5414	RC	4.1	14	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5416	RC	4.1	16	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.5418	RC	4.1	18	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7408	RC	4.8	8	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7410	RC	4.8	10	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7412	RC	4.8	12	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7414	RC	4.8	14	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7416	RC	4.8	16	100485	Type 54	New
	Bone Level Tapered Implants Titanium SLA	021.7418	RC	4.8	18	100485	Type 54	New
	Extra Oral					100353	Type 01	
	Narrow Neck					100390	Type 17	
	Ortho					100368	Type 09	
	Palatal					100368	Type 09	
	Standard Implants Roxolid SLActive		RN	3.3		100350	Type 04	
	Standard Implants Roxolid SLActive		RN	4.1		100350	Type 04	
	Standard Implants Roxolid SLActive		RN	4.8		100350	Type 04	
	Standard Implants Roxolid SLActive		WN	4.8		100350	Type 04	
	Standard Implants Titanium SLA		RN	3.3		100350	Type 04	
	Standard Implants Titanium SLA		RN	4.1		100350	Type 04	
	Standard Implants Titanium SLA		RN	4.8		100350	Type 04	
	Standard Implants Titanium SLA		WN	4.8		100350	Type 04	
	Standard Plus Implants Roxolid SLActive		NNC	3.3		100560	Type 68	
	Standard Plus Implants Roxolid SLActive		RN	3.3		100350	Type 04	
	Standard Plus Implants Roxolid SLActive		RN	4.1		100350	Type 04	
	Standard Plus Implants Roxolid SLActive		RN	4.8		100350	Type 04	
	Standard Plus Implants Roxolid SLActive		WN	4.8		100350	Type 04	
	Standard Plus Implants Titanium SLA		RN	3.3		100350	Type 04	
	Standard Plus Implants Titanium SLA		RN	4.1		100350	Type 04	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
STRAUMANN	Standard Plus Implants Titanium SLA		RN		4.8		100350	Type 04
	Standard Plus Implants Titanium SLA		WN		4.8		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN		3.3		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		RN		4.1		100350	Type 04
	Tapered Effect Implants Roxolid SLActive		WN		4.8		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN		3.3		100350	Type 04
	Tapered Effect Implants Titanium SLA		RN		4.1		100350	Type 04
	Tapered Effect Implants Titanium SLA		WN		4.8		100350	Type 04
SWEDEN & MARTINA	EZ						100380	Type 07
	Global				3.8		100479	Type 48
	Global				4.3		100479	Type 48
	Global				4.8		100479	Type 48
	Global				5.5		100479	Type 48
	Kohno				3.8		100440	Type 32
	Kohno				4.25		100380	Type 07
	Kohno				5.0		100380	Type 07
	Premium - new				3.3		100440	Type 32
	Premium - new				3.8		100440	Type 32
	Premium				4.25		100380	Type 07
	Premium				5.0		100380	Type 07
	Premium - old				3.3		100402	Type 21
SYBRON	Endopore External				3.5		100436	Type 30
	Endopore External				4.1		100353	Type 01
	Endopore External				5.0		100353	Type 01
	Entegra External				3.25		100353	Type 01
	Entegra External				4.0		100353	Type 01
	Pitt-Easy				4.0		100412	Type 24
	Pitt-Easy				4.9		100412	Type 24
	SybronPRO XRT				3.3		100378	Type 06
	SybronPRO XRT				3.5		100378	Type 06
	SybronPRO XRT				4.1		100378	Type 06
TA-DENT					3.3		100455	Type 38

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
TA-DENT				4.8		100378	Type 06	
				5.6		100480	Type 49	
TAG DENTAL	Axis			3.3		100431	Type 27	New
	Axis			3.75		100431	Type 27	New
	Axis			4.2		100431	Type 27	New
	Axis			5.0		100431	Type 27	New
	Axis			6.0		100431	Type 27	New
	Massif			3.75		100431	Type 27	New
	Massif			4.2		100431	Type 27	New
	Massif			5.0		100431	Type 27	New
	Massif			6.0		100431	Type 27	New
	TBR	External Hexagon					100353	Type 01
Internal Octagon						100398	Type 18	
M-implant				3.2		100402	Type 21	
M-implant				3.9		100402	Type 21	
M-implant				4.7		100402	Type 21	
Z1-M implant				3.2		100402	Type 21	
Z1-M implant				3.9		100402	Type 21	
Z1-M implant				4.7		100402	Type 21	
TEKKA	Inkone			3.5		100480	Type 49	
	Inkone			4.5		100480	Type 49	
	Krestal			3.5		100378	Type 06	
	Krestal			4.0		100378	Type 06	
	Krestal			4.2		100378	Type 06	
	Krestal			5.1		100378	Type 06	
	Progress			3.4		100378	Type 06	
	Progress			4.0		100378	Type 06	
	T-quest			4.0		100480	Type 49	
	T-quest			4.5		100480	Type 49	
Twinkone			3.5		100390	Type 17		
TFI SYSTEM	TFI			3.3		100425	Type 26	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New
TFI SYSTEM	TFI			3.75		100425	Type 26	
	TFI			4.2		100425	Type 26	
	TFI			5.0		100425	Type 26	
THOMMEN MEDICAL	SPI System Platform			3.5		100404	Type 22	
	SPI System Platform			4.0		100402	Type 21	
	SPI System Platform			4.5		100386	Type 15	
	SPI System Platform			5.0		100386	Type 15	
	SPI System Platform			6.0		100386	Type 15	
TIDAL SPIRAL				3.75		100431	Type 27	
				4.2		100431	Type 27	
				5.0		100431	Type 27	
				6.0		100431	Type 27	
TRATE	Rootform			4.0		100497	Type 57	New
TRI DENTAL IMPLANTS	TRI-Octa 4.8			3.75		100544	Type 61	
	TRI-Octa 4.8			4.1		100544	Type 61	
	TRI-Octa 4.8			4.7		100544	Type 61	
	TRI-Vent			3.5		100440	Type 32	
U-IMPL	External Hex					100353	Type 01	
	Internal Hex					100440	Type 32	
WARANTEC	Oneplant Implant			4.3		100412	Type 24	
	Oneplant Internal Tapered					100378	Type 06	
WEGO				3.4	9	100484	Type 53	
				4.3	7	100455	Type 38	
				4.3	9	100411	Type 23	
WIELAND IMPLANT SYSTEM				3.5		100472	Type 45	
				4.3		100472	Type 45	
				5.0		100472	Type 45	

Implant Company	Model	Catalog No	Platform	Diameter (mm)	Length (mm)	Art. No	SmartPeg	New	
YAMAHACHI DENTAL	μ-one HA Implant 2 piece				3.4		100455	Type 38	
	μ-one HA Implant 2 piece				3.8		100455	Type 38	
	μ-one HA Implant 2 piece				4.3		100431	Type 27	
	μ-one HA Implant 2 piece				5.0		100398	Type 18	
ZIMMER DENTAL	AdVent				4.5		100446	Type 35	
	AdVent				5.7		100444	Type 34	
	Eztetic Implant System CT31				3.1		100543	Type 60	
	Screw-vent		3.5		3.3		100440	Type 32	
	Screw-vent		3.5		3.7		100440	Type 32	
	Screw-vent		4.5		4.7		100431	Type 27	
	Spline				3.25		100436	Type 30	
	Spline				3.75		100436	Type 30	
	Spline				5.0		100436	Type 30	
	SwissPlus				3.8		100425	Type 26	
	SwissPlus				4.8		100350	Type 04	
	Tapered Screw-vent		3.5		3.7		100440	Type 32	
	Tapered Screw-vent		3.5		4.1		100440	Type 32	
	Tapered Screw-vent		4.5		4.7		100431	Type 27	
Tapered Screw-vent		5.7		6.0		100442	Type 33		
ZUGA MEDICAL	Bone Level				3.5		100554	Type 63	
	Bone Level				5.0		100557	Type 65	

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	Art. No	SmartPeg	New
CAMLOG	Comfour Bar abutment		Yellow		3.8		100579	Type 73	New
	Comfour Bar abutment		Purple		4.3		100579	Type 73	New
	Comfour Bar abutment		Blue		5.0		100580	Type 74	New
COCHLEAR	Baha BA 210					5.5	100487	Type 55	
	Baha BA 210					8.5	100487	Type 55	
	Baha BA 300, BIA 300					6.0	100487	Type 55	
	Baha BA 300, BIA 300					9.0	100487	Type 55	
	Baha BA 400, BIA 400					10.0	100487	Type 55	
	Baha BA 400, BIA 400					12.0	100487	Type 55	
	Baha BA 400, BIA 400					14.0	100487	Type 55	
	Baha BA 400, BIA 400					6.0	100487	Type 55	
	Baha BA 400, BIA 400					8.0	100487	Type 55	
	Vistafix Standard abutment					3.0	100363	Type A3	
	Vistafix Standard abutment					4.0	100363	Type A3	
	Vistafix Standard abutment					5.5	100363	Type A3	
	Vistafix Standard abutment					7.0	100363	Type A3	
	Vistafix VXA 300					3.5	100363	Type A3	
	Vistafix VXA 300					4.5	100363	Type A3	
	Vistafix VXA 300					6.0	100363	Type A3	
Vistafix VXA 300					7.5	100363	Type A3		
NEOSS	Access abutment						100483	Type 52	
NOBEL BIOCARE	Brånemark system MUA		NP				100363	Type A3	
	Brånemark system MUA		RP				100363	Type A3	
	On1 Base CC		NP				100579	Type 73	New
	On1 Base CC		RP				100579	Type 73	New
	On1 Base CC		WP				100579	Type 73	New
	Mua		NP				100415	Type 25	New
	Mua		RP				100415	Type 25	New
	Mua		WP				100415	Type 25	New
OTICON MEDICAL	Ponto Implant abutment					4.0	100487	Type 55	

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	Art. No	SmartPeg	New
STRAUMANN	Coping for screw-retained abutment, TI Bridge	023.2749	NC		3.5		100415	Type 25	New
	Coping for screw-retained abutment, TI Bridge	023.4751	NC/RC		4.6		100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2749	NC	17	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2755	NC	17	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0010	NC	17	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2751	NC	30	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.2757	NC	30	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0012	NC	30	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4747	RC	17	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4752	RC	17	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4749	RC	30	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.4754	RC	30	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0014	RC	17	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type A	022.0016	RC	30	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2750	NC	17	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2756	NC	17	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0011	NC	17	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2752	NC	30	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.2758	NC	30	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0013	NC	30	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4748	RC	17	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4753	RC	17	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4750	RC	30	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.4755	RC	30	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0017	RC	30	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, angled, Type B	022.0015	RC	17	4.6	5.5	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2745	NC	0	3.5	1	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2746	NC	0	3.5	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2753	NC	0	3.5	4	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2747	NC	0	4.6	1	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2748	NC	0	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.2754	NC	0	4.6	4	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4745	RC	0	4.6	1	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4746	RC	0	4.6	2.5	100415	Type 25	New
	Screw-retained abutment, TAN, straight	022.4751	RC	0	4.6	4	100415	Type 25	New

Implant Company	Abutment Model	Catalog No	Platform	Angle	Diameter (mm)	Length or GH (mm)	Art. No	SmartPeg	New
THOMMEN MEDICAL	SPI System Bar abutment						100447	Type 36	
	VARIOMulti abutment						100415	Type 25	



A fresh start every time!

To avoid damage to the implant, the SmartPeg material is quite soft, which limits the lifetime of the threads. Therefore, SmartPegs are single use items.

Please visit
www.osstell.com/smartpegguide/
to get the latest SmartPeg
Reference Guide, as it is regularly
updated.

