




Zirconia Abutments

Special two-part design consists of a Titanium base and a Zirconia abutment in various designs, for maximal restoration flexibility and best results.

- **Esthetic** - due to its colour and texture the light absorbance and reflection are identical to a natural tooth. Special anatomical design provides long-term soft tissue support.
- **Strong** - superb mechanical strength, stable and durable connection between the titanium base and the implant.
- **High restoration flexibility** - two part design facilitates the customization process.

Straight Zirconia Abutment		Angled Zirconia Abutments	
Code	ZHBZ	ZHBZ-15	ZHBZ-25
Ref. No.	6054	6058	6057

Accessories for Zirconia Abutments



Material	Titanium	Titanium	Plastic	Titanium
Code	HBZ	HBZ-R	PHBZ*	SHBZ
Ref. No.	6043	6044	6080	6053
Instructions	For restoration on single implant	Bridge restoration on multiple implants. Provides easy solution for divergence of up to 30° between implants	Used for fabrication of customized Zirconia abutment using Pantograph technique	Used with Titanium bases and customized Zirconia abutment or bridges