

# Cytoplast™ Titanium-Reinforced

Titanium-reinforced, high-density PTFE membrane



## Anterior Narrow

12 mm x 24 mm

Designed for narrow single-tooth extraction sites, especially where one bony wall is missing

## Anterior Singles

14 mm x 24 mm

Designed for single-tooth extraction sites, especially where one or more bony walls are missing

## Buccal

17 mm x 25 mm

Designed for large buccal defects

## Posterior Singles

20 mm x 25 mm

Designed for posterior extraction sites and limited ridge augmentation

## Posterior Singles T<sup>2</sup>

25 mm x 36 mm

Designed for posterior extraction sites and limited ridge augmentation

## Posterior Large

25 mm x 30 mm

Designed for large bony defects, including ridge augmentation

### Ti-250 (250 µm thick)

Ti250ANL-1

Ti250ANL-2

Ti250AS-1

Ti250AS-2

Ti250BL-1

Ti250BL-2

Ti250PS-1

Ti250PS-2

Ti250PST-1

Ti250PST-2

Ti250PL-1

Ti250PL-2

### Ti-150 (150 µm thick)

Ti150ANL-1

Ti150ANL-2

Ti150AS-1

Ti150AS-2

Ti150BL-1

Ti150BL-2

Ti150PS-1

Ti150PS-2

Ti150PST-1

Ti150PST-2

Ti150PL-1

Ti150PL-2

shown actual size.

(1 membrane per box)

(2 membranes per box)

(1 membrane per box)

(2 membranes per box)

(1 membrane per box)

(2 membranes per box)

(1 membrane per box)

(2 membranes per box)

(1 membrane per box)

(2 membranes per box)

(1 membrane per box)

(2 membranes per box)

