



**Ti-150 membranes are 40% thinner than Ti-250 membranes, providing clinicians another handling option in Cytoplast™ Titanium-Reinforced Membranes.*

Posterior Large T²

30 mm x 41 mm

Designed for large bony defects, including ridge augmentation

XL

30 mm x 40 mm

Designed for very large bony defects, including ridge augmentation

XLK

30 mm x 40 mm

Designed for very large bony defects, including ridge augmentation

Anterior Perio

13 mm x 19 mm

Designed to fit periodontal defects in the anterior

Posterior Perio

13 mm x 18 mm

Designed to fit periodontal defects in the posterior

	Ti-250 <i>(250 µm thick)</i>	Ti-150 <i>(150 µm thick)</i>
	Ti250PLT-1	Ti150PLT-1
	Ti250PLT-2	Ti150PLT-2
	Ti250XL-1	Ti150XL-1
	Ti250XL-2	Ti150XL-2
	Ti250XLK-1	Ti150XLK-1
	Ti250XLK-2	Ti150XLK-2
	Ti250AP-1	Ti150AP-1
	Ti250AP-2	Ti150AP-2
	Ti250PP-1	Ti150PP-1
	Ti250PP-2	Ti150PP-2

shown actual size.

(1 membrane per box)

(2 membranes per box)



(1 membrane per box)

(2 membranes per box)



(1 membrane per box)

(2 membranes per box)



(1 membrane per box)

(2 membranes per box)



(1 membrane per box)

(2 membranes per box)

