

HU-FRIEDY NATIONAL PROMOTIONS

HUNGARY
SPRING 2019

Valid from 1st March to 31st May 2019

For more information about Hu-Friedy's Total Solutions visit [HU-FRIEDY.EU](https://www.hu-friedy.eu)

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved.

How the best perform






PRODUCTS & SERVICES

Since 1908, Hu-Friedy has been putting the highest quality and most innovative dental products in the hands of dental professionals worldwide. And today as one of the world's leading dental instrument and product companies we continue our commitment to advance dental performance by supporting these products with world class service, education and peer communities.

We're committed to making sure the best perform that way.





Visit our new **IMS INSTRUMENT
MANAGEMENT SYSTEM** site
**Organisation, Efficiency and Safety
in a simple click!**
Use the QrCode

INDEX

Cutting-Edge clinical solutions
Check out what's new in Hu-Friedy Family!

02

NEW PRODUCTS

- 02 • 3ssential Kit by Werestore.it
- 03 • Titanium implant Scalers Kit

NEW

**Our most loved instruments
for your best outcomes**

04

FOCUS ON

ENDODONTIC BLACK LINE

- 04 • Explorers and Excavators
- 05 • Pluggers, Spreaders and Rubber Dam Clamps

DIAGNOSTIC

- 06 • Colorvue Biotype Probes, Colorvue Probes and Expro
- 07 • Explorers, Dressing Pliers and Mirrors

PERIODONTAL

- 08 • EverEdge 2.0 and Sharpening Instruments

ORTHODONTIC

- 09 • Orthodontic Line

SURGICAL

- 10 • Needle Holders
- 11 • Scissors
- 12 • **Black Line** Surgical
- 13 • Scalpel Blades

RESTORATIVE

- 14 • Composite/Plastic Filling Instruments
- 15 • Carvers, Excavators and Cord Packers

**Instrument Management System
and Infection Prevention**

16

IMS PRODUCTS

- 16 • IMS Infinity Cassettes

Your daily recharge

18

IMS PRODUCTS

- 18 • Enzymax Detergent and Sporeclear Surfaces Disinfectant

Innovation by our Key Opinion Leaders

19

OUR KOL

- 19 • Swiss Perio Kit by Dr. Rino Burkhardt
- 21 • Opinion Leaders Education
Adhesive Dentistry Procedures - Werestore.it online Workshop



CUTTING-EDGE CLINICAL SOLUTIONS

CHECK OUT WHAT'S NEW
IN THE HU-FRIEDY FAMILY!

For more than one hundred years, Hu-Friedy has led the dental industry through **technological innovation**.

We consistently collaborate with **research centers, universities, and key opinion leaders** to understand the emerging needs of practitioners, to stay at the forefront of the latest dental technologies and methodologies, and to develop new and ever-better products

NEW 3SESSENTIAL KIT BY WERESTORE.IT

NEW


THE ESSENTIAL RESTORATIVE KIT INCLUDING "PLUG&SCULPT" INSTRUMENT

The Hu-Friedy **3SESSENTIAL KIT** by Werestore.it is the most straightforward kit of restorative ever made, with just **3 essential instruments at once: Anterior, Posterior and Spatula**, from the Hu-Friedy Black Line Collection.

FEATURES AND BENEFITS

- **Essential method:** having only three instruments makes restorative procedures straightforward
- **Plug&Sculpt concept:** The **Posterior** can plug and sculpt the composite by using its two different tips
- **Solid brush:** The **Anterior** models and spreads the composite easily, just like a brush
- **Different colours:** each colour helps to easily identify the instrument: RED for Anterior, BLUE for Posterior and GREY for Spatula
- **Black Line Collection instruments:** smooth surface reduces stickiness to composite material and the black colour has more contrast to the dental hard and soft tissues.
- **Ergonomics and visibility improved:** allow dental clinicians to perform at their best.

TNCIPAOX Anterior



TNEXBKR1X Posterior



CVKOT1X Spatula



3ssential Kit by Werestore.it
IMS Infinity Cassette 1/4 DIN 5 instruments Black
+ 3 instruments (Spatula, Posterior, Anterior)

RESTOREITDIN

Sugg. Retail Price

97.231 Ft

3SESSENTIAL instruments: CODE

XTS PFI Essential WeRestore.it Anterior **TNCIPAOX**

XTS PFI Essential WeRestore.it Posterior **TNEXBKR1X**

XTS PFI Essential WeRestore.it Spatula **CVKOT1X**

Sugg. Retail Price (per instrument)

19.907 Ft



WeRestore .it

Would you like to discover
some tips and tricks to
achieve a perfect isolation?

Do you know which are the
preparation tools and how to
use them?

Preheated composite, flowable
composite or
dual-cure cement for overlay
composite luting?

**CHECK THIS OUT AND MUCH
MORE ON THE WERESTORE.IT
CHANNEL, see page 21.**




NEW TITANIUM IMPLANT SCALERS

PERFORM AT YOUR BEST with Hu-Friedy's new, state of the art line of Titanium Implant Scalers. Made from the same titanium alloy as implants and abutments, these scalers are expertly designed and manufactured for implant maintenance, debridement, biofilm removal and can be used both supra and subgingival.

- **Improved visual acuity** with Hu-Friedy's signature Teal colored Anodized Titanium
- **Optimized sharpness** with cutting edges that are finely honed and sharpened post Anodization process
- **Enhanced instrument characteristics** for Implant maintenance with the ability to remove calculus and biofilm
- **Increased instrument value** with the ability to be sharpened at any time to extend the useful life of the instrument
- **Ergonomic sensitivity** featuring our Ergo handle that offers a large diameter and cross/hatch design for maximum tactile grip and support



IMP204SDT IMPLG1/2T IMPM1/2T IMPM11/12T IMPM13/14T



IME14DIN53

DESCRIPTION	CODE	Sugg. Retail Price
Complete Kit IMS Infinity Cassette 1/4 DIN Orange + 5 instruments	IMEDIMPLM	153.600 Ft
DESCRIPTION	CODE	Sugg. Retail Price (per instrument)
204SD	IMP204SDT	26.194 Ft
Langer 1/2	IMPLG1/2T	
Mini Five 1/2	IMPM1/2T	
Mini Five 11/12	IMPM11/12T	
Mini Five 13/14	IMPM13/14T	

OUR MOST LOVED INSTRUMENTS FOR YOUR BEST OUTCOMES

With more than **10,000 products** in a wide range of specialties that include orthodontics, pedodontics, restoratives and more, we lead the industry in innovation and technology development.

Our products are meticulously produced by our own in-house, **highly skilled artisans** and we provide **the highest-level of full-service** support to our customers before, during and especially after the sale.

ENDODONTIC



FOCUS ON BLACK LINE ENDODONTIC INSTRUMENTS

EXPLORERS

EXCAVATORS



#DG16 DE*,
BLACK LINE
EXDG16/23X

Sugg.
Retail Price
15.926 Ft

#DG16/17 DE*,
BLACK LINE
EXDG16/17X

Sugg.
Retail Price
15.926 Ft



#31W DE*,
Oval Spoon,
BLACK LINE
EXC31WX

Sugg.
Retail Price
20.955 Ft

#31L DE*,
BLACK LINE
EXC31LX

Sugg.
Retail Price
20.955 Ft



#33L DE*,
BLACK LINE
EXC33LX

Sugg.
Retail Price
20.955 Ft

BLACK LINE

Featuring a performance engineered coating to enhance lubricity, reduce light reflection under harsh lighting, and optimize edge retention, the **Black Line Endo instruments** are designed to provide exceptional clinical outcomes and improved efficiency during endodontic procedures.



FOCUS ON BLACK LINE ENDODONTIC INSTRUMENTS

PLUGGERS



#5/7 Root Canal DE**,
BLACK LINE RCP5/7X

Sugg. Retail Price
20.326 Ft

#1/3 Root Canal DE*,
BLACK LINE RCP1/3X

Sugg. Retail Price
20.326 Ft

SPREADERS



#D11TS SE**,
Root Canal,
BLACK LINE RCSD11TSX

Sugg. Retail Price
16.135 Ft

RUBBER DAM CLAMPS

Hu-Friedy's 48 clamps, made in the USA from proprietary stainless steel, provide full-mouth coverage.

POINTS OF PERFORMANCE

- **Satin Steel®** which resists corrosion while aiding in flexibility and memory
- **A dark coating** allowing for improved intraoral contrast and visual acuity
- **Matte finish** for decreased light reflection and glare



Sugg. Retail Price 7.963 Ft - Valid for the complete Black Line RDC range



*DE=Double Ended **SE=Single Ended



COLORVUE BIOTYPE PROBES BY DR. RASPERINI & DR. TESTORI

White probe tips

Green probe tips

Blue probe tips

Starter Kit with 3 handles and 12 tips
(4 green, 4 blue, 4 white)
PBTPKIT12

Sugg. Retail Price 76.486 Ft
Promo Price 53.540 Ft

COLORVUE PROBES

Markings
1-2-3-4-5-
6-7-8-9-10-
11-12 mm

UNC12
6 Pack,
1 handle + 7 tips
PCVNCKIT6

Sugg. Retail Price
33.738 Ft
**Promo Price
23.616 Ft**

EXPRO HANDLES

TU17 Expro
and Handle 6*
XPTU17/PH6

Sugg. Retail Price
10.268 Ft
**Promo Price
7.188 Ft**

*Colorvue PerioScreen tips
sold separately.



EXPLORERS



Nr. 23 Shepherd's Hook, Handle 30
EXS23

Sugg. Retail Price
6.915 Ft


Promo Price
4.841 Ft



Nr. 9, Handle 6 SE**
EXS96

Sugg. Retail Price
7.963 Ft

Promo Price
5.574 Ft



Nr. 5, Handle 6 DE*
EXD56

Sugg. Retail Price
10.687 Ft

Promo Price
7.481 Ft

DRESSING PLIERS



Nr. 2 College
DP2

Sugg. Retail Price
8.801 Ft

Promo Price
6.601 Ft



Nr. 17
DP17

Sugg. Retail Price
9.220 Ft

Promo Price
6.915 Ft



MIRRORS



ø 22 mm

Nr. 4, European Simple Stem, Single-Sided HD Mirror, Single
MIR4HDE

Sugg. Retail Price
5.239 Ft

Promo Price
4.715 Ft

Valid also for **MIR3HDE** and **MIR5HDE**

Nr. 6 European Simple Stem Mirror Handle
MHE6

Sugg. Retail Price
9.011 Ft

Promo Price
6.307 Ft





EVEREDGE® 2.0 CURETTES & SCALERS



EverEdge 2.0
Curettes &
Scalers,
Handle 9

Sugg.
Retail Price
22.212 Ft

Promo Price
16.659 Ft

Valid for the complete
EverEdge 2.0 Range



EverEdge 2.0
Curettes &
Scalers,
Handle Resin 8

Sugg.
Retail Price
20.955 Ft

Promo Price
15.716 Ft

Valid for the complete
EverEdge 2.0 Range

FOCUS ON: AUTOMATIC OR MANUAL? SHARPENING INSTRUMENTS

► AUTOMATIC

The Hu-Friedy Sidekick sharpener is designed for routine maintenance sharpening of scalers and curettes. Correct routine sharpening helps extend instrument life.



Sidekick Sharpening Kit Includes:
Sidekick Unit, User Manual & CD-ROM,
Ceramic Sharpening Stone, Plastic Test Stick,
Magnifying Lens and 2 AA Batteries.
SDKKIT

Sugg. Retail Price 191.529 Ft

Promo Price 162.799 Ft

SHARPENING STONES

► MANUAL



NO OIL NEEDED

Dual Grit Ceramic Stone
SSDGRIT

Sugg. Retail Price
14.878 Ft

Promo Price
10.415 Ft



OIL NEEDED

Nr. 4 Flat Arkansas
Hard Fine Grit
SS4

Sugg. Retail Price
17.183 Ft

Promo Price
14.606 Ft



Sharpen-Ez,
Sharpening Stone Oil
SSO

Sugg. Retail Price
8.172 Ft

Promo Price
6.947 Ft



ORTHODONTIC LINE

VALID FOR ALL PART CODES OF ORTHO STANDARD PLIERS. See our General Catalogue at page H10-H13.



UTILITY PLIERS

Strong, yet lightweight and perfectly balanced, the dimensions of these instruments are perfectly proportioned, they fit your grasp with unsurpassed ease.

Designed to provide superior control with effortless versatility, Hu-Friedy's Utility Pliers assure comfort to both you and your patients.

Slim Weingart Plier,
wire size \leq .030"
(.76 mm)
678-202

Sugg. Retail Price 61.189 Ft
Promo Price 48.951 Ft



WIRE FORMING PLIERS

These instruments perform as natural extensions of your professional artistry.

Our impeccably smooth finished tips allow you to bend the wire without risk of damage. Whether you are shaping stainless steel or nickel titanium, Hu-Friedy has a hand-finished, precision instrument to help bring your brilliant workmanship to life.

Tweed Loop Forming Plier,
wire size \leq .022" x .025"
(.56 x .64 mm)
678-306

Sugg. Retail Price 61.189 Ft
Promo Price 48.951 Ft

20% OFF FROM RETAIL PRICES



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.** Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS) are tested on Hu-Friedy instruments.

ILS
Instrument Lubricant Spray 236 ml

Sugg. Retail Price 10.897 Ft
Promo Price 7.628 Ft

IPS
Penetrating Oil 236 ml

Sugg. Retail Price 10.897 Ft
Promo Price 7.628 Ft






NEEDLE HOLDERS



Perma Sharp Hu-Friedy Perma Sharp needle holders feature tungsten carbide inserts with a matched diamond-cut pyramidal profile for enhanced longevity, a secure, positive suture needle grasp and reduced risk of cutting the suture material. With a precision designed hinge, the smooth surface reduces the risk of catching and or cutting the suture material during the knot tying process, enhancing clinical efficiencies. For optimal grip and control, Hu-Friedy offers a variety of design options, including standard, finger and palm controlled closure.



Castroviejo
Microsurgical,
14 cm
Recommended suture
sizes: 5-0, 6-0, 7-0, 8-0
NHCV

Sugg.
Retail Price
104.356 Ft

Promo Price
73.049 Ft



Castroviejo
Microsurgical,
18 cm
Recommended suture
sizes: 5-0, 6-0, 7-0, 8-0
NHCV18

Sugg.
Retail Price
113.576 Ft

Promo Price
79.503 Ft



Simion Straight
Castroviejo
Perma Sharp,
round handle,
18 cm
Recommended suture
sizes: 3-0, 4-0
NH5024SIM

Sugg.
Retail Price
167.850 Ft

Promo Price
117.495 Ft



Curved beak
Castroviejo
Perma Sharp,
round handle,
18 cm
Recommended suture
sizes: 5-0, 6-0, 7-0,
8-0, 9-0, 10-0
NH5024RC

Sugg.
Retail Price
167.850 Ft


Promo Price
117.495 Ft



Micro Straight
Castro
Diamond Dusted,
Perma Sharp,
round handle,
18 cm
Recommended suture
sizes: 6-0, 7-0, 8-0, 9-0
NHM-5026R

Sugg.
Retail Price
152.762 Ft

Promo Price
106.933 Ft



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilized in the open position.** Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS) are tested on Hu-Friedy instruments.

ILS
Instrument Lubricant Spray 236 ml
Sugg. Retail Price 10.897 Ft
Promo Price 7.628 Ft

IPS
Penetrating Oil 236 ml
Sugg. Retail Price 10.897 Ft
Promo Price 7.628 Ft



SCISSORS

SURGICAL

Super-Cut cutting blades are honed/serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

Perma Sharp

Scissors feature tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.



Nr. 14 LaGrange, Double-Curved, 1 serrated blade, 11.5 cm

S14

Sugg. Retail Price

25.775 Ft

Promo Price

20.620 Ft



Nr. 16 Goldman-Fox Curved, 1 serrated blade. Tapered beaks for interproximal use. 12.5 cm

S16

Sugg. Retail Price

19.069 Ft

Promo Price

15.255 Ft



Cortellini, Curved, 12 cm

SCORT

Sugg. Retail Price

23.260 Ft

Promo Price

18.608 Ft



Cortellini Curved, **Super-Cut** with 1 serrated blade. 12 cm


SCORTSC

Sugg. Retail Price

41.910 Ft

Promo Price

33.528 Ft



Straight Goldman-Fox **Perma Sharp**, 1 serrated blade, 13 cm

S5080

Sugg. Retail Price

75.857 Ft

Promo Price

60.686 Ft




BLACK LINE SURGICAL PERIOSTEALS



A full thickness mucoperiosteal pouch raised by blunt subperiosteal dissection.



DESCRIPTION	CODE	Sugg. Retail Price	Promo Price
Nr. 3 Prichard	PPR3X	42.120 Ft	33.696 Ft
Nr. 9 Molt	P9X	33.528 Ft	26.822 Ft
Allen, Anterior	PPAELAX	46.101 Ft	36.881 Ft
Allen, Anterior, short handle	PPAELASX	46.101 Ft	36.881 Ft
Buser	PPBUSERX	37.929 Ft	30.343 Ft
Buser, short handle	PPBUSERSX	37.929 Ft	30.343 Ft
Nr. 152 Kramer-Nevis	PKN152X	33.528 Ft	26.822 Ft
Nr. 152 Kramer-Nevis, short handle	PKN152SX	33.528 Ft	26.822 Ft



PERIODONTAL KNIVES



Incision made within the sulcus to sever the attachment from the base of the sulcus to the osseous crest.



DESCRIPTION	CODE	Sugg. Retail Price	Promo Price
Nr. 1 Quincey Tunneling	FOUINCEY1X	44.425 Ft	35.540 Ft
Nr. 2 Quincey Tunneling	FOUINCEY2X	44.425 Ft	35.540 Ft
Nr. 1 Tunneling	TKN1X	44.844 Ft	35.875 Ft
Nr. 1 Tunneling, short handle	TKN1SX	44.844 Ft	35.875 Ft
Nr. 2 Tunneling	TKN2X	44.844 Ft	35.875 Ft
Nr. 2 Tunneling, short handle	TKN2SX	44.844 Ft	35.875 Ft
Nr. 1/2 Allen Orban, Round	K012KP03RX	33.528 Ft	26.822 Ft
Allen End- Cutting Intrасulcular	KPAX	33.528 Ft	26.822 Ft
Nr. 15/16 Kirkland DE	KK15/16X	34.785 Ft	27.828 Ft
Nr. 15/16 Kirkland DE, short handle	KK15/16SX	34.785 Ft	27.828 Ft
1/2 Orban DE	K01/2X	34.785 Ft	27.828 Ft
11 Goldman-Fox DE	KGF11X	34.785 Ft	27.828 Ft

20% OFF FROM RETAIL PRICES

PERIODONTAL BONE CHISELS



Buser Chisel recontouring the alveolar crest prior to implant placement.



DESCRIPTION	CODE	Sugg. Retail Price	Promo Price
Nr. 1/2 Buser Modified 4 mm/5 mm	PPBUSER12X	38.767 Ft	31.013 Ft
Nr. 2/3 Buser Modified 5 mm/6 mm	PPBUSER23X	38.767 Ft	31.013 Ft
Nr. 1/3 Buser Modified 4 mm/6 mm	PPBUSER13X	38.767 Ft	31.013 Ft

WHY CHOOSE BLACK LINE SURGICAL INSTRUMENTS?

- Performance engineered coating for a harder, smoother surface for optimal edge retention and enhanced lubricity
- Reduced light reflection afforded by a matte finished handle and black working ends
- Unique, smooth, large diameter lightweight handle for maximum comfort, reduced hand fatigue and increased control



20% OFF FROM RETAIL PRICES

SINUS LIFT INSTRUMENTS



Kramer-Nevins Sinus Lift instrument delicately separating the Schneiderian Membrane from the sinus floor.



DESCRIPTION	CODE	Sugg. Retail Price	Promo Price
Kramer-Nevins, Obtuse	IMP6578SX	45.263 Ft	36.210 Ft
Kramer-Nevins, Acute	IMP6577SCX	45.263 Ft	36.210 Ft
Nr. 1 Sinus Curette	SINC1X	44.425 Ft	35.540 Ft
Nr. 2 Sinus Curette	SINC2X	44.425 Ft	35.540 Ft
Nr. 2 Sinus Curette, short handle	SINC2SX	44.425 Ft	35.540 Ft
Nr. 1 Sinus Curette, short handle	SINC1SX	44.425 Ft	35.540 Ft

SURGICAL CURETTE



Molt curette debriding the osteotomy.



DESCRIPTION	CODE	Sugg. Retail Price	Promo Price
Nr. 2/4 Molt	CM2/4X	36.043 Ft	28.834 Ft
Nr. 86 Lucas	CL86X	35.833 Ft	28.666 Ft
Nr. 87 Lucas	CL87X	35.833 Ft	28.666 Ft

SCALPELS AND BLADES



Nr. 12 Sterile Blade, Box of 100, Carbon Steel, Expiry Date **10-255-12**

Sugg. Retail Price
25.356 Ft

Promo Price
17.749 Ft



Nr. 12D* Sterile Blade, Box of 100, Carbon Steel, Expiry Date **10-256-12**

Sugg. Retail Price
30.594 Ft

Promo Price
21.416 Ft



Nr. 15 Sterile Blade, Box of 100, Carbon Steel, Expiry Date **10-255-15**

Sugg. Retail Price
25.356 Ft

Promo Price
17.749 Ft



Nr. 15C Sterile Blade, Box of 100, Carbon Steel, Expiry Date **10-256-15**

Sugg. Retail Price
30.594 Ft

Promo Price
21.416 Ft



Nr. 5E, European Style Round Scalpel Handle, Matte Finish
10-130-05E

Sugg. Retail Price
30.385 Ft

Promo Price
21.269 Ft





COMPOSITE/ PLASTIC FILLING INSTRUMENTS



Nr. 1, Goldstein Flexi-Thin, Handle 6
CIGFT16

Sugg. Retail Price
17.183 Ft


Promo Price
11.169 Ft



Nr. 2, Goldstein Flexi-Thin, Handle 6
CIGFT26

Sugg. Retail Price
17.183 Ft

Promo Price
11.169 Ft



Nr. 4, Goldstein Flexi-Thin, Handle 6
CIGFT46

Sugg. Retail Price
17.183 Ft

Promo Price
11.169 Ft



Nr. 3, Goldstein Flexi-Thin, Handle 6
CIGFT36

Sugg. Retail Price
17.183 Ft

Promo Price
11.169 Ft



Nr. 5, Goldstein Flexi-Thin, Handle 41
CIGFT5

Sugg. Retail Price
15.088 Ft


Promo Price
9.807 Ft



Nr. 2, Goldstein Flexi-Thin, Handle 41
CIGFT6

Sugg. Retail Price
15.088 Ft

Promo Price
9.807 Ft





CARVERS



Nr. 8 Wiland,
Handle 41
CVW18

Sugg.
Retail Price
14.669 Ft

Promo Price
9.535 Ft



IPC Interproximal,
Handle 41
CVIPC

Sugg.
Retail Price
14.669 Ft

Promo Price
9.535 Ft

EXCAVATORS



Nr. 17WH Heavy
Oval Spoon,
Handle 41
EXC17WH

Sugg.
Retail Price
15.088 Ft

Promo Price
9.807 Ft



Nr. 18WH Heavy
Oval Spoon,
Handle 41
EXC18WH

Sugg.
Retail Price
15.088 Ft

Promo Price
9.807 Ft


CORD PACKERS




Nr. S6 Gingival,
Handle 41
GCPS6

Sugg.
Retail Price
17.183 Ft

Promo Price
11.169 Ft





INSTRUMENT MANAGEMENT SYSTEM AND INFECTION PREVENTION



Hu-Friedy's **Instrument Management System** is the most comprehensive, reliable, and complete instrument processing and infection control system on the market today. IMS combines **sterilisation, storage and organization** of instruments in one **integrated system**.

Hu-Friedy's IMS system allows you to focus on your patients rather than wasting valuable time looking for the instruments to complete a procedure.

WORK SMARTER, NOT HARDER: CHOOSE THE IMS SOLUTION

Why using IMS Solution?

- Enhanced chairside efficiency and easy training with standardised procedural set-ups incrementing the productivity
- Reduced set-up time with organised storage of instruments within procedural labelled cassettes
- Entire cassette can be placed into Thermo-disinfector, ultrasonic bath and autoclave
- Reduced need for manual sorting of instruments into pouches
- Cleaning time reduced, allowing for a more productive team.














IME14DIN5_

IMS Infinity Cassette 1/4 DIN,
5 Instruments
Dimensions: 72 x 184 x 34 mm.


Choose the color.

Just add the color coding number to the last part of the partcode (replace "_" with the number)


 0 = Red	 6 = Lavander
 1 = Grey	 7 = Black
 2 = White	 8 = Blue
 3 = Orange	 9 = Green
 4 = Purple	 11 = Ocean Blue
 5 = Yellow	

Sugg. Retail Price 49.035 Ft

Promo Price 36.776 Ft



Visit our new IMS INSTRUMENT MANAGEMENT SYSTEM site
imsuser.hu-friedy.eu/en/
ORGANISATION, EFFICIENCY and SAFETY in a simple click!




Sugg. Retail Price 55.112 Ft

Promo Price 41.334 Ft

IME12DIN1_

IMS Infinity Cassette 1/2 DIN, 8 Instruments with an accessory area. Ideal for diagnostics and prophylaxis.
Dimensioni: 152 x 184 x 34 mm



IMEDIN18_

IMS Infinity Cassette DIN, 18 Instruments
Dimensioni: 290 x 184 x 34 mm

Sugg. Retail Price 74.181 Ft

Promo Price 55.636 Ft



IMEDIN18_S

IMS Infinity Cassette DIN, 18 Instruments, Surgical, with acc. area and clips
Dimensioni: 290 x 184 x 34 mm

Sugg. Retail Price 79.629 Ft

Promo Price 59.722 Ft

Hu-Friedy



YOUR DAILY RECHARGE

We strive to support your day by day profession by providing a wide range of **dental consumables** from hygiene to surgical products.

Our **commitment** remains the same: supporting and enabling the **most advanced clinical procedures** and the **most efficient way to work.**

IMS ENZYMAX DETERGENT



Enzymax Spray Gel: Dual enzyme gel with corrosion inhibitor acts as a **ready-to-use instrument pre-cleaner**. Maintains moisture to help prevent debris from drying on instruments before cleaning and **helps dissolve blood, tissue, and other debris**. Can also be used as a linen spot remover.

Sugg. Retail Price

9.639 Ft

Promo Price

6.266 Ft

IMS Enzymax Spray Gel, 709 ml, up to 450 Sprays
IMS-1229

Enzymax Dual Enzyme Ultrasonic Detergent is specifically formulated for efficient ultrasonic cleaning or pre-soaking of cassettes and instruments. Contains Steelgard agents to protect instruments from harmful minerals naturally found in water that contribute to buildup, spotting and corrosion. Can also be used as an evacuation system cleaner.



IMS Enzymax Liquid, Quart Bottle, Makes 121 lts.
IMS-1224

Sugg. Retail Price

20.955 Ft

Promo Price

13.621 Ft



IMS Enzymax Liquid, Gallon, Makes 485 lts.
IMS-1226

Sugg. Retail Price

42.329 Ft

Promo Price

27.514 Ft

IMS SPORE CLEAR SURFACE DISINFECTANTS

This one-step formula **cleans and disinfects** critical workspaces **within 1 minute of contact**. It offers a fast, simple and effective answer to creating a safe work environment protecting dental staff and their patients from harmful pathogens. **Sporeclear** is **non-alcohol based, free of aldehydes, chlorine, phenols, peroxide** and **nonylphenol ethoxylates**, (NPEs), making them **compatible with all washable surfaces**.

Fast solutions for effective disinfection in LESS THAN 1 MINUTE
SPORICIDAL - TUBERCULOCIDAL - BACTERICIDAL - VIRUCIDAL - FUNGICIDAL



Sporeclear® Disinfectant for non-invasive Medical Devices, 160 wipes canister
MSCD354160

Sugg. Retail Price

4.610 Ft

Promo Price

3.688 Ft



Sporeclear® Disinfectant for non-invasive Medical Devices, 160 wipes refill pack
MSCD354161

Sugg. Retail Price

3.981 Ft

Promo Price

3.185 Ft



Sporeclear® Disinfectant for non-invasive Medical Devices, Liquid Spray, 750 ml bottle (ready-to-use)
MSCD356750

Sugg. Retail Price

5.448 Ft

Promo Price

4.359 Ft



INNOVATION BY OUR KEY OPINION LEADERS

The collaboration of the most **influential dental professionals** along with our world class team of engineers working together to develop the products of the future is the secret of our **continuous innovation** that, combined with the talents of our team, allows us to **anticipate the needs of our customers** in every country and build a customised solution, unique to their individual requirements.

OUR KEY OPINION LEADERS SWISS PERIO KIT - DR. RINO BURKHARDT


THE SWISS PERIO KIT INCLUDES:

IMS Infinity DIN Surgical cassette for 18 instruments, black	IMEDIN187S
Swiss Perio Microsurgical Needle Holder Stainless Steel	SPNHDPVN
Swiss Perio Microsurgical Plier Straight Stainless Steel	SPTPDSSPV
Swiss Perio Microsurgical Plier Straight Stainless Steel, anatomical	SPTPDAPV
Swiss Perio Microsurgical Scissor Stainless Steel	SPSPV
Swiss Perio Microsurgical Scissors Stainless Steel, 45 degree	SPSPV45
Swiss Perio Periosteal Elevators Mod. P24G, Black Line	SPP24GSPX
Swiss Perio Periosteal Elevators Mod. P24GSIC, Black Line	SPP24GSICX
Swiss Perio Tissue Dissector Mod. PH2, Black Line	SPPH2X
Swiss Perio Papilla Elevator Mod. PPSCHLEE, Black Line	SPPSCHLEEX
84 Lucas Double Ended Surgical Curette, Black Line	CL84X
Swiss Perio Periosteal Elevators Mod. P24GSP, Black Line	SPP24GSP2X
YG7/8 Younger-Good, EverEdge 2.0, Smooth Handle Stainless Steel	SYG7/8CHE2
Colorvue Probe Handle, Smooth Handle Stainless Steel	PH6CH
Mini Scalpel Blade / Mirror Handle Stainless Steel	SHDPV
5 European Style Round Scalp Handle	10-130-05E
N.6 Simple Steam Mirror Handle, Smooth Handle	MHE6CH
17 Dressing Pliers	DP17
Q2N Color-Coded Nabers Probe, Smooth Handle	PQ2NCH
12 Colorvue Probe tips, 12 Replacements Tips	PCV12PT
European Simple Steamn.5 Double-Sided HD Mirror, Single	MIR5DSHDE

Complete Swiss Perio Kit
SPKITDIN

Sugg. Retail Price
1.369.828 Ft

Promo Price
1.095.863 Ft



Rino Burkhardt, one of the Board Members of SwissPerio Education, is responsible for the training programs. He is an EFP (European Federation of Periodontology) certified specialist in periodontology and received his Masters degree from the Medical Faculty of the University of Berne, Switzerland. Since 1996 he runs a private clinic in Zurich, limited to periodontology and implantology. Additionally, he holds an Honorary Professorship at the University of Hong Kong and a lectureship at the University of Zurich. He is an active member of the Swiss Society of Periodontology (SSP), the European Academy of Aesthetic Dentistry (EAED), the European Association for Osseointegration (EAO), and Board Member of the Swiss Society of Oral Implantology (SGI). He published numerous book chapters and articles and won a research prize from the EFP.



OUR KEY OPINION LEADERS SWISS PERIO KIT - DR. RINO BURKHARDT

BEST SOLD SWISS PERIO INSTRUMENTS



Swiss Perio
Microsurgical
Needle Holder
SPNHDPVN

Sugg.
Retail Price
277.654 Ft

Promo Price
194.358 Ft



Swiss Perio
Microsurgical
Plier Straight
SPTPDSSPV

Sugg.
Retail Price
189.433 Ft

Promo Price
132.603 Ft



Swiss Perio
Microsurgical
Plier Straight,
anatomical
SPTDAPV

Sugg.
Retail Price
148.571 Ft

Promo Price
104.000 Ft



Swiss Perio
Microsurgical
Scissor
SPSPV

Sugg.
Retail Price
221.285 Ft

Promo Price
154.899 Ft



Swiss Perio
Microsurgical
Scissor, 45°
SPSPV45

Sugg.
Retail Price
265.500 Ft

Promo Price
185.850 Ft



WeRestore .it

WeRestore.it is a project founded in 2016 by Gaetano Paolone and Salvatore Scolavino.

WeRestore.it aims for excellence in adhesive direct and indirect restorative dentistry. Both Salvatore and Gaetano are active members of the Italian Academy of Restorative Dentistry. They wrote a book about direct posterior restorations edited by Quintessence Int. Publishing.

Gaetano is adjunct professor of Restorative Dentistry at Vita e Salute University, San Raffaele, Milan. They lecture all around the world in restorative dentistry and run master courses in the WeRestore.it course center in Rome.


OPINION LEADERS EDUCATION **GAETANO PAOLONE & SALVATORE SCOLAVINO**

Would you like to discover some tips and tricks to achieve a perfect isolation?

Do you know which are the preparation tools and how to use them?

Preheated composite, flowable composite or dual-cure cement for overlay composite luting?

CHECK THIS OUT AND MUCH MORE ON WERESTORE.IT CHANNEL.



FOR MORE INFORMATION:
www.adhesive-dentistry.com

- Online video course
- Full HD Quality
- PC, MAC, Tablet e Mobile

• Hosted on:  **osteocom**
THE DENTAL LEARNING PLATFORM

The Adhesive Dentistry Procedures online channel aims to teach on a wide variety of innovative procedures concerning Adhesive Dentistry, in order to be aligned on techniques, materials and clinical advances in both direct and indirect restorations.

Hu-Friedy

CONTACTS

Riccardo Boschetto

Junior Territory Manager Eastern Europe,

Via Mauro Macchi, 27 - 20124 Milan, Italy

Tel. +39 02 36 58 96 00 ext. 202

Mob. +39 340 885 96 43

Mail rboschetto@hu-friedy.com

OUR NEXT GIGS...



IDS 2019 - HALL 10.1 - Booth D040/E042

12th - 16th March 2019, Cologne



DOWNLOAD OUR CATALOGUE

For all of us at Hu-Friedy, this is much more than a book that contains the highest quality dental instruments and products. It's a manifestation of the dedication, innovative spirit and hard work of a company that has served and supported the dental community for over 110 years.

For more information about
Hu-Friedy's Total Solutions visit

HU-FRIEDY.EU


PREVENTIVE




RESTORATIVE



ORTHODONTIC




INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



POWER



Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA
E-Mail: info@hu-friedy.eu | Website: Hu-Friedy.eu

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved. The contents of this Product Catalog & Reference Guide including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of Hu-Friedy Mfg. Co., LLC and may not be copied, reproduced or distributed without prior written permission. Hu-Friedy Mfg. Co., LLC. reserves the right to change features of products, images and text at any time without notice. HFL-070/0219

How the best perform

Hu-Friedy