

# HU-FRIEDY NATIONAL PROMOTIONS

---

## HUNGARY

Valid from 2<sup>nd</sup> March until 31<sup>st</sup> May 2020

All prices are VAT included

For more information about Hu-Friedy's Total Solutions visit [HU-FRIEDY.EU](https://hu-friedy.eu)

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved.

How the best perform





# NEW AKRO-FLEX™

| Create Beautiful Restorations

Akro-Flex™ is an extremely flexible composite instrument which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations.

**The inclusion of a ductile material allows Akro-Flex to perform as a solid brush.**



**Akro-Flex**  
by Werestore.it

Watch our tutorial video for Akro-Flex  
on YouTube

Video courtesy of werestore.it



**Akro-Flex**  
by Dr. Alessandro Conti

Watch our tutorial video for Akro-Flex  
on YouTube

Video courtesy of Dr. Alessandro Conti



## UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex acts as a solid brush. The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.



## HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease. The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



## AKRO-FLEX

Akro-Flex Composite Instrument

Sugg. Retail Price

**31.090 FT**

## ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.



Patent pending

Patent pending





## CONSUMABLES

### SURFACE DISINFECTANT



One-step formula, cleans and disinfects critical workspaces within 1 minute of contact, offers a fast, simple and effective answer to creating a safe work environment. Sporeclear is non-alcohol based, free of aldehydes, chlorine, phenols, peroxide and nonylphenol ethoxylates, (NPEs), making them compatible with all washable surfaces and eco-friendly.

#### MSCD354200

IMS Sporeclear Disinfectant for non-invasive Medical Devices, Wipes, 200 wipes canister

Sugg. Retail Price 5.959 FT  
**Promo 3.873 FT**

### IMS WRAPS



Disposable paper for sterile packing of the cassettes. After the packed cassette has been opened, the sterilisation paper may be used as a base, thus forming a predominantly sterile environment. Additionally, contaminated cotton wool rolls, gauze and similar items can be disposed of hygienically in this paper.

CODE	DESCRIPTION	Sugg. Retail Price	Promo
IMS-1210M	IMS Sterilisation Paper DIN and Container Cassettes 60x60cm 500pcs./pack	129.152 FT	<b>77.491 FT</b>
IMS-1210MH	IMS Sterilisation Paper DIN and Container Cassettes 60x60cm 100pcs./pack.	26.298 FT	<b>15.779 FT</b>
IMS-1211M	IMS Sterilisation Paper for 1/4 DIN Cassette 37,5x37,5cm 1000pcs./pack.	113.373 FT	<b>68.024 FT</b>
IMS-1211MH	IMS Sterilisation Paper for 1/4 DIN Cassette 37,5x37,5cm 100pcs./pack.	11.688 FT	<b>7.013 FT</b>
IMS-1212M	IMS Sterilisation Paper for 1/2 DIN Cassette 50x50cm 500pcs./pack.	101.393 FT	<b>60.836 FT</b>
IMS-1212MH	IMS Sterilisation Paper for 1/2 DIN Cassette 50x50cm 100pcs./pack.	20.746 FT	<b>12.448 FT</b>

### CLEANING & CARE



Keeps hinged instruments moving freely. Do not contains silicone.

#### ILS

IMS Lubricant Spray for instruments 8oz/237 ml Pump

Sugg. Retail Price 11.486 FT  
**Promo 6.892 FT**



Dual enzyme ultrasonic detergent/ presoak. A fast-acting, concentrated, liquid formula.

#### IMS-1226

IMS Detergent Enzymax liquid, 3,8 l Gallon Bottle

Sugg. Retail Price 45.037 FT  
**Promo 22.518 FT**



Dual enzyme gel with corrosion inhibitor acts as a ready-to-use instrument pre-cleaner. Maintains moisture to help prevent debris from drying on instruments before cleaning and helps dissolve blood, tissue, and other debris. Can also be used as a linen spot remover.


#### IMS-1229

IMS Enzymax Spray Gel, 710 ml bottle

Sugg. Retail Price 10.320 FT  
**Promo 6.192 FT**

## HD MIRRORS

#### HD VS RHODIUM COATING



HD Mirror provides superior brilliance and color for quicker and more accurate visibility of the mouth.\*

- **113 %\*** reflection factor for exceptional image clarity
- **38,5 %\*** brighter than rhodium coated mirror glass
- **50 %\*** brighter than other front surface mirror glass

\* Data on file and available upon request.



#4 **single-sided**, 22 mm Ø

VERSION	CODE	Sugg. Retail Price	Promo
European Stem, 3pcs./pack	MIR4HD3E	16.020 FT	<b>10.413 FT</b>
European Stem, 12pcs./pack	MIR4HD12E	37.739 FT	<b>26.418 FT</b>

#5 **single-sided**, 24 mm Ø

VERSION	CODE	Sugg. Retail Price	Promo
European Stem, 3pcs./pack	MIR5HD3E	16.020 FT	<b>10.413 FT</b>
European Stem, 12pcs./pack	MIR5HD12E	37.739 FT	<b>26.418 FT</b>



#4 **double-sided**, 22 mm Ø

VERSION	CODE	Sugg. Retail Price	Promo
European Stem, 1pc./pack	MIR4DSHDE	6.995 FT	<b>5.596 FT</b>

#5 **double-sided**, 24 mm Ø

VERSION	CODE	Sugg. Retail Price	Promo
European Stem, 1pc./pack	MIR5DSHDE	6.995 FT	<b>5.596 FT</b>



#### MHE6


#6 European Simple Stem Mirror Handle

Sugg. Retail Price 9.672 FT  
**Promo 5.610 FT**

## OXFORD PROBES



Colorvue probes feature a vivid yellow tip with black markings for enhanced contrast to intraoral structures in both natural and implant dentition. Optimal flexibility and a rounded tip provide greater patient comfort and acceptability. The probe tips are safe for use with implants and feature a twist-on design for easy replacement of worn tips.



#### PCVOKIT6

Colorvue Kit #UNC15 Oxford, 6 tips + 1 handle PH6, double-ended probe

Sugg. Retail Price 36.573 FT  
**Promo 21.213 FT**

#### PCVO12PT

Colorvue Tips #UNC15 Oxford, 12 tips/pack.

Sugg. Retail Price 48.966 FT  
**Promo 28.400 FT**

## PROBES AND EXPLORERS

1-2-3-4-5-6-7-8-9-10-11-12-13-14-15



#### PCPUNC156

Probe UNC #UNC15 handle #6

Sugg. Retail Price 14.854 FT



#### EXS66

Explorer #6 handle #6

Sugg. Retail Price 8.506 FT



#### EXS236

Explorer Shepherds Hook #23 handle #6

Sugg. Retail Price 8.506 FT



#### EXS96

Explorer #9 handle #6

Sugg. Retail Price 8.506 FT

**BUNDLE OFFER**  
**BUY 1 PCPUNC156 AND GET 1 EXPLORER OF YOUR CHOICE INCLUDED IN THE PRICE**



## EVEREDGE® 2.0 THE MOST ADVANCED SCALER EVER

**BUNDLE OFFER**  
**BUY 6 EVEREDGE 2.0 CURETTES**  
**+ SIDEKICK SHARPENER FOR ONLY**  
**190.039 FT (INSTEAD OF 345.526 FT)**

### SIDEKICK SHARPENER

The Hu-Friedy Sidekick sharpener is designed for routine maintenance sharpening of scalers and curettes. Correct routine sharpening helps extend instrument life.

- The Sidekick sharpener is cordless, for use anywhere in the dental office.
- Instrument Guide Channels and Vertical Backstops are “template-like” guides that help control blade angulation, providing consistent results.
- The Sidekick sharpener is compact and requires minimal space.

### SDKKIT

Sidekick Sharpening Kit

Kit Includes:  
Sidekick Unit, User Manual  
& CD-ROM, Ceramic  
Sharpening Stone, Plastic  
Test Stick, Magnifying Lens  
and 2 AA Batteries.

Sugg. Retail Price 204.069 FT



Watch our tutorial  
video for Sidekick  
on YouTube

### STANDARD GRACEY CURETTES

Description	Handle 9
1/2 Gracey Standard	SG1/29IE2
3/4 Gracey Standard	SG3/499E2
5/6 Gracey Standard	SG5/695E2
7/8 Gracey Standard	SG7/897E2
9/10 Gracey Standard	SG9/1096E2
11/12 Gracey Standard	SG11/1293E2
13/14 Gracey Standard	SG13/1498E2
15/16 Gracey Standard	SG15/1693E2
17/18 Gracey Standard	SG17/1898E2
11/14 Gracey Standard	SG11/149E2
12/13 Gracey Standard	SG12/139E2
1/2 Gracey Standard, Rigid	SG1/2R9E2
3/4 Gracey Standard, Rigid	SG3/4R9E2
5/6 Gracey Standard, Rigid	SG5/6R9E2
7/8 Gracey Standard, Rigid	SG7/8R9E2
9/10 Gracey Standard, Rigid	SG9/10R9E2
11/12 Gracey Standard, Rigid	SG11/12R9E2
13/14 Gracey Standard, Rigid	SG13/14R9E2
15/16 Gracey Standard, Rigid	SG15/16R9E2
17/18 Gracey Standard, Rigid	SG17/18R9E2

### AFTER FIVE GRACEY CURETTES

Description	Handle 9
1/2 After Five Gracey	SRPG1/29IE2
3/4 After Five Gracey	SRPG3/499E2
5/6 After Five Gracey	SRPG5/695E2
7/8 After Five Gracey	SRPG7/897E2
11/12 After Five Gracey	SRP11/1293E2
13/14 After Five Gracey	SRP13/1498E2
15/16 After Five Gracey	SRP15/1693E2
11/14 After Five Gracey	SRPG11/149E2
12/13 After Five Gracey	SRPG12/139E2
1/2 After Five Gracey, Rigid	SRP1/2R9E2
3/4 After Five Gracey, Rigid	SRP3/4R9E2
5/6 After Five Gracey, Rigid	SRP5/6R9E2
7/8 After Five Gracey, Rigid	SRP7/8R9E2
11/12 After Five Gracey, Rigid	SRP11/12R9E2
13/14 After Five Gracey, Rigid	SRP13/14R9E2
15/16 After Five Gracey, Rigid	SRP15/16R9E2

### MINI FIVE GRACEY CURETTES

Description	Handle 9
1/2 Mini Five Gracey	SAS1/29IE2
3/4 Mini Five Gracey	SAS3/49E2
5/6 Mini Five Gracey	SAS5/695E2
7/8 Mini Five Gracey	SAS7/897E2

### MINI FIVE GRACEY CURETTES

Description	Handle 9
11/12 Mini Five Gracey	SAS11/1293E2
13/14 Mini Five Gracey	SAS13/1498E2
15/16 Mini Five Gracey	SAS15/1693E2
17/18 Mini Five Gracey	SAS17/1898E2
11/14 Mini Five Gracey	SAS11/149E2
12/13 Mini Five Gracey	SAS12/139E2
1/2 Mini Five Gracey, Rigid	SAS1/2R9E2
3/4 Mini Five Gracey, Rigid	SAS3/4R9E2
5/6 Mini Five Gracey, Rigid	SAS5/6R9E2
7/8 Mini Five Gracey, Rigid	SAS7/8R9E2
11/12 Mini Five Gracey, Rigid	SAS11/12R9E2
13/14 Mini Five Gracey, Rigid	SAS13/14R9E2
15/16 Mini Five Gracey, Rigid	SAS15/16R9E2

### MICRO MINI FIVE GRACEY CURETTES

Description	Handle 9
1/2 Micro Mini Five Gracey	SM1/29IE2
7/8 Micro Mini Five Gracey	SMS7/897E2
11/12 Micro Mini Five Gracey	SMS11/1293E2
13/14 Micro Mini Five Gracey	SMS13/1498E2

### BIOGENT CURETTES

Description	Handle 9
1/2 Biogent	SG1/29LCE2
7/8 Biogent	SG7/89LCE2
11/12 Biogent	SG11/129LCE2
13/14 Biogent	SG13/149LCE2

### UNIVERSAL CURETTES

Description	Handle 9
1/2 Barnhart	SBH1/29E2
5/6 Barnhart	SBH5/69E2
1/2 Barnhart, USC/Suter	SBH1/2ZUSC9E2
5/6 Barnhart, USC/Suter	SBH5/6ZUSC9E2
1/2 Barnhart, Rigid	SBH1/2R9E2
5/6 Barnhart, Rigid	SBH5/6R9E2
2R/2L Columbia	SC2R/2L9E2
4R/4L Columbia	SC4R/4L9E2
13/14 Columbia	SC13/149E2
13/14 McCall	SM13/149E2
YG7/8 Younger-Good	SYG7/89E2
1/2 Langer	SL1/29E2
3/4 Langer	SL3/49E2
5/6 Langer	SL5/69E2
17/18 Langer	SL17/189E2

### UNIVERSAL CURETTES

Description	Handle 9
1/2 After Five Langer	SL1/2AF9E2
3/4 After Five Langer	SL3/4AF9E2
5/6 After Five Langer	SL5/6AF9E2
1/2 Mini Five Langer	SL1/2MF9E2
3/4 Mini Five Langer	SL3/4MF9E2
5/6 Mini Five Langer	SL5/6MF9E2

### SICKLE SCALERS

Description	Handle 9
204S Scaler	S204S9E2
204SD Scaler	S204SD9E2
H5/33 Scaler	SH5/339E2
H6/H7 Scaler	SH6/79E2
204 IUFW Scaler	SIUFW2049E2
30/33 Jacquette	SJ30/339E2
31/32 Jacquette	SJ31/329E2
34/35 Jacquette	SJ34/359E2
0/00 Morse	SM0/009E2
13S/14S McCall	SM13/14S9E2
17/18 McCall	SM17/189E2
17/18 McCall's, Ptd Tip	SM17/18PT9E2
23 Scaler	SM239E2
23A Scaler	SM23A9E2
128 Nebraska	SM128DE9E2
135 Scaler	SM1359E2
137 Scaler	SM1379E2
U15/30 Scaler	SU15/309E2
U15/33 Scaler	SU15/339E2
1 Nevi Anterior	SCNEV19E2
2 Nevi Posterior	SCNEV29E2
3 Nevi Scaler Posterior	SCNEV39E2
4 Nevi Scaler Posterior	SCNEV49E2
1 Nevi H5, Anterior	SNEV1/H59E2

Sugg. Retail Price 23.576 FT



## ULTRASONIC TIPS SERIE S

**BUY 1 AND  
GET 1 FOR FREE**

Compatible with  
Satelec/Acteon\* generators.





\* Our offered items are compatible with but not original products of device manufacturing companies Satelec and/or EMS and are not produced, licensed or recommended from them. All used brand names and definitions are registered trademarks and brands of Acteon, Inc. and Electro Medical Systems respectively. They are intended only to show compatibility of items with products of different manufacturers.

## ULTRASONIC TIPS SERIE E

**BUY 1 AND  
GET 1 FOR FREE**

Compatible with all  
EMS\* generators





## NEEDLE HOLDERS



**Perma Sharp** Hu-Friedy Perma Sharp needle holders feature tungsten carbide inserts with a matched diamond-cut pyramidal profile for enhanced longevity, a secure, positive suture needle grasp and reduced risk of cutting the suture material. With a precision designed hinge, the smooth surface reduces the risk of catching and or cutting the suture material during the knot tying process, enhancing clinical efficiencies. For optimal grip and control, Hu-Friedy offers a variety of design options, including standard, finger and palm controlled closure.



Castroviejo  
Microsurgical,  
14 cm  
**Recommended suture  
sizes: 5-0, 6-0, 7-0, 8-0**  
**NHCV**

Sugg. Retail Price  
111.102 FT

**Promo**  
**77.771 FT**



Castroviejo  
Microsurgical,  
18 cm  
**Recommended suture  
sizes: 5-0, 6-0, 7-0, 8-0**  
**NHCV18**

Sugg. Retail Price  
121.077 FT

**Promo**  
**84.754 FT**



Simion Straight  
Castroviejo  
Perma Sharp,  
round handle,  
18 cm  
**Recommended suture  
sizes: 3-0, 4-0**  
**NH5024SIM**

Sugg. Retail Price  
178.981 FT

**Promo**  
**125.287 FT**



Curved beak  
Castroviejo  
Perma Sharp,  
round handle,  
18 cm  
**Recommended suture  
sizes: 5-0, 6-0, 7-0,  
8-0, 9-0, 10-0**  
**NH5024RC**

Sugg. Retail Price  
178.981 FT

**Promo**  
**125.287 FT**



Micro Straight  
Castro  
Diamond Dusted,  
Perma Sharp,  
round handle,  
18 cm  
**Recommended suture  
sizes: 6-0, 7-0, 8-0, 9-0**  
**NHM-5026R**

Sugg. Retail Price  
162.961 FT

**Promo**  
**114.073 FT**




### DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Instrument Lubricant Spray (**ILS**), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.** Hu-Friedy Instrument Lubricant Spray (**ILS**) is tested on Hu-Friedy instruments.

**ILS**  
Instrument Lubricant Spray 236 ml

Sugg. Retail Price 11.486 FT  
**Promo 6.892 FT**

## SCISSORS



#14 LaGrange  
Double-Curved Super-Cut,  
1 serrated blade, 11.5 cm  
**S14SC**

Sugg. Retail Price  
45.900 FT

**Promo**  
**28.917 FT**



#16 Goldman-Fox  
Curved, Super-Cut,  
1 serrated blade, 12.5 cm  
**S16SC**

Sugg. Retail Price  
45.900 FT

**Promo**  
**26.622 FT**



#17 Straight Iris, Super-Cut,  
Delicate, 11.5 cm  
**S17SC**

Sugg. Retail Price  
35.926 FT

**Promo**  
**23.352 FT**



#18 Curved Iris, Super-Cut,  
Delicate, 1 serrated blade,  
11.5 cm  
**S18SC**

Sugg. Retail Price  
35.926 FT

**Promo**  
**23.352 FT**




Cortellini Curved, Super-Cut  
with 1 serrated blade, 12 cm  
**SCORTSC**

Sugg. Retail Price  
44.734 FT

**Promo**  
**31.314 FT**


**Super-Cut** cutting blades are honed/serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.







## SURGICAL CURETTES

**CL84**

#84 Lucas Spoon Shape, DE, handle #522

Sugg. Retail Price 34.760 FT  
**Promo 19.118 FT**

**CL85**

#85 Lucas Spoon Shape, DE, handle #522

Sugg. Retail Price 34.760 FT  
**Promo 19.118 FT**

**CL87**

#87 Lucas Spoon Shape, DE, handle #522

Sugg. Retail Price 34.760 FT  
**Promo 19.118 FT**

**CL86**

#86 Lucas Spoon Shape, DE, handle #522

Sugg. Retail Price 34.760 FT  
**Promo 19.118 FT**

**CM2/4**

#2/4 Molt Surgical Curette, DE, handle #522

Sugg. Retail Price 35.062 FT  
**Promo 19.284 FT**

BLADES AND  
BLADE HANDLES**10-130-05E**

#5E European Round Scalpel Handle, Matte Finish

Sugg. Retail Price 32.385 FT  
**Promo 21.050 FT**

**10-130-03E**

#3E European Scalpel Handle #3 graduated, Matte Finish

Sugg. Retail Price 6.391 FT  
**Promo 5.304 FT**

**10-255-12**

#12 Sterile Blade, Box of 100, Carbon Steel

Sugg. Retail Price 17.704 FT  
**Promo 10.976 FT**

**10-255-15**

#15 Sterile Blade, Box of 100, Carbon Steel

Sugg. Retail Price 17.704 FT  
**Promo 10.976 FT**

**10-256-12**

#12D Sterile Blade, Box of 100, Carbon Steel

Sugg. Retail Price 21.158 FT  
**Promo 13.118 FT**

**10-256-15**

#15C Sterile Blade, Box of 100, Carbon Steel

Sugg. Retail Price 21.158 FT  
**Promo 13.118 FT**

## SUTURES

SPECIAL OFFER

1+1

50%\* PROMOTION  
ON ALL HU-FRIEDY  
SUTURES

DISCOVER MORE ABOUT  
PERMA SHARP SUTURES  
ON OUR DIGITAL CATALOGUE  
SURGICAL CHAPTER PAGE E130

## NYLON NON ABSORBABLE

NEEDLE		Suture Size	CODE
C-3		13 mm	5-0 <b>PSN698N</b>
C-6		18.7 mm	4-0 <b>PSN662N</b>
			5-0 <b>PSN661N</b>

## SILK NON ABSORBABLE

NEEDLE		Suture Size	CODE
C-3		13 mm	4-0 <b>PSN641S</b>
			5-0 <b>PSN640S</b>
			3-0 <b>PSN7772S</b>
C-6		18.7 mm	4-0 <b>PSN683S</b>
			5-0 <b>PSN682S</b>
C-7		24.3 mm	3-0 <b>PSN684S</b>
			2-0 <b>PSN737S</b>
C-9		23 mm	3-0 <b>PSN632S</b>
			4-0 <b>PSN631S</b>
C-17		12 mm	4-0 <b>PSN736S</b>
C-19		18.3 mm	3-0 <b>PSN733S</b>
			4-0 <b>PSN734S</b>
C-22		16.2 mm	5-0 <b>PSN1628S</b>
			6-0 <b>PSN681S</b>
			3-0 <b>PSN7762S</b>
C-43		16.2 mm	4-0 <b>PSN7763S</b>
			5-0 <b>PSN7764S</b>
D-14		17.5 mm	3-0 <b>PSN953S</b>
			4-0 <b>PSN952S</b>
			3-0 <b>PSN872S</b>
T-28		17.5 mm	4-0 <b>PSN871S</b>
			5-0 <b>PSN870S</b>
T-4		22.2 mm	3-0 <b>PSN7343S</b>
			4-0 <b>PSN7342S</b>

## POLYPROPYLENE NON ABSORBABLE

NEEDLE		Suture Size	CODE
C-1		10.7 mm	6-0 <b>PSN8697P</b>
			7-0 <b>PSN8696P</b>
C-3		13 mm	5-0 <b>PSN8698P</b>
			6-0 <b>PSN8695P</b>
C-6		18.7 mm	4-0 <b>PSN8683P</b>
			5-0 <b>PSN8661P</b>
C-4		12.7 mm	5-0 <b>PSN8245P</b>
C-17		12 mm	5-0 <b>PSN8385P</b>
			6-0 <b>PSN8384P</b>
C-22		16.2 mm	5-0 <b>PSN8614P</b>
			6-0 <b>PSN8660P</b>
D-14		17.5 mm	5-0 <b>PSN7773P</b>
T-28		17.5 mm	5-0 <b>PSN8756P</b>

## PGA FAST ABSORPTION

NEEDLE		Suture Size	CODE
C-3		13 mm	4-0 <b>PSN510FA</b>
			5-0 <b>PSN515FA</b>
C-6		18.7 mm	4-0 <b>PSN525FA</b>
			5-0 <b>PSN530FA</b>
C-6		Regular Cut 18.7 mm	3-0 <b>PSN540FA</b>
			4-0 <b>PSN535FA</b>
C-9		23 mm	3-0 <b>PSN585FA</b>


## PGA ABSORBABLE

NEEDLE		Suture Size	CODE
C-1		10.7 mm	5-0 <b>PSN490V</b>
			6-0 <b>PSN489V</b>
			4-0 <b>PSN494V</b>
C-3		13 mm	5-0 <b>PSN493V</b>
			6-0 <b>PSN492V</b>
			3-0 <b>PSN393V</b>
C-6		18.7 mm	4-0 <b>PSN392V</b>
			4-0 <b>PSN397V</b>
C-17		12 mm	4-0 <b>PSN386V</b>
			5-0 <b>PSN385V</b>
C-19		18.3 mm	4-0 <b>PSN724V</b>
			5-0 <b>PSN723V</b>
			4-0 <b>PSN389V</b>
C-22		16.2 mm	5-0 <b>PSN388V</b>
			6-0 <b>PSN387V</b>
D-14		17.5 mm	4-0 <b>PSN271V</b>

## POLYESTER NON ABSORBABLE

NEEDLE		Suture Size	CODE
C-3		13 mm	5-0 <b>PSNR698L</b>
C-6		18.7 mm	3-0 <b>PSNR7772L</b>
			4-0 <b>PSNR683L</b>
			3-0 <b>PSN6953L</b>
D-14		17.5 mm	5-0 <b>PSN7773L</b>





## PERIOSTEALS

SPECIAL OFFER

# -35%

APPLY 35% DISCOUNT  
ON YOUR FAVOURITE  
PERIOSTEALS\*

From 10.413 FT\*\*




DISCOVER THE COMPLETE  
PERIOSTEALS RANGE ON OUR  
DIGITAL CATALOGUE  
SURGICAL CHAPTER - PAGE E92



\* Valid for all range of periosteals, with any kind of handle, except **BLACK LINE**.

\*\* Discount already applied on the cheapest instrument.

## KNIVES



SPECIAL OFFER

# -35%

APPLY 35% DISCOUNT  
ON YOUR FAVOURITE  
PERIODONTAL KNIVES\*

From 13.388 FT\*\*



DISCOVER THE COMPLETE  
PERIODONTAL KNIVES RANGE ON  
OUR DIGITAL CATALOGUE  
SURGICAL CHAPTER - PAGE E71



\* Valid for all range of periodontal knives, with any kind of handle, except **BLACK LINE**.

\*\* Discount already applied on the cheapest instrument.





# COMPOSCULP BY DR. DIETSCHI COMPOSITE INSTRUMENT KIT

## POINTS OF PERFORMANCE

- **Remarkable Restorations:** Highly polished stainless steel working ends resist composite material, ensuring a stunning final restoration.
- **Reliable Hu-Friedy Quality:** Working ends crafted with Immunity Steel®, an optimal blend of carbon and chromium providing strength and resisting corrosion.
- **Color-Coding Indicators:** Each CompoSculp instrument exhibits a color identifier, ensuring the desired instrument is at hand.
- **Lightweight Handles:** Durable ResinEight® and sleek Smooth handle options allow for maximum comfort, reduced hand fatigue and an increased level of tactile control.

### KOL CORNER




Instruments made in collaboration with **Dr. Dietschi**, Opinion Leader Hu-Friedy



**PFIDDDIN**  
Kit contains 5 instruments and 1 IMS Infinity cassette 1/4 DIN Orange

The Kit includes:		CODE
Dietschi Composite DIN Cassette of 5 instruments		<b>IME14DIN53</b>
Composculp # 1/2		<b>PFIDD1/2</b>
Composculp # 3/4		<b>PFIDD3/4</b>
Composculp # 5/6		<b>PFIDD5/6</b>
Composculp # 7/8		<b>PFIDD7/8</b>
Composculp # 9/10		<b>PFIDD9/10</b>


VERSION	CODE	Sugg. Retail Price	Promo
Dr. Dietschi Kit	<b>PFIDDDIN</b>	136.708 FT	<b>102.531 FT</b>



# EXTENDING YOUR INSTRUMENTS LIFE

Our newly upgraded **technical care center**, based in our facilities in Tuttlingen, further increases our goal of **maximising the useful life of our instruments**. This post-sale service allows dental professionals to ensure their instruments always perform at their best and for the longest time possible.

## HU-FRIEDY TECHNICAL CARE CENTER



### EXTENDING YOUR INSTRUMENTS LIFE!

#### SHARPENING

Instruments should be kept sharp and true to their original design. Dental procedures are most effective when using sharp instruments because they reduce fatigue, improve deposit removal, save time, enhance tactile sensitivity, and minimize patient discomfort.

#### PREPARATION

Our experienced service team in Tuttlingen checks, repairs and upgrades your Hu-Friedy Instruments.

#### MAINTENANCE


Regular maintenance extends your instruments life, improves reliability and ensures highest comfort and safety for patient and user.

#### CUSTOMIZATION

You can choose to customize your instrument, engraving the name of your office, its logo, or more simply, by procedural kit, making them easily identifiable for you and your team. Our laser marking technology allows our Technical Care Center Team to offer this.

Contact our experts for more information, please send us your enquiry:  
[HFETechService@hu-friedy.com](mailto:HFETechService@hu-friedy.com)

### SURGICAL







## DOWNLOAD OUR CATALOGUE

For all of us at Hu-Friedy, this is much more than a book that contains the highest quality dental instruments and products. It's a manifestation of the dedication, innovative spirit and hard work of a company that has served and supported the dental community for over 110 years.

[catalog.hu-friedy.eu](http://catalog.hu-friedy.eu)

## FOLLOW US ON:



@HuFriedyEU



@hu\_friedy\_europe



@Hu-Friedy Mfg. Co., LLC.EU

## OUR NEXT GIGS

**EFP Perio Master Clinic**  
6-7/03/2020 - Dublin

**EXPODENTAL**  
12-14/03/2020 - Madrid

**DENTAL EXPO AMSTERDAM 2020**  
19-21/03/2020 - Amsterdam

For more information about  
Hu-Friedy's Total Solutions visit  
**HU-FRIEDY.EU**

PREVENTIVE



RESTORATIVE



ORTHODONTIC




INSTRUMENT  
MANAGEMENT



ENDODONTIC



SURGICAL



POWER



ASK YOUR LOCAL REPRESENTATIVE OR DEALER FOR THE AVAILABILITY IN YOUR COUNTRY

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

All prices indicated are the recommended prices in EU and don't include VAT and other taxes.

Chicago | Milan | Frankfurt | Shanghai | Tokyo | Tuttlingen

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell St. | Chicago, IL 60618 | USA

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. The contents of this promotions including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of Hu-Friedy Mfg. Co., LLC and may not be copied, reproduced or distributed without prior written permission. Hu-Friedy Mfg. Co., LLC. reserves the right to change features of products, images and text at any time without notice. SporeClear is a registered trademark of Altos Medical, LLC. HFL-070/0220

How the best perform

**Hu-Friedy**