



NATIONAL PROMOTION HUNGARY

AUTUMN
2020

Valid from 15th of October
until 31th of December 2020

Introducing...

HuFriedyGroup

The Best In Practice

Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

Service & Support


Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

Proud members of HuFriedyGroup:



RESTORATIVE INSTRUMENTS

CARVERS HU-FRIEDY



CVWI8

Carver Willand
#8 hdl #41

Sugg. Retail

Ft 16.180

Promo

Ft 10.517

CVIPC

Carver
hdl #41 Interproximal


Sugg. Retail

Ft 16.180

Promo

Ft 10.517

EXCAVATORS HU-FRIEDY



EXC17WH

Excavator #17WH
hdl #41

Sugg. Retail

Ft 16.491

Promo

Ft 10.719

EXC18WH

Excavator #18WH
hdl #41 heavy

Sugg. Retail

Ft 16.491

Promo

Ft 10.719

PLASTIC FILLING HU-FRIEDY



CIGFT16

Plastic Filling #1
Instrument Goldstein
hdl #6 Flexi-Thin

Sugg. Retail

Ft 19.025

Promo

Ft 11.415

CIGFT26

Plastic Filling #2
Instrument Goldstein
hdl #6 Flexi-Thin

Sugg. Retail

Ft 19.025

Promo

Ft 11.415



GCPS6


Gingival Cord Packer
#S6 hdl #41, smooth

Sugg. Retail

Ft 18.669

Promo

Ft 11.201



BB216

Burnisher #21
hdl #6

Sugg. Retail

Ft 16.847

Promo

Ft 10.108

BB186

Burnisher #18
hdl #6

Sugg. Retail

Ft 16.847

Promo


Ft 10.108

RESTORATIVE INSTRUMENTS

AKRO-FLEX™ COMPOSITE INSTRUMENT HU-FRIEDY


CREATE BEAUTIFUL RESTORATIONS

Akro-Flex™ is an extremely flexible composite instrument which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations. **The inclusion of a ductile material allows Akro-Flex™ to perform as a solid brush.**



HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease. The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex™ acts as a solid brush. The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.



ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.

AKRO-FLEX

Akro-Flex™
Composite Instrument

Sugg. Retail


Ft 32.004

Promo

Ft 28.803



WATCH OUR
TUTORIAL VIDEO
FOR AKRO-FLEX™
ON YOUTUBE



Patent pending.

RESTORATIVE INSTRUMENTS BY WERESTORE.IT[†]

CALIBRA™ COMPOSITE INSTRUMENT HU-FRIEDY

FOR CALIBRATED RESTORATIONS

Calibra™ is an **XTS (AlTiN)^{††}** coated composite instrument which helps to provide the proper anatomical curvature to the dentin material layer of a restoration and leaves the ideal space for the final layer of enamel material. The **Calibra™ “IN”** point is used to define the dentinal skeleton, while the **Calibra™ “OUT”** point helps to reduce the chances of **excess composite** material being placed during the final layer of the restoration. Calibra™ acts as a **guide**, helping to prevent over-contouring at the preparation margin while assisting in reducing excess composite thickness build up as well as **modeling** the final layer of composite material.

TNCALIBRA

XTS PFI WeRestore.it
Calibra™

NEW

Sugg. Retail

Promo

Ft 28.092

Ft 19.665

HOW CALIBRA™ IN WORKS



HOW CALIBRA™ OUT WORKS



Calibra™ consists of two terminals: “IN” which defines the dentinal skeleton. “OUT” which allows to verify the thickness and shape of the last layer without excess. This instrument is in fact both the first caliper ever produced to check the thickness (a practical alternative to the sagittal silicone index) and the first modeling tool ever produced for the last layer.

[†] WeRestore.it is not owned by or affiliated with Hu-Friedy. ^{††} Data on file, available upon request.



KOL CORNER

Instruments made
in collaboration with
Dr. Gaetano Paolone
and **Dr. Salvatore**
Scolavino, Opinion
leaders Hu-Friedy.

CALIBRA™
OUT

CALIBRA™
IN



CALIBRA™ BY
WERESTORE.IT
WATCH OUR
PRACTICAL
VIDEO

DIAGNOSTIC INSTRUMENTS HU-FRIEDY



EXS3A6

Explorer #3A
hdl #6

Sugg. Retail Ft 8.757

Promo Ft 6.567



EXS6A6

Explorer #6A
hdl #6

Sugg. Retail Ft 8.757

Promo Ft 6.567



EXS236

Explorer Shepherds
Hook #23 hdl #6

Sugg. Retail Ft 8.757

Promo Ft 6.567



EXTU176

Explorer Wilkins/Tufts
#TU17 hdl #6

Sugg. Retail Ft 8.757

Promo Ft 6.567



EXD11/12A6

Explorer #11/12 hdl #6
After Five

Sugg. Retail Ft 11.824

Promo Ft 8.868

1-2-3....15



PCPUNC156

Probe UNC #UNC15
hdl #6

Sugg. Retail Ft 15.291

Promo Ft 11.468

3.5....11.5



PCP11.5B6

Probe #11.5B
hdl #6

Sugg. Retail Ft 15.291

Promo Ft 11.468



EXTU17/236

Explorer Wilkins/Tufts
#TU17/23 hdl #6

Sugg. Retail Ft 11.824

Promo Ft 8.868

3-6-9-12



PCP126

Probe #12
hdl #6

Sugg. Retail Ft 15.291

Promo Ft 11.468

1-2-3-5-7-8-9-10



PQOW6

Probe Michigan/Williams
hdl #6

Sugg. Retail Ft 13.735

Promo Ft 10.301



P15/11.5B6

Probe #15/11.5B
hdl #6

Sugg. Retail Ft 24.892

Promo Ft 18.669

3-6-9-12



PQ2N6

Birfurcation Probe
Nabers #2N hdl #6

Sugg. Retail Ft 24.892

Promo Ft 18.669

EVEREDGE™ 2.0 INSTRUMENTS

EVEREDGE™ 2.0 SCALERS AND CURETTES

NO SHARPENING... DO YOU REALLY BELIEVE IT?



PROPHY MAINTENANCE KIT DIN HU-FRIEDY


IMS™ Infinity Cassette with 5 EverEdge 2.0 Instruments:

IMEDPROPHYM

DESCRIPTION	CODE
Scaler Neville #1 anterior Sickle	SCNEVI19E2
Scaler #204SD Sickle	S204SD9E2
Curette Columbia #4R/4L	SC4R/4L9E2
Curette Gracey Labriola-Cortellini #13/14 DE	SG13/149LCE2
Curette Gracey Labriola-Cortellini #11/12	SG11/129LCE2

Sugg. Retail **Promo**

Ft 160.331 **Ft 104.215**



For more information
on **EVEREDGE™ 2.0**
visit our site hu-friedy.eu



STANDARD GRACEY CURETTE HU-FRIEDY

Description	Handle 9	Resin 8 Color
1/2 Standard Gracey	SG1/29IE2	SG1/2C8E2
7/8 Standard Gracey	SG7/897E2	SG7/8C8E2
11/12 Standard Gracey	SG11/1293E2	SG11/12C8E2
13/14 Standard Gracey	SG13/1498E2	SG13/14C8E2

MINI FIVE GRACEY CURETTES HU-FRIEDY

Description	Handle 9	Resin 8 Color
1/2 Mini Five Gracey	SAS1/29IE2	SAS1/2C8E2
7/8 Mini Five Gracey	SAS7/897E2	SAS7/8C8E2
11/12 Mini Five Gracey	SAS11/1293E2	SAS11/12C8E2
13/14 Mini Five Gracey	SAS13/1498E2	SAS13/14C8E2

UNIVERSAL CURETTES HU-FRIEDY

Description	Handle 9	Resin 8 Color
1/2 Langer	SL1/29E2	SL1/2C8E2
3/4 Langer	SL3/49E2	SL3/4C8E2
5/6 Langer	SL5/69E2	SL5/6C8E2
17/18 Langer	SL17/189E2	SL17/18C8E2

SICKLE SCALERS HU-FRIEDY

Description	Handle 9	Resin 8 Color
H6/H7 Scaler	SH6/79E2	SH6/7C8E2
204S Scaler	S204S9E2	S204SC8E2
204SD Scaler	S204SD9E2	S204SDC8E2
1 Nevi Anterior	SCNEVI19E2	SCNEVI1C8E2
2 Nevi Posterior	SCNEVI29E2	SCNEVI2C8E2
3 Nevi Scaler Posterior	SCNEVI39E2	SCNEVI3C8E2
4 Nevi Scaler Posterior	SCNEVI49E2	SCNEVI4C8E2

Valid for all **EVEREDGE 2.0** range



Handle 9

Sugg. Retail

Ft 24.270

Promo

Ft 18.202




Handle 8 Color

Sugg. Retail

Ft 23.070

Promo

Ft 17.302



*Products tested: Hu-Friedy EverEdge™ 2.0 H7, American Eagle XP H6-7, PDT H6-7, LMEM/EN H6-7, and G M Hartzell H6-7. For more information, contact David Platt, Hu-Friedy Mfg. Co., LLC, 3232 N. Rockwell St., Chicago, IL 60618.

SCISSORS & NEEDLE HOLDERS HU-FRIEDY

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. All hinged instruments should be sterilised in the open position.



NH5038

PermaSharp™ Needle Holder Crile-Wood #5038, 15 cm

Sugg. Retail Ft 63.430

Promo Ft 44.401



NHCW

Needle Holder Crile-Wood grooved, 3-0 - 6-0, 15cm

Sugg. Retail Ft 19.025

Promo Ft 13.317



NMH

Needle Holder Mayo-Hegar grooved, 1-0 - 4-0 Needle, 16cm

Sugg. Retail Ft 20.892

Promo Ft 14.624



NH5024CORT

Needle Holder Cortellini Micro Castro #5024 round handle, straight

Sugg. Retail Ft 184.245

Promo Ft 128.972



NH5024SIM

PermaSharp™ Needle Holder Simion #5024R round handle, straight, 18cm

Sugg. Retail Ft 184.245

Promo Ft 128.972



S17SC

SuperCut™ Scissors Iris, #17, straight, 11,5 cm

Sugg. Retail Ft 36.982

Promo Ft 25.888



S18SC

SuperCut™ Scissors Iris #18 fine angled, 11,5cm

Sugg. Retail Ft 36.982

Promo Ft 25.888



S14SC

SuperCut™ Scissors LaGrange #14 double angled, 11,5cm

Sugg. Retail Ft 47.250

Promo Ft 33.075



S16

SuperCut™ Scissors Iris #18 fine angled, 11,5 cm **Black Line**

Sugg. Retail Ft 21.203

Promo Ft 14.842



SCORTSC

SuperCut™ Scissors Cortellini curved, 12cm

Sugg. Retail Ft 46.050

Promo Ft 32.235



S5083

Scissors Iris #5083 angled, PermaSharp™, 11,5 cm

Sugg. Retail Ft 83.344

Promo Ft 58.341

SUPER CUT™

Cutting blade honed/serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

PERIOSTEALS


PERIOSTEAL HU-FRIEDY



P24G
Periosteal #24
long handle

Sugg. Retail
Ft 33.604

Promo
Ft 25.203



PPR3S
Periosteal Prichard
#3 hdl #522 short

Sugg. Retail
Ft 39.827

Promo
Ft 29.870



PPBUSER6
Periosteal Buser
hdl #6

Sugg. Retail
Ft 38.583


Promo
Ft 28.937



P146
Periosteal
Goldman-Fox
#14 hdl #6

Sugg. Retail
Ft 34.358


Promo
Ft 25.903



CL84
Surgical Curette
Lucas #84, hdl #522
2,2mm

Sugg. Retail
Ft 35.782


Promo
Ft 25.048



CL85
Surgical Curette
Lucas #85, hdl #522
2,5mm

Sugg. Retail
Ft 35.782


Promo
Ft 25.048



CL86
Surgical Curette
Lucas #86, hdl #522
2,8mm

Sugg. Retail
Ft 35.782

Promo
Ft 25.048



CL87
Surgical Curette
Lucas #87, hdl #522
3,5mm

Sugg. Retail
Ft 16.847

Promo
Ft 10.108



CM2/4
Surgical Curette Molt
#2/4, hdl #522

Sugg. Retail
Ft 36.093

Promo
Ft 25.265

SURGICAL CURETTES HU-FRIEDY

ENDODONTIC INSTRUMENTS

EXCEPTIONAL RESTORATION STARTS WITH EXCELLENT INSULATION



RDPA
Rubber Dam Punch
Ainsworth
17 cm

Sugg. Retail Ft 95.745

Promo Ft 76.596



RDF4E
Rubber Dam Clamp
Forceps #4
European

Sugg. Retail Ft 56.896

Promo Ft 45.517



RDCSET7
Rubber Dam Clamps
Starter Kit with 7 clamps

Sugg. Retail Ft 46.361

Promo Ft 37.089



RDCOB
IMS™ Rubber Dam Clamp
Organisation Board
for 8 clamps

Sugg. Retail Ft 32.937


Promo Ft 26.350



RDAF6
Rubber Dam Frame
Adult, 15cm

Sugg. Retail Ft 15.602

Promo Ft 12.482



RDCF5
Rubber Dam Frame
Child
13cm

Sugg. Retail Ft 15.602

Promo Ft 12.482



IMS-1273
IMS™ Accessory Box
large

Sugg. Retail Ft 23.959


Promo Ft 19.167



IMS-1271
IMS™ Accessory Box
small

Sugg. Retail Ft 17.736

Promo Ft 14.188



IMEDRDAM
IMS™ Infinity DIN Collection Rubber Dam Kit


Sugg. Retail Ft 326.219

Promo Ft 260.975

INSTRUMENT MANAGEMENT SYSTEM

IMS™ INFINITY SERIES CASSETTES

IMS Infinity Series Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, cassettes can be specialized to meet your needs.



1/4 DIN CASSETTE HU-FRIEDY IME14DIN5_

IMS Cassette 1/4 DIN for 5 instruments
Dimensions: 72 x 184 x 34 mm
(without instruments)

Sugg. Retail

Ft 54.451

Promo

Ft 43.561



1/2 DIN CASSETTE HU-FRIEDY IME12DIN1_

IMS Cassette 1/2 DIN for 10 instruments
Dimensions: 152 x 184 x 34 mm
(without instruments)

Sugg. Retail

Ft 60.941

Promo

Ft 48.753



DIN SURGERY CASSETTE HU-FRIEDY IMEDIN18_S

DIN Cassette for 18 instruments
with an accessory area + 4 clips
Dimensions: 290 x 184 x 34 mm
(without instruments)

Sugg. Retail

Ft 88.589

Promo

Ft 70.871

IMS SIGNATURE SERIES™ TUBS HU-FRIEDY IMS-140_

Together, procedure tubs and cassettes create a more efficient and organised system for managing instruments and consumable products.

Includes:

- Signature Series Tub
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover



Sugg. Retail

Ft 34.671

Promo

Ft 27.737

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace " _ " with the number)

Red	●	0	Lavender	●	6
Grey	●	1	Black	●	7
White	●	2	Blue	●	8
Orange	●	3	Green	●	9
Purple	●	4	Ocean Blue	●	11
Yellow	●	5			

INFECTION CONTROL

DENTAPURE™ CARTRIDGE CROSSTEX**

REVOLUTIONIZE YOUR WATERLINE MAINTENANCE! DAILY WATERLINE MAINTENANCE BECOMES A SIMPLE, ANNUAL ROUTINE**

Benefits for both independent water bottle systems & municipal systems plumbed to city water:

- NO routine shocking
- NO distilled water required
- NO harsh chemicals
- NO silver
- NO special disposal requirements
- NO allergenic iodine proteins¹
- NO concerns with dental materials and bonding²
- NO concerns with dental unit corrosion or etching³
- An effective method of treating municipally plumbed dental units

For more information on how the multi-award-winning DentaPure™ Cartridge can help simplify your DUWL maintenance protocols, visit: DentaPure.Crosstex.com/EU

DESCRIPTION	QUANTITY	CODE	Sugg. Retail	Promo
Independent Water Bottle Cartridge	1 Each	DPI365B10	Ft 111.125	Ft 66.675
Municipal Cartridge	1 Each	DPI365M10	Ft 111.125	Ft 66.675
Iodine Test Strips	50 Test Strips/Bottle	DPTTEST	Ft 15.558	Ft 12.290

Municipal Cartridge DPI365M10

For waterlines plumbed directly to city water. Installs onto incoming lines with included fittings.



ONE
DENTAPURE™
CARTRIDGE
= **365 DAYS**
COMPLIANT,[†]
DENTAL UNIT
WATER*

Independent Water Bottle Cartridge DPI365B10

For independent water bottles. Easily installed by office staff onto existing pickup tube with included fitting.

[†]Compliant = Compliant with EPA potable water standard <500 CFU/mL. ^{††}Or, 240L of water if usage records are kept. ¹Schabelman E., Witting M. The relationship of radiocontrast, iodine, and seafood allergies: a medical myth exposed, J Emerg Med. 2010 Nov;39(5):701-7. doi: 10.1016/j.jemermed.2009.10.014. Epub 2010 Jan 4. ²Puttaiah R. Effects of Low Grade Iodine in Dental Unit Waterlines On Shear Bond Strength of a Dentin Bonding Agent, Baylor College of Dentistry. ³Data on file at Crosstex.

BACTRYL™ DISINFECTANT CROSSTEX**

READY-TO-USE SPRAY AND WIPES FOR THE CLEANING AND DISINFECTION OF SURFACES AND DENTAL MEDICAL DEVICES

- Ready-to-use, broad spectrum detergent and disinfectant (One Step)
- Fully tested and compliant to most rigid EU Standards (EN14885)
- Wipes are Compliant to new compulsory EN 16615 in clean and dirty conditions.
- VAH-listed Disinfectant
- Molecule effective against gram positive and gram negative bacteria, mycobacteria, bacteria, viruses (HIV, HBV, HCV) and fungi. Effective on all enveloped viruses.
- Fully compatible with surfaces such as worktops, lamps, trolleys, dental units, trays, sinks, instruments and furnishings of dental surgery with long-lasting residual bacteriostatic action
- Safe for dental professionals and patients
- Available as a spray or wipe

500046

4 x 1000 ml Bottles
+ 2 x Spray
Triggers/Case

Sugg. Retail Ft 28.893

Promo Ft 17.336



500047

6 x 110 Wipes/
Container

Sugg. Retail Ft 45.339

Promo Ft 36.271



Cantel Medical Italy Srl is the legal manufacturer

PROCEDURAL KITS

BASIC IMPLANT SET OMNIA*

The Basic Implant Set has been designed for all routine surgical activities where a surgical protocol ensuring an adequate degree of sterility and guaranteeing the safety of operators and patients is needed. Due to its reasonable price, it is especially fit for all surgical procedures with immediate loading or for operations, such as re-opening of implants, flap surgery, operations on minor surgical wounds.



12.S5621.00

Basic Implant Set

DESCRIPTION

1	Surgically folded Special Plus model gown (size L) and 2 hand towels wrapped in medical paper	SMS Microfibre	Light blue
1	Gown with cotton cuffs	NWF	Light blue
2	Surgical caps with fixed laces	Absorbent NWF	Light blue
1	Round surgical cap for the patient	Absorbent NWF	Light blue
1	Absorbent/liquid proof drape 50x75 cm with self-adhesive [†] side	NWF + Polyethylene	Light blue
1	Absorbent/liquid proof drape 75x90 cm with adhesive [†] U-shaped cut 11x9 cm	NWF + Polyethylene	Light blue
1	Omnisleeve sheath 120x7 cm	PE	Transparent
1	Disposal bag for contaminated material	PVC	Red

[†] **Bold numbers** indicate the **adhesive** side

Sugg. Retail

Promo

Ft 9.779

Ft 8.801

PERIODONTAL SET OMNIA*

Sterile set for periodontal surgery or small oral surgery. It contains patient drapes and surface drapes. It includes sheaths for hand pieces, adhesive films for handles, keypads and folded compresses. The set is packed in a double packaging. The practical peel-off label allows easy traceability of the operation.



12.T0013.00

Periodontal Set

DESCRIPTION

1	Adhesive [†] U shaped 11 x 9 cm absorbent/water-proof drape cm 75 x 90 cm	Non woven +Polyethylene	Blue
1	Absorbent/water proof drape 75 x 90 cm	Non woven +Polyethylene	Blue
2	Absorbent/water proof drapes 50 x 75 cm	Non woven +Polyethylene	Blue
2	Tubing sleeves 120 x 7 cm with fixing adhesives	LDPE	Transparent
2	Folded compresses 5 x 5 cm, 4-ply	Non woven	White

[†] **Bold numbers** indicate the **adhesive** side

Sugg. Retail

Promo

Ft 6.223

Ft 5.601

IMPLANTOLOGY SET OMNIA*

Complete sterile set conceived for all Implantology, Parodontology and Maxillo Facial surgeries. The set is composed of patient drape, surface drapes and tubing sleeves. The set is packed in sterile wrapping.

The box contains 10 sets.



12.T0202.00

Implantology Set

DESCRIPTION

1	Surgical drape 60 x 60 cm	SMS	Light blue
1	Surgical Drape 150 x 150 cm with adhesive hole	NWF/TNT	Light blue
2	Surgical Drape 50 x 75 cm	PP+PE	Light blue
2	Adesivi medicali 9 x 51 cm	SPPB	Light blue
2	Adhesive Velcro Tape 14 x 2,5 cm	Velcro	White
2	2 Tube Cover 6 x 120 cm	PE	Transparent

Sugg. Retail

Promo


Ft 68.898

Ft 62.008

SURGERY: THE ESSENTIALS

SUCTION TIPS OMNIA*

Omnia disposable suction tips for surgical aspirators are made of medical grade plastic material. Suction tips ensure a targeted and safe aspiration during surgical operations. They are available with differently shaped hand-grips and terminal ends according to the particular needs of the dental surgeon, and they are supplied without connection tubing.



CODE	DESCRIPTION	Sugg. Retail	Promo
32.F5051.00	Aspirator, 2,50 m, with suction tip with ergonomic hand-grip, and that is also enabled for the autologous bone harvesting trap (each pack contains 1 adapter Adat. 05), PVC, 5 pcs.	Ft 17.780	Ft 14.224
32.F6565.00	Omniasurg model sterile suction tip (each pack contains 1 adapter Adat. 08) (tip length: 20 cm), PVC, 20 pcs.	Ft 11.113	Ft 8.890
32.F2065.00	Surgical aspirator, 2,20 m, with Omniasurg suction tip (each pack contains 1 adapter Adat. 05), PVC, 10 pcs.	Ft 20.003	Ft 16.002



VISUAL OMNIA*

The Visual is an individual protection device available to ensure your safety. Its wide polyethylene fog preventing shield and the polyurethane insert make this device also fit for use by people wearing glasses, 12 pcs.

10.M0200.00

Visual

Sugg. Retail	Promo
Ft 13.335	Ft 12.002



SHELL COVER OMNIA*

The SHELL COVER gown cover is perfect for routine practice that create high aerosol dispersion levels. Material: Medical grade LDPE waterproof, with high insulation standard. The additional Thumb-Up system helps keeping the gown fixed while worn.

30.D1002.00

Shell Cover - Not Sterile 100 pcs

Sugg. Retail	Promo
Ft 106.680	Ft 96.012

22.D1002.00

Shell Cover - Sterile 50 pcs

Sugg. Retail	Promo
Ft 83.344	Ft 75.009

IRRIGATION SETS OMNIA*

PHYSIODISPENSER AND PIEZOELECTRIC SURGERY

- Proper irrigation during the creation of the implant site is a fundamental element for successful bone integration of the implant and the overall outcome of the operation.
- Wide range of mechanical irrigation lines which are compatible with most implant motor models. Irrigation systems include a vented spike, a roller or clamp for flow control, fittings, and an integrated peristaltic pump section.
- Convenient terminal cap makes the connection of the irrigation line to the hand piece easy to achieve. The cap allows you to connect tubes to the wide variety of cooling needle models available.
- The mechanical irrigation systems offer the following options:
 - a. Internal irrigation only
 - b. External irrigation only
 - c. Internal and external irrigation



CODE (10 Pcs, par boîte)	DESCRIPTION	Sugg. Retail	Promo
32.F1201.00	Irrigation Tubing for NSK Surgic XT, Surgic XT Plus, Osseocision, Vario Surg, Kyocera Implantor Neo, Aseptico Aeu 6000 series, Aseptico Aeu 7000 series, Bonart Piezo surgical system	Ft 48.895	Ft 39.116
32.F0134.00	NSK Surgic XT, Surgic XT Plus, Osseocision®, VarioSurg, SurgicPro, Kyocera Implantor Neo, Bonart Piezo surgical system, Dentium iCTmotor, NSK 3i Osseocision, Delfine Intra-Lock, 1/2Y	Ft 34.227	Ft 27.381
32.F0139.00	Nobel Biocare80/ W&H Implantmed/ Implant Med new design 2005/ Implantmed SI 923/ Astratech Implant SI 923/ Surgical Motor System/Ace/ Straumann/InnovaCorp. 3i/Centerpulse/Lifecore/ DU 900/FRIOS® UNIT S/ FRIOS® UNIT S/I, Osseoset 100/Osseoset 200/Osseocare, Elcomed SA 310/SA2006, Piezomed, W&H SI 1023, 1/2Y	Ft 34.227	Ft 27.381
32.F0170.00	Bien Air Chiropro 980 / EMS PIEZON® MASTER surgery / Esacrom Surgisonic I / UBS-Ultrasonic bone surgery/ Surgybone - Acrobone - Easybone Silfradent, BioSAFin Easy Surgery, MyTUTOR, Heracle Sweden&Martina, Cefla dental Chair, Osseodoc BienAir, 1/2Y	Ft 34.227	Ft 27.381
32.F1204.00	Nobel Biocare80/ W&H Implantmed/ Implant Med new design 2005/ Implantmed SI 923/ Astratech Implant SI 923/ Surgical Motor System/Ace/ Straumann/InnovaCorp. 3i/Centerpulse/Lifecore/ DU 900/FRIOS® UNIT S/ FRIOS® UNIT S/I, Drilltech/Osseoset 100/Osseoset 200/Osseocare, 1/2Y	Ft 48.895	Ft 39.116



**FOR MORE
INFORMATION
ON OUR
IRRIGATION
SET.**



For more information on all available irrigation lines, see :

<https://www.omniaspa.eu/en/products/oral-irrigation-kits.html>



@HuFriedyEU



@hu_friedy_europe



@Hu-Friedy Mfg. Co., LLC.EU



@OmniaSPA



@Omnia S.p.A.

*Products provided by Omnia. ** Products provided by Crosstex. All other products provided by Hu-Friedy.

For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit Hu-Friedy.eu

For more information on Crosstex products, call +39 05 24 52 74 53 or visit crosstex.com

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit omniaspa.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. All prices indicated are intended VAT included. All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.

©2020 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/1020