

NATIONAL PROMOTION HUNGARY

SPRING/SUMMER 2021

Valid from 15th of May
until 31st of July 2021

Introducing...

HuFriedyGroup

The Best In Practice

Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

Service & Support

Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

Proud members of HuFriedyGroup:



DIAGNOSTIC INSTRUMENTS HU-FRIEDY®



EXD5
Explorer #5 hdl #31

Sugg. Retail Ft 12.264

Promo Ft 7.972



EXD3CH
Explorer #3 cow horn

Sugg. Retail Ft 12.264

Promo Ft 7.972



EXD11/126
Explorer Dominion #11/12
hdl #6

Sugg. Retail Ft 12.875

Promo Ft 8.369



1-2-3...15
PCPUNC157
Probe UNC #UNC15 hdl
#7 Qulix 1-2-3-...15

Sugg. Retail Ft 17.974

Promo Ft 11.683



3-6-8-11
PCP11
Probe #11 hdl #30 Qulix
3-6-8-11

Sugg. Retail Ft 15.248

Promo Ft 9.911



3-6-9-12
PQ2N6
Bifurcation Probe
Nabers #2N hdl #6

Sugg. Retail Ft 27.113

Promo Ft 17.624



DPU17
Dressing Pliers Universal
#17

Sugg. Retail Ft 8.811

Promo Ft 5.286



DP2
Dressing Pliers College
#2

Sugg. Retail Ft 10.573

Promo Ft 6.872



DP18
Dressing Pliers #18

Sugg. Retail Ft 11.583

Promo Ft 6.950



MIR5HD12E
High Definition Mirror #5
single-sided, European,
12 pcs./pack

Sugg. Retail Ft 42.314

Promo Ft 25.389



MIR4HD12E
High Definition Mirror #4
single-sided, European,
12 pcs./pack

Sugg. Retail Ft 42.314

Promo Ft 25.389




MHE6
Mirror Handle hdl #6
European Thread

Sugg. Retail Ft 10.855

Promo Ft 6.513

NO SHARPENING... DO YOU REALLY BELIEVE IT?



SG5/695E2

Curette Gracey #5/6 hdl #9 oral/labial, yellow

Sugg. Retail

Ft 26.432

Promo

Ft 17.181

Valid for all #9 HANDLE




EVEREDGE™ 2.0 INSTRUMENTS SCALERS AND CURETTES

SCALERS EVEREDGE™ WERE SHARP..

SCALERS EVEREDGE™ 2.0 ARE EVEN MORE SHARPENED.

Using an enhanced finishing process, Hu-Friedy's professional artisans are creating edges that are scientifically measured to be over 60% sharper than the next closest competitor. Less force means easier, more comfortable scaling for the clinician and patient.

*Products tested: Hu-Friedy EverEdge™ 2.0 H7, American Eagle XP H6-7, PDT H6-7, LM EM/EN H6-7, and G M Hartzell H6-7. Hu-Friedy cooperated with CATRA, an independent specialized blade laboratory, to conduct this testing. For more information, contact David Platt, Hu-Friedy Mfg. Co., LLC, 3232 N. Rockwell St, Chicago, IL 60618.



For more information on **EVEREDGE™ 2.0** Instruments, visit our site [HuFriedy.eu](https://www.HuFriedy.eu)





IMP204SDT
Implant Scaler Titanium
#204SD

Sugg. Retail Ft 31.178

Promo Ft 20.266



IMPLG1/2T
Implant Scaler Titanium
#Langer 1/2

Sugg. Retail Ft 31.178

Promo Ft 20.266



IMPM1/2T
Implant Scaler Titanium
#Mini-Five 1/2

Sugg. Retail Ft 31.178

Promo Ft 20.266



IMPM13/14T
Implant Scaler Titanium
#Mini-Five 13/14

Sugg. Retail Ft 31.178

Promo Ft 20.266



IMPM11/12T
Implant Scaler Titanium
#Mini-Five 11/12

Sugg. Retail Ft 31.178

Promo Ft 20.266



IMEDIMPLM
Titanium Implant Maintenance Kit
(IMS™ Infinity Series™ Cassette 1/4 DIN, purple + 5 Instruments)

Sugg. Retail Ft 186.104

Promo Ft 120.968

RESTORATIVE INSTRUMENTS HU-FRIEDY®



EXC17H
Excavator #17 hdl #41 heavy

Sugg. Retail Ft 17.974

Promo Ft 10.784



EXC18H
Excavator #18H hdl #41 heavy

Sugg. Retail Ft 17.974

Promo Ft 10.784



EXC19H
Excavator #19 hdl #41 heavy

Sugg. Retail Ft 17.974

Promo Ft 10.784



CIGFT16
Plastic Filling Instrument
Goldstein #1, hdl #6 Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFT26
Plastic Filling Instrument
Goldstein #2, hdl #6 Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFT36
Plastic Filling Instrument
Goldstein #3, hdl #6 Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFT46
Plastic Filling Instrument
Goldstein #4, hdl #6 Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFT5
Plastic Filling Instrument
Goldstein #5, hdl #41 Flexi-Thin

Sugg. Retail Ft 17.974

Promo Ft 10.784



CIGFT6
Plastic Filling Instrument
Goldstein #6, hdl #41 Flexi-Thin

Sugg. Retail Ft 17.974

Promo Ft 10.784



CIGFTMIN16
Plastic Filling Instrument
Goldstein #1, hdl #6
mini Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFTMIN36
Plastic Filling Instrument
Goldstein #3, hdl #6
mini Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CIGFTMIN46
Plastic Filling Instrument
Goldstein #4, hdl #6
mini Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434


RESTORATIVE INSTRUMENTS HU-FRIEDY®

AKRO-FLEX™ COMPOSITE INSTRUMENT

CREATE BEAUTIFUL RESTORATIONS


Akro-Flex™ Composite Instrument is extremely flexible which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations.

The inclusion of a ductile material allows Akro-Flex™ Composite Instrument to perform as a solid brush.



HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease. The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex™ Composite Instrument acts as a solid brush. The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.



ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.

AKRO-FLEX

Akro-Flex™
Composite Instrument

Sugg. Retail

Ft 34.867

Promo

Ft 24.407



WATCH OUR
TUTORIAL VIDEO FOR AKRO-FLEX™
COMPOSITE INSTRUMENT ON YOUTUBE



Patent pending.

RESTORATIVE INSTRUMENTS BY WERESTORE.IT[†]



KOL CORNER

Instruments made in collaboration with **Dr. Gaetano Paolone and Dr. Salvatore Scolavino**, Opinion leaders Hu-Friedy.

CALIBRA™ COMPOSITE INSTRUMENT HU-FRIEDY®

FOR CALIBRATED RESTORATIONS

Calibra™ instrument is an **XTS (AlTiN)^{††}** coated composite instrument which helps to provide the proper anatomical curvature to the dentin material layer of a restoration and leaves the ideal space for the final layer of enamel material. The **Calibra™ “IN”** point is used to define the dentinal skeleton, while the **Calibra™ “OUT”** point helps to reduce the chances of excess composite material being placed during the final layer of the restoration. Calibra Composite Instrument acts as a **guide**, helping to prevent over-contouring at the preparation margin while assisting in reducing excess composite thickness build up as well as **modeling** the final layer of composite material.

TNCALIBRA

XTS PFI WeRestore.it
Calibra™ Instrument

NEW


Sugg. Retail

Promo


Ft 30.168

Ft 21.117

HOW
CALIBRA™ IN
POINT WORKS



HOW
CALIBRA™ OUT
POINT WORKS



Calibra™ Composite Instrument consists of two terminals: “IN” which defines the dentinal skeleton. “OUT” which allows to verify the thickness and shape of the last layer without excess. This instrument is a caliper to check the thickness (a practical alternative to the sagittal silicone index) and a modelling tool for the last layer.












CALIBRA™
COMPOSITE
INSTRUMENT BY
WERESTORE.IT
WATCH OUR
PRACTICAL VIDEO

[†] WeRestore.it is not owned by or affiliated with Hu-Friedy. ^{††} Data on file, available upon request.

RESTORATIVE INSTRUMENTS BY WERESTORE.IT

HU-FRIEDY 3SEENTIAL KIT HU-FRIEDY®

 <p>ANTERIOR</p>   <p>TNCIPAOX XTS PFI Essential WeRestore.it Anterior</p> <p>Sugg. Retail Ft 23.730 Promo Ft 15.424</p>	 <p>POSTERIOR</p>   <p>TNEXBKR1X XTS PFI Essential WeRestore.it Posterior</p> <p>Sugg. Retail Ft 23.730 Promo Ft 15.424</p>	 <p>SPATULA</p>   <p>CVKOT1X XTS PFI Essential WeRestore.it Spatula</p> <p>Sugg. Retail Ft 23.730 Promo Ft 15.424</p>
--	---	---


THE ESSENTIAL RESTORATIVE KIT INCLUDING "PLUG&SCULPT" INSTRUMENT

The Hu-Friedy 3SEENTIAL KIT by Werestore.it is the most straightforward kit of restorative ever made, with just 3 essential instruments at once: Anterior, Posterior and Spatula, from the Hu-Friedy Black Line Collection. The three instruments are synergic and interchangeable and the word "Essential" refers to the easiness of all restorative procedures focused on the essential steps to create high esthetic composite restorations while using this kit.

The revolutionary concept of "Plug&Sculpt" of the Posterior is very simple: one tip is used to plug the composite into the cavity and the other one can sculpt composite in additive and subtractive modelling techniques. The Anterior, instead, is called by the doctors "solid brush" since it models and spreads the composite very easily, just like a brush.

RESTOREITDIN

Set of 3 instruments and cassette



Ce kit comprend:

IMS™ Infinity Series™ Cassette 1/4 DIN, 5 instruments, black	IME14DIN57
XTS PFI Essential WeRestore.it Spatula	CVKOT1X
XTS PFI Essential WeRestore.it Posterior	TNEXBKR1X
XTS PFI Essential WeRestore.it Anterior	TNCIPAOX

DESCRIPTION	CODE	Sugg. Retail	Promo
3SEENTIAL Kit	RESTOREITDIN	Ft 116.065	Ft 75.442

LUXATING ELEVATORS HU-FRIEDY®



EL2S

Luxating Elevator hdl #510,
straight 2 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942



EL2C

Luxating Elevator hdl #510,
curved 2 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942



EL3S

Luxating Elevator hdl #510,
straight 3 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942



EL3C

Luxating Elevator hdl #510,
curved 3 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942



EL4S

Luxating Elevator hdl #510,
straight 4 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942



EL4C

Luxating Elevator hdl #510,
curved 4 mm

Sugg. Retail Ft 47.060

Promo Ft 32.942

SURGICAL CURETTES HU-FRIEDY®



CL846

Surgical Curette Lucas #84,
hdl #6 2,2 mm

Sugg. Retail Ft 39.988

Promo Ft 23.993



CL856

Surgical Curette Lucas #85,
hdl #6 2,8 mm

Sugg. Retail Ft 39.988

Promo Ft 23.993



CL866

Surgical Curette Lucas #86,
hdl #6 2,8 mm

Sugg. Retail Ft 39.988

Promo Ft 23.993



CL876

Surgical Curette Lucas #87,
hdl #6 3,5 mm

Sugg. Retail Ft 39.988

Promo Ft 23.993



CM2/4S6

Surgical Curette Molt #2/4,
short hdl #6

Sugg. Retail Ft 40.270

Promo Ft 24.162

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY® Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



NEEDLE HOLDERS, SCISSORS & TISSUE PLIERS HU-FRIEDY®

NEEDLE HOLDERS HU-FRIEDY®



NHMH
Needle Holder Mayo-Hegar,
grooved, 1-0 - 4-0 Needle, 16 cm

Sugg. Retail Ft 22.767

Promo Ft 14.798



NHCW
Needle Holder Crile-Wood,
grooved, 3-0 - 6-0, 15 cm

Sugg. Retail Ft 20.7233

Promo Ft 13.470



NH5038
Needle Holder Crile-Wood
#5038, PermaSharp™, 15 cm

Sugg. Retail Ft 69.075

Promo Ft 44.899



NH5024R
Needle Holder Castroviejo
#5024R, round handle,
PermaSharp™, straight, 18 cm

Sugg. Retail Ft 200.624

Promo Ft 130.405



NH5024SIM
Needle Holder Simion #5024R,
round handle, PermaSharp™,
straight, 18 cm

Sugg. Retail Ft 200.624

Promo Ft 130.405

SUPER CUT™
Cutting blade honed/serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

PERMA SHARP™
Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

SCISSORS HU-FRIEDY®



SCORT
Scissors Cortellini curved 12 cm

Sugg. Retail Ft 27.842

Promo Ft 18.097



S16
Scissors Goldman-Fox #16
curved serrated 12,5 cm

Sugg. Retail Ft 23.096

Promo Ft 15.012



S14
Scissors LaGrange #14 2-x
angled 11,5 cm

Sugg. Retail Ft 30.849

Promo Ft 20.052



S18SC
Scissors Iris #18 fine angled
SuperCut™ 11,5 cm

Sugg. Retail Ft 40.270

Promo Ft 26.176

TISSUE PLIERS HU-FRIEDY®



TP41
Tissue Pliers Adson #41
anatomical straight 12 cm

Sugg. Retail Ft 13.956

Promo Ft 9.071



TP42
Tissue Pliers Adson #42
surgical 1x2-straight 12 cm

Sugg. Retail Ft 13.956

Promo Ft 9.071

PERIODONTAL KNIVES AND SINUS LIFT INSTRUMENTS HU-FRIEDY®

PERIODONTAL KNIVES HU-FRIEDY®



TKN1
Tunneling Knife #1 hdl #6

Sugg. Retail Ft 48.752

Promo Ft 29.251



TKN2
Tunneling Knife #2 hdl #6

Sugg. Retail Ft 48.752

Promo Ft 29.251



PFIWDS1MK
Tunneling knife hdl #6

Sugg. Retail Ft 26.056

Promo Ft 15.634



IMP65786
Sinuslift Instrument
Kramer-Nevins #LG/SM
hdl #6, 2x angled

Sugg. Retail Ft 50.162


Promo Ft 35.113



IMP6524S6
Sinuslift Instrument short hdl #6

Sugg. Retail Ft 50.162

Promo Ft 35.113



IMP6577SC6
Sinuslift Instrument
Kramer-Nevins #LG/SM
hdl #6, strongly curved

Sugg. Retail Ft 50.162

Promo Ft 35.113

ENDODONTICS HU-FRIEDY®

RUBBER DAM PUNCH HU-FRIEDY®



RDPA
Rubber Dam Punch Ainsworth, 17 cm

Sugg. Retail Ft 104.271

Promo Ft 62.526

RUBBER DAM FORCEPS HU-FRIEDY®



RDF4E
Rubber Dam Clamp Forceps #4 European

Sugg. Retail Ft 61.956

Promo Ft 37.174



EXDG16
Explorer Endo #DG16 hdl #31

Sugg. Retail Ft 15.977

Promo Ft 9.586

RUBBER DAM CLAMPS HU-FRIEDY®



For more information
on **RUBBER DAM CLAMPS**
visit [HuFriedy.eu](https://www.hufriedy.eu)



*Valid for all the stainless steel
Rubber Dam Clamps Line

RDCM00
Rubber Dam Clamps #00 UK anterior

Sugg. Retail Ft 8.200

Promo Ft 5.740*

IMS™ INSTRUMENT MANAGEMENT SYSTEM HU-FRIEDY®

IMS™ INFINITY SERIES™ CASSETTES



SIGNATURE SERIES™ TUBS

IMS-140_

Together, IMS™ Signature Series™ Tubs and Cassettes create the most efficient and organised system for managing instruments and consumable products.

Include:

- Tub Signature
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover



Sugg. Retail

Ft 37.765

Promo

Ft 24.542


Description

	CODE	Sugg. Retail	Promo
1 IMS™ Infinity Series™ Cassette DIN, 18 instruments, Surgical, Dimensions: 290 x 184 x 34 mm	IMEDIN18_S	Ft 96.470	Ft 62.706
2 IMS™ Infinity Series™ Cassette 1/4 DIN, 5 instruments, Dimensions: 72 x 184 x 34 mm	IME14DIN5_	Ft 59.301	Ft 38.546
3 IMS™ Infinity Series™ Cassette 1/2 DIN, 10 instruments, Dimensions: 152 x 184 x 34 mm	IME12DIN1_	Ft 66.373	Ft 43.143

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace "_" with the number)

Red	●	0	Lavender	●	6
Grey	●	1	Black	●	7
White	●	2	Blue	●	8
Orange	●	3	Green	●	9
Purple	●	4	Ocean Blue	●	11
Yellow	●	5			



For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit [HuFriedy.eu](https://www.HuFriedy.eu)
For more information on Crosstex products, call +39 05 24 52 74 53 or visit [crosstex.com](https://www.crosstex.com)
For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit [omniaspa.eu](https://www.omniaspa.eu)

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.
All prices indicated are intended. VAT included.
Hu-Friedy Mfg. Co., LLC. • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedy.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2021 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/0521