



NATIONAL PROMOTION HUNGARY

FALL/WINTER 2021

Valid from 22th of November 2021
until 15th of January 2022

Introducing...

HuFriedyGroup

The Best In Practice

Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

Service & Support

Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

Proud members of HuFriedyGroup:



IRRIGATION LINES OMNIA™

Omnia offers a wide range of mechanical irrigation lines, which are compatible with most physiodispensers available for sale.

Irrigation systems include a perforator, a roller or clamp for flow control, the fittings and an integrated peristaltic pump section.

OMNIA mechanical irrigation systems offer the following options:


a. Internal irrigation only

b. External irrigation only

c. Internal and external irrigation

The practical terminal cap makes the connection of the irrigation line to the hand piece easier and allows connecting our tubes with

all the cooling needles of the different hand pieces, whose diameters often differ.



PART CODE	DESCRIPTION	Sugg. Retail	Promo
32.F0134.00	Irrigation set 1/2YNSK Surgic XT, Surgic XT Plus, Osseocision®, VarioSurg, SurgicPro	Ft 36.182	Ft 27.137
32.F0139.00	Irrigation set 1/2Y Nobel Biocare80/ W&H Implantmed/ Implant Med new design 2005	Ft 36.182	Ft 27.137
32.F0014.00	Irrigation set 1/2Y Micro Dispenser 8000 con llave, MD10S, MD20 con llave	Ft 36.182	Ft 27.137
32.F0105.00	Irrigation set 1/2Y Aseptico Aeu 6000 series, Aseptico Aeu 7000 series, Saeyang Ki-20	Ft 36.182	Ft 27.137
32.F0108.00	Irrigation set Piezosurgery Touch Mectron	Ft 36.182	Ft 27.137
32.F0112.00	Irrigation set 1/2Y De Giorgi Steri Implant, INTRAMATIC 2, INTRAMAX PLUS, Solujet	Ft 36.182	Ft 27.137
32.F0147.00	Irrigation set 1/2Y Saeshin Traus SIP10, X-Cube, De Giorgi Traus sip 10	Ft 36.182	Ft 27.137
32.F0170.00	Irrigation set 1/2Y Bien Air-Hersel, Chiropro 980 / PIEZON® MASTER surgery/ Piezosurgery II	Ft 36.182	Ft 27.137
30.E0008.00	Hand Fuser Irrigation system	Ft 47.009	Ft 35.256

PROCEDURAL KIT OMNIA™

BASIC IMPLANT SET OMNIA*

The Basic Implant Set has been designed for all routine surgical activities where a surgical protocol ensuring an adequate degree of sterility and guaranteeing the safety of operators and patients is needed. Due to its reasonable price, it is especially fit for all surgical procedures with immediate loading or for operations, such as re-opening of implants, flap surgery, operations on minor surgical wounds

12.S5621.00

Basic Implant Set (sterile)



Description	Material	Color
1 Surgically folded Special Plus model gown (size L) and 2 hand towels wrapped in medical paper	SMS Microfibre	Light blue
1 Gown with cotton cuffs	NWF	Light blue
2 Surgical caps with fixed laces	Absorbent NWF	Light blue
1 Round surgical cap for the patient	Absorbent NWF	Light blue
1 Absorbent/liquid proof drape 50 x 75 cm with self-adhesive* side	NWF + PE	Light blue
2 Absorbent/liquid proof drape 75 x 90 cm with adhesive* U-shaped cut 11x9 cm	NWF + PE	Light blue
1 Omnisleeve sheath 120 x 7 cm	PE	Transparent
1 Disposal bag for contaminated material	PVC	Red

Sugg. Retail	Promo
Ft 10.338	Ft 8.270

IMPLANT SET OMNIA*

12.T3055.00

Implant Set



Description	Material	Color
3 Fluid Repellent Gowns with cuffs	NWF	Blue
2 Round surgical caps with elastics	NWF	Blue
1 Surgical cap with fixed laces	NWF	Blue
3 Surgical Masks with elastics	NWF	Blue
1 Absorbent/Waterproof drape 50x75cm with hole Ø 6 cm	NWF	Blue
1 Absorbent/Waterproof drape 75x90 cm	NWF	Blue
2 Absorbent/Waterproof drape 50x50cm	NWF	Blue
3 Pair of surgical gloves size available between 6 and 8 1/2	NWF	White

Sugg. Retail	Promo
Ft 7.049	Ft 5.639

PROCEDURAL SET OMNIA*

Procedural Set essential line, 8 Uds.

12.S6007.00

Procedural Set



Description	Material	Color
1 Liquid proof operator gown with cuffs	SPPB*	Light blue
2 Liquid proof assistant gowns with elastics	SPPB*	Light blue
1 Surgeon cap with laces	NWF absorbent	green
2 Assistant caps with elastic	NWF absorbent	Light blue
1 Surgeon face mask with laces	NWF	Light blue
2 Assistant face masks with elastic	NWF	Light blue
2 Plain drapes 50x75 cm	SPPB*	Light blue
1 Patient drape 75x90 cm with adhesive U shaped 11x9 cm	SPPB*	Light blue
1 Tubing sleeve 120x7 cm with elastics	PE	transparent

* Spunbonded Polypropylene

Sugg. Retail	Promo
Ft 54.508	Ft 43.607

DIAGNOSTIC INSTRUMENTS HU-FRIEDY®



PCP11

Probe #11 hdl #30 Qulix 3-6-8-11

Sugg. Retail Ft 15.248

Promo Ft 9.911



EXD5

Explorer #5 hdl #31

Sugg. Retail Ft 12.264

Promo Ft 7.972



DP2

Dressing Pliers College #2

Sugg. Retail Ft 10.573

Promo Ft 6.872



MHE6

Mirror Handle hdl #6 European Thread

Sugg. Retail Ft 10.855

Promo Ft 7.056



MIR5HDE

High Definition Mirror #5 single-sided, European, 1pc./pack

Sugg. Retail Ft 6.438

Promo Ft 4.184



CONTACT
YOUR DEALER TO
FIND OUT **ALL THE**
CODES ON
OFFER

EVEREDGE™ 2.0 INSTRUMENTS



EVEREDGE™ SCALERS ARE SHARP.
EVEREDGE™ 2.0 SCALERS ARE EVEN SHARPER.

Using an enhanced finishing process, Hu-Friedy's professional artisans are creating edges that are over 60 % sharper than the next closest competitor†. Less force means easier, more comfortable scaling for the clinician, which enables them to deliver the best clinical outcome every time.

†Products tested: Hu-Friedy EverEdge™ 2.0 H7, American Eagle XP H6-7, PDT H6-7, LM EM/EN H6-7, and G M Hartzell H6-7. For more information, contact David Platt, Hu-Friedy Mfg. Co., LLC, 3232 N. Rockwell St, Chicago, IL 60618.



AVAILABLE WITH DIFFERENT HANDLES:

HANDLE #9 EVEREDGE™ 2.0



HANDLE HARMONY™



NEW

SG13/1498E2

Curette Gracey #13/14 hdl #9 distal, blue

Sugg. Retail	Promo
Ft 26.432	Ft 18.502

SG5/6XE2

Curette Gracey #5/6 hdl #Harmony oral/ labial, purple

Sugg. Retail	Promo
Ft 32.423	Ft 27.560

For more information on EVEREDGE™ 2.0 Instruments, visit our site [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)



EVEREDGE™ 2.0 INSTRUMENTS

NEW

NEW HARMONY™ ERGONOMIC SCALERS AND CURETTES

REDUCED PINCH FORCE


REDUCES PINCH FORCE UP TO 65%⁴

LESS PRESSURE

REDUCES PRESSURE ON TOOTH BY 37%⁴

BALANCE

A PERFECTLY BALANCE WEIGHT



COLOR CODE

SECURE AND AGILE HANDLING



1) Hayes MJ, Cockrell D, Smith DR. A systematic review of musculoskeletal disorders among dental professionals. Int J Dent Hygiene. 2009;7:159-165. 2) Data on file. Available upon request 3) Compared to the Hu-Friedy #9 metal handle, these are the nominal values. Data on file. Available upon request. 4) Data on file. Available upon request. 5) Int J Dent Hygiene 7, 2009; 159-165 DOI: 10.1111/1601-5037.2009.00395.x. Hayes MJ, Cockrell D, Smith DR. A systematic review of musculoskeletal disorders among dental professionals. 6) Rempel, David, et al. "The Effects of Periodontal Curette Handle Weight and Diameter on Arm Pain." The Journal of the American Dental Association, vol. 143, no. 10, 2012, pp. 1105-1113. doi:10.14219/jada.archive.2012.0041. 7) Lalumandier, James A, and Scott D McPhee. "Prevalence and Risk Factors of Hand Problems and Carpal Tunnel Syndrome among Dental Hygienists." Journal of Dental Hygiene, vol. 75, no. 11, 2001, pp. 130-134. 8) Mulimani P, Hoe VCW, Hayes MJ, Idiculla JJ, Abas ABL, Karanth L. Ergonomic interventions for preventing musculoskeletal disorders in dental care practitioners. Cochrane Database of Systematic Reviews 2018, Issue 10. Art. No.: CD011261. DOI: 10.1002/14651858.CD011261.pub2. 9) Visit us online at www.hufriedy.eu/Harmony-Scaler to know the general conditions.



CONTACT
YOUR DEALER TO
FIND OUT **ALL THE**
CODES ON
OFFER

RESTORATIVE INSTRUMENTS



EXC18H

Excavator #18H hdl #41 heavy

Sugg. Retail Ft 17.974

Promo Ft 11.683



EXC19H

Excavator #19 hdl #41 heavy

Sugg. Retail Ft 17.974

Promo Ft 11.683



CIGFT36

Plastic Filling Instrument Goldstein #3,
hdl #6 Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434



CONTACT
YOUR DEALER TO
FIND OUT **ALL THE**
CODES ON
OFFER



CIGFT16

Plastic Filling Instrument Goldstein #1, hdl #6
Flexi-Thin

Sugg. Retail Ft 20.723

Promo Ft 12.434

SURGICAL INSTRUMENTS

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



CL856

Surgical Curette Lucas #85,
hdl #6 2,8 mm

Sugg. Retail Ft 39.988

Promo Ft 23.993



TPG1

Tissue Pliers Gerald hdl #
straight 18 cm

Sugg. Retail Ft 21.686

Promo Ft 14.096



PERMA SHARP™ INSTRUMENT

Featuring tungsten carbide
inserts bonded to surgical-
grade stainless steel on both
cutting edges for smooth,
controlled cutting and
maximum durability.

NH5024R

Needle Holder Castroviejo
#5024R, round handle,
PermaSharp, straight, 18 cm

Sugg. Retail Ft 200.624

Promo Ft 130.405



SUPER-CUT

Cutting blade honed/
serrated to „razor“
sharpness with opposing
blade finished with fine
serrations to avoid
tissue slippage.

SCORTSC

Scissors Cortellini curved,
SuperCut 12 cm

Sugg. Retail Ft 50.162

Promo Ft 32.605



**CONTACT
YOUR DEALER TO
FIND OUT ALL THE
CODES ON
OFFER**



TKN1

Tunneling Knife #1 hdl #6

Sugg. Retail Ft 48.752

Promo Ft 31.689



PPBUSER6

Periosteal Buser hdl #6

Sugg. Retail Ft 42.033

Promo Ft 27.321

ENDODONTIC INSTRUMENTS HU-FRIEDY®

HU-FRIEDY'S ENDODONTIC INSTRUMENTATION

Hu-Friedy's comprehensive line of endodontic instruments and accessories give you versatile options to meet all of your clinical needs. If you are looking for exceptional instrument design, quality and reliability, and unmatched customer service, choose Hu-Friedy's Endodontic Instruments.

THE CHALLENGE

The sensitive nature of endodontic procedures necessitates instruments with exceptional response, integrity and durability to accommodate a wide variety of endodontic procedures — from routine root canals to complex microsurgery.

THE SOLUTION

Hu-Friedy instruments are individually handcrafted to precise specifications from quality materials to ensure superior performance, longevity and reliability. Instrument design combined with proprietary materials and incomparable workmanship result in instruments with unique balance, working leverage and control. They function as a seamless extension of the practitioner's skill.

THE RESULT

high-quality Hu-Friedy instrumentation provides unsurpassed clinical efficiency, so practitioners can deliver superior outcomes and focus on the patient. The result is outstanding patient care.

RUBBER DAM PUNCH



RDPA

Rubber Dam Punch Ainsworth, 17 cm

Sugg. Retail Ft 104.271

Promo Ft 72.990

RUBBER DAM CLAMPS



RDCM00

Rubber Dam Clamps #00
UK anterior

Sugg. Retail Ft 8.200

Promo Ft 5.740

EACH



CONTACT
YOUR DEALER TO
FIND OUT **ALL THE**
CODES ON
OFFER

IMS™ INSTRUMENT MANAGEMENT SYSTEM HU-FRIEDY

IMS™ INFINITY SERIES™ CASSETTES

IMS™ Infinity Series™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, cassettes can be specialized to meet your needs.



IME14DIN50

IMS™ Infinity Series™ Cassette 1/4 DIN, 5 instruments, red.
Dimensions: 72 x 184 x 34 mm

Sugg. Retail Ft 59.301

Promo Ft 38.546

Main Benefits



Customized kits for surgical and non-surgical procedures



Better Safety



Streamlined Reprocessing



Improved Organization



Instrument Protection



Return on Investment

SIGNATURE SERIES™ TUBS

IMS-140_

Together, IMS™ Signature Series™ Tubs and Cassettes create the most efficient and organised system for managing instruments and consumable products.

Include:

- Tub Signature
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover



Sugg. Retail

Ft 37.756

Promo

Ft 24.542

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace "_" with the number)

Red	●	0	Lavender	●	6
Grey	●	1	Black	●	7
White	●	2	Blue	●	8
Orange	●	3	Green	●	9
Purple	●	4	Ocean Blue	●	11
Yellow	●	5			



For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit HuFriedy.eu

For more information on Crosstex products, call +39 05 24 52 74 53 or visit crosstex.com

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit omniaspa.eu

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

All prices indicated are intended VAT included.

All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.

©2021 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/1121