



BLACK DAYS promotions HUNGARY

.....

Valid from 21st until 25th
November 2022

Introducing...

HuFriedyGroup

The Best In Practice

Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

Service & Support

Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

Proud members of HuFriedyGroup:



BLACK DAYS promotions HUNGARY

.....

Valid from 21st until 25th
November 2022

DIAGNOSTIC

EXPLORERS P.4

SURGICAL

TISSUE PLIERS P.4

SCISSORS P.5

HEMOSTATS P.6

SCALPEL HANDLES P.6

CHEEK RETRACTORS P.6

RONGEURS P.6

NEEDLE HOLDERS P.7

TISSUE NIPPERS P.8

SURGICAL CURETTES P.8

SYRINGES P.8

RETRACTORS P.8

PERIOSTEALS P.9

HAND SCALING

CURETTES AND SCALERS EVEREDGE™ 2.0 P.9

RESTORATIVE


COMPOSITE INSTRUMENTS P.10

GINGIVAL CORD PACKERS P.10

IMS™ INSTRUMENT MANAGEMENT SYSTEM

IMS™ INFINITY SERIES™ CASSETTES P.11

EXPLORERS




EXD3CH
#3CH DE Cow Horn,
#31 Handle

Retail 16.408,40 Ft
Promo 9.024,62 Ft

EXS23
#23 Shepherd's Hook SE,
#30 Handle

Retail 11.415,71 Ft
Promo 6.278,64 Ft

TISSUE PLIERS



PERMA SHARP™
Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

TP5042
Adson 1x2,
Perma Sharp™, 12 cm

Retail 100.744,34 Ft
Promo 55.409,39 Ft

TP43
#43 Adson-Brown,
Straight, 12 cm

Retail 33.923,29 Ft
Promo 18.657,81 Ft


TPG1
Gerald,
Straight/Serrated, 18 cm

Retail 27.473,28 Ft
Promo 15.110,30 Ft

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**




SCISSORS



S1
#1 Kelly, Curved,
16 cm

Retail 33.626,43 Ft
Promo 18.494,53 Ft



S2
#2 Kelly, Straight,
16 cm

Retail 32.304,04 Ft
Promo 17.767,22 Ft




SCORTSC
Cortellini, Super-Cut,
Curved, 12 cm

Retail 62.934,85 Ft
Promo 37.760,91 Ft



S5080
Goldman-Fox
Perma Sharp™, Straight,
13 cm

Retail 114.912,78 Ft
Promo 68.947,67 Ft



S19SCX
Metzenbaum, Super-Cut,
14.25 cm, Black Line

Retail 121.956,51 Ft
Promo 73.173,91 Ft



S18SCX
#18 Iris, Curved,
Super-Cut, 11.5 cm,
Black Line

Retail 68.224,40 Ft
Promo 34.112,20 Ft

SUPER-CUT
Cutting blade honed/ serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

PERMA SHARP™
Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY Instrument Lubricant Spray (**ILS**), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



HEMOSTATS



HKRC

Kelly-Rankin, Curved
16 cm

Retail 21.913,85 Ft

Promo 13.148,31 Ft

SCALPEL HANDLES



10-130-03E

#3E European Scalpel
Handle

Retail 9.013,83 Ft

Promo 5.858,99 Ft

CHEEK RETRACTORS



CRW2

Weider Cheek & Tongue
Retractor, Large

Retail 46.769,34 Ft

Promo 23.384,67 Ft

RONGEURS



R1

#1 Hu-Friedy™
16,5 cm

Retail 149.186,90 Ft

Promo 82.052,80 Ft

RMF

Mini 90° Friedman,
14,5 cm

Retail 132.076,83 Ft

Promo 72.642,25 Ft

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



NEEDLE HOLDERS


1-0, 2-0, 3-0, 4-0



NH5074E
Mathieu-Kocher Perma Sharp™, 14 cm

Retail 123.063,00 Ft
Promo 73.837,80 Ft

6-0, 7-0, 8-0



NHDPV
Diamond Dusted
Castroviejo Microsurgical,
Straight, 18 cm

Retail 428.453,55 Ft
Promo 257.072,13 Ft


5-0, 6-0, 7-0, 8-0



NH5052
Fine Swedish Pattern,
Perma Sharp™, 15 cm

Retail 94.753,11 Ft
Promo 56.851,87 Ft

3-0, 4-0, 5-0, 6-0




NH5068
Olsen-Hegar,
Perma Sharp™, 17 cm

Retail 114.912,78 Ft
Promo 68.947,67 Ft

PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.


3-0, 4-0, 5-0, 6-0



NHOH
Olsen-Hegar,
17 cm

Retail 76.212,70 Ft
Promo 45.727,62 Ft

6-0, 7-0, 8-0



NH5084
Micro-Vascular Pattern
Perma Sharp™, 18 cm

Retail 130.295,65 Ft
Promo 78.177,39 Ft

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



TISSUE NIPPERS



NIPS
Goldman-Fox, 10 cm

Retail 88.734,90 Ft
Promo 26.620,47 Ft

SURGICAL CURETTES



CL866
#86 Lucas DE, 2.8 mm,
#6 Short Handle

Retail 50.169,76 Ft
Promo 30.101,86 Ft



CM2/4
#2/4 Molt DE, 7-3.5 mm
#522 Handle

Retail 49.306,16 Ft
Promo 29.583,70 Ft


SYRINGES



SYBG
Bone Grafting,
7 mm

Retail 100.258,56 Ft
Promo 60.155,14 Ft

RETRACTORS



IR
Implant,
23 cm


Retail 159.766,00 Ft
Promo 95.859,60 Ft

PERIOSTEALS



PPBUSER6
Buser Elevator,
#6 Handle

Retail 52.220,81 Ft
Promo 31.332,49 Ft



P14
#14 Goldman-Fox,
#2 Short Handle

Retail 45.473,94 Ft
Promo 27.284,36 Ft



PFIWDS1MK
Urban Mini Me
Microsurgical,
#6 Handle

Retail 32.385,00 Ft
Promo 19.431,00 Ft



P24G
24G, #2 Handle

Retail 45.473,94 Ft
Promo 27.284,36 Ft

CURETTES AND SCALERS EVEREDGE™ 2.0



SAS5/695E2
Gracey #5/6 Mini Five™,
yellow, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SAS7/897E2
Gracey #7/8 Mini Five™,
green, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SG7/897E2
Standard Gracey #7/8
green, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SG5/695E2
Standard Gracey #5/6,
yellow, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SRPG7/897E2
Gracey #7/8 After Five™,
green, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



S204SD9E2
#204SD Posterior Sickle
Scaler, purple, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SRP11/1293E2
Gracey 11/12 After Five™
orange, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SG13/1498E2
Standard Gracey #13/14,
blue, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft



SH6/79E2
#H6/7 Sickle Scaler,
purple, #9 Handle

Retail 33.788,35 Ft
Promo 18.583,59 Ft

COMPOSITE INSTRUMENTS



TNCVIPC
Interproximal Carver,
XTS Handle

Retail 29.794,20 Ft

Promo 17.876,52 Ft



TNPFIW3
#W3,
XTS Handle

Retail 29.794,20 Ft

Promo 16.386,81 Ft



PFIA6
#A6,
#41 Handle

Retail 22.750,46 Ft

Promo 13.650,28 Ft



CVIPC
Interproximal Carver,
#41 Handle

Retail 22.534,56 Ft

Promo 13.520,74 Ft

GINGIVAL CORD PACKERS



GCP1136
#113, Serrated, #6 Handle

Retail 27.122,44 Ft

Promo 16.273,46 Ft



APF2E
Articulating Paper
Forceps

Retail 25.314,28 Ft












Promo 13.922,85 Ft

IMS™ INSTRUMENT MANAGEMENT SYSTEM HU-FRIEDY

IMS™ INFINITY SERIES™ CASSETTES

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace “_” with the number)

Red		0
Grey		1
White		2
Orange		3
Purple		4
Yellow		5
Lavender		6
Black		7
Blue		8
Green		9
Ocean Blue		11

IMEDIN13_

IMS™ Infinity Series™
1 DIN Cassette,
13 Instruments,
Dimensions:
290 x 184 x 34 mm

Retail 113.968,21 Ft

Promo 68.380,93 Ft



IME12DIN1_

IMS™ Infinity Series™ 1/2 DIN
Cassette with accessory area,
10 Instruments.
Dimensions:
152 x 184 x 34 mm

Ideal for diagnostics and prophylaxis.

Retail 84.039,08 Ft

Promo 50.423,45 Ft

IME14DIN5_

IMS™ Infinity Series™
1/4 DIN Cassette - 5
Instruments.
Dimensions:
72 x 184 x 34 mm

Retail 75.079,23 Ft

Promo 45.047,54 Ft



For more information on HuFriedyGroup products or to place an order,
call at 00800 48 37 43 39 or visit [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.
©2022 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/1122