



-45%

-50%

-60%

**BLACK
DAYS**
promotions
H U N G A R Y

.....

Valid from 20th until 24th
November 2023

INDEX

Personal Protection Equipment

ULTRA™ EARLOOP BLACK FACEMASK
with Secure Fit™ Mask Technology P.1

Diagnostic

DRESSING PLIERS P.2

MIRRORS P.2

PROBES P.2

EXPLORERS P.2

Hand Scaling

CURETTES AND SCALERS
EVEREDGE™ 2.0 P.3

Restorative

COMPOSCULP INSTRUMENT KIT P.3

PAPER FORCEPS P.3

EXCAVATORS P.3

CROWN REMOVERS P.4

PRIMARY CROWNS P.4

Endodontic

MEASURING DEVICES P.5

ENDODONTIC MIRRORS P.5

LOCKING PLIERS P.5

POST & SILVER POINT
REMOVAL FORCEPS P.5

ENDODONTIC EXPLORERS P.5

RUBBER DAM FRAME P.5

RUBBER DAM PUNCH P.5

Surgical

SCALPEL BLADES AND HANDLES P.6

CHEEK RETRACTORS P.6

PERIOSTEALS P.6

SCISSORS P.7

SURGICAL CURETTES P.8

TISSUE PLIERS P.8

NEEDLE HOLDERS P.9

PERSONAL PROTECTION EQUIPMENT CROSSTEX

-60%

ULTRA™ EARLOOP BLACK FACEMASK with Secure Fit™ Mask Technology

FEATURE & BENEFIT

Secure Fit™ Technology

The innovative Crosstex™ Ultra™ Earloop Facemask with patented Secure Fit™ Technology creates a snug, custom fit that adapts to most size or shape face to significantly minimizes gapping and provides high level of protection.

Fluid resistant outer layer

Spunbond outer layer material for maximum protection against fluids traveling through the mask.

Fluid resistant white inner layer

Unique layer barrier for maximum protection against fluids.

Extra-long aluminium nose piece

Significantly reduces gapping for maximum protection.

Not made with natural rubber latex

Prevents latex allergy

Fluid resistance of 160 mmHg

Ideal for procedures where high amounts of fluid, spray and/or aerosols are produced.



GCFBKSF

Ultra™ Earloop Facemask Secure Fit™ Technology, Black, (EN14683 Type IIR Enhanced), 50 pcs/box

Retail 8.852 Ft

Promo 3.541 Ft

-50%

DRESSING PLIERS



DP18
#18

Retail 14.722 Ft
Promo 7.361 Ft

DP2
#2 College

Retail 13.435 Ft
Promo 6.717 Ft

MIRRORS



MIR5DSHDE
#5 HD Double-Sided Mirror, EU Simple Stem, single

Retail 10.190 Ft
Promo 5.095 Ft

MHE6
#6 European Simple Stem Mirror Handle

Retail 13.408 Ft
Promo 6.704 Ft

PROBES

1-2-3-5-7-8-9-10 mm

3-6-8-11 mm



PQW
Qulix SE Williams Probe, handle #30

Retail 15.875 Ft
Promo 7.938 Ft



PCP116
Qulix SE Probe, handle #6

Retail 20.970 Ft
Promo 10.485 Ft

EXPLORERS



EXD11/126
#11/12 ODU DE Explorer handle #6

Retail 17.109 Ft
Promo 8.554 Ft

HAND SCALING HU-FRIEDY

RESTORATIVE HU-FRIEDY

CURETTES AND SCALERS EVEREDGE™ 2.0



-45%

Valid for the whole portfolio of EverEdge™ 2.0 Scalers and Curettes handle #9

SG5/695E2

Gracey 5/6 Curette, Standard, handle #9

Retail 33.574 Ft

Promo 18.466 Ft

SG13/1498E2

Gracey 13/14 Curette, Standard, handle #9

Retail 33.574 Ft


Promo 18.466 Ft

COMPOSCULP INSTRUMENT KIT

SMOOTH HANDLES

-45%

CompoSculp instruments were designed by Dr. Didier Dietschi with the intention of aiding clinicians in achieving superior esthetic restorations. The CompoSculp kit is comprised of 5 complimentary composite instruments for both anterior and posterior restorations.



Contains:

Dietschi CompoSculp 1/2, #CCH Handle	PFIDD1/2
Dietschi CompoSculp 3/4, #CCH Handle	PFIDD3/4
Dietschi CompoSculp 5/6, #CCH Handle	PFIDD5/6
Dietschi CompoSculp 7/8, #CCH Handle	PFIDD7/8
Dietschi CompoSculp 9/10, #CCH Handle	PFIDD9/10
IMS™ Infinity Series™ 1/4 DIN Cassette, 5 Instruments, Orange	IME14DIN53

PFIDDDIN

CompoSculp Instrument Kit, Smooth Handles #CCH

Retail 192.754 Ft

Promo 106.015 Ft

PAPER FORCEPS

-60%



APF2E

Miller Articulating Paper Forceps

Retail 25.153 Ft

Promo 106.015 Ft

EXCAVATORS



-50%

EXC18H

18 Heavy Excavator, handle #41

Retail 23.035 Ft

Promo 11.517 Ft

RESTORATIVE HU-FRIEDY

CROWN REMOVERS



-40%

CRU
Trial Crown
Remove

Retail 142.661 Ft
Promo 85.597 Ft

CRL
Trial Crown
Remover, Lower

Retail 142.661 Ft
Promo 85.597 Ft

CRCH1
Christensen Crown
Remover Straight,
Anterior

Retail 51.728 Ft
Promo 31.037 Ft

PRIMARY CROWNS

Stainless Steel Pedo Crowns come in 48 sizes for use on the first and second primary molars. These crowns feature a more narrow mesial-distal maxillary first primary molar width for quicker adjacent crown placement.

-45%



Upper Right Primary 1st #5, 5 pcs/box	SSC-URD5-CE
Upper Right Primary 2nd #3, 5 pcs/box	SSC-URE3-CE
Upper Left Primary 1st #3, 5 pcs/box	SSC-ULD3-CE
Upper Left Primary 1st #4 5 pcs/box	SSC-URD4-CE

Retail 13.301 Ft each Promo 7.315 Ft each

ENDODONTIC HU-FRIEDY

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



MEASURING DEVICES

-45%



CLR6

Gracey 13/14 Curette, Standard, handle #9

Retail 19.629 Ft

Promo 10.796 Ft

ENDODONTIC MIRRORS

-45%



MM5

Micro Mirror Round, 5 mm

Retail 53.981 Ft

Promo 29.689 Ft

-45%



MM8

Micro Mirror, 3 mm x 9 mm

Retail 66.450 Ft

Promo 36.548 Ft

LOCKING PLIERS

-50%



EPL1


1 Endo Locking Pliers, 160-5 mm

Retail 24.751 Ft

Promo 12.376 Ft

POST & SILVER POINT REMOVAL FORCEPS

-50%



RF45

Steiglitz (Archer) 45°, 130-5 mm

Retail 67.308 Ft

Promo 33.654 Ft

ENDODONTIC EXPLORERS

-50%



EXDG16

#DG16 DE, Endodontic Explorer, handle #31

Retail 19.737 Ft

Promo 9.868 Ft

RUBBER DAM FRAME

-50%

RDAF6

Rubber Dam Frame, Adult - 15 cm

Retail 20.997 Ft

Promo 10.498 Ft

RUBBER DAM PUNCH

-50%



RDPA

Rubber Dam Punch, Ainsworth, 17,5 cm

Retail 131.184 Ft

Promo 65.592 Ft

SCALPEL BLADES AND HANDLES

-60%



10-255-11
#11 Sterile Scalpel
Blades, 100 pcs/box,
Carbon Steel

Retail 24.751 Ft

Promo 9.900 Ft

-45%



10-130-05E
#5 Scalpel Handle

Retail 45.265 Ft

Promo 24.896 Ft

-60%



10-256-12
#12D Sterile Scalpel
Blades, 100 pcs/box,
Carbon Steel

Retail 29.578 Ft

Promo 11.831 Ft

CHEEK RETRACTORS

-45%



CRM2
Cawood-Minnesota
Retractor

Retail 20.782 Ft

Promo 11.430 Ft

CRO2
Orringer Retractor
Medium

Retail 27.031 Ft

Promo 14.867 Ft

CRW2
Weider Cheek & Tongue
Retractor, Large

CRW1
Weider Cheek & Tongue
Retractor, Small

Retail 46.472 Ft each

Promo 25.560 Ft each

PERIOSTEALS

-45%



PWILL
Williger Periosteal,
Large

Retail 43.120 Ft

Promo 23.716 Ft



PP25
#25 Palti Periosteal
Large

Retail 54.410 Ft

Promo 29.925 Ft

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



SCISSORS

SUPER-CUT

Cutting blade honed/ serrated to „razor“ sharpness with opposing blade finished with fine serrations for reduced tissue slippage.

PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

-60%



S14

#14 LaGrange Scissors,
11,5 cm, 1 serrated blade,
Standard

Retail 38.454 Ft

Promo 15.382 Ft

-60%



S16

#16 Iris Scissors, 12,5
cm, 1 serrated blade,
Standard

Retail 28.800 Ft

Promo 11.520 Ft

-45%



S5082

Iris Perma Sharp™
Scissors, Straight,
11,5 cm

Retail 103.430 Ft

Promo 56.886 Ft

-45%



SCORT

Scissors Cortellini,
Curved, 12 cm

Retail 34.700 Ft

Promo 19.085 Ft

-45%



SCORTSC

Scissors Cortellini,
Super-Cut, Curved,
12 cm, 1 serrated blade

Retail 62.535 Ft

Promo 34.394 Ft

SURGICAL CURETTES

-45%



CL88

#88 Lucas DE Surgical Curette, handle #522

Retail 50.629 Ft

Promo 27.846 Ft

-45%



CM46

#4 Molt SE Curette, handle #6

Retail 43.925 Ft

Promo 24.159 Ft

TISSUE PLIERS

-45%



TP42

#42 Adson 1x2 Tissue Pliers, 12 cm

Retail 17.565 Ft

Promo 9.660 Ft



TP5043

Adson-Brown Perma Sharp™ Tissue Plier, 12 cm

Retail 105.736 Ft

Promo 58.155 Ft



TPG1

Gerald Tissue Pliers, Straight/Serrated, 17,78 cm

Retail 27.299 Ft

Promo 15.014 Ft



TPDAPV

Anatomic Microsurgical Tissue Plier, 18 cm

Retail 204.982 Ft

Promo 112.740 Ft



TPDSSPV

Microsurgical Tissue Plier, Straight, 18 cm

Retail 261.269 Ft

Promo 143.698 Ft

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



NEEDLE HOLDERS

1-0, 2-0, 3-0, 4-0



NH5072

Olsen-Hegar Perma Sharp™ Needle Holder, 14 cm

Retail 108.203 Ft

Promo 59.512 Ft

5-0, 6-0, 7-0, 8-0, 9-0



SPNHDPVN

Swiss Perio Microsurgical Needle Holder, 18 cm

Retail 424.793 Ft

Promo 233.636 Ft

6-0, 7-0, 8-0, 9-0, 10-0




NHTMBH

Titanium Diamond Dusted Castroviejo Microsurgical Needle Holder, 18 cm

Retail 671.876 Ft

Promo 369.532 Ft

6-0, 7-0, 8-0



NHDPV

Diamond Dusted Castroviejo Microsurgical Needle Holder, Straight, 18 cm

Retail 425.732 Ft

Promo 234.152 Ft

1-0, 2-0, 3-0, 4-0



NHMH

Mayo-Hegar Needle Holder, 16 cm

Retail 28.130 Ft

Promo 15.472 Ft

6-0, 7-0, 8-0, 9-0, 10-0, 11-0




NHM-5027

Castroviejo Microsurgical Needle Holder, Curved, 13 cm

Retail 224.611 Ft

Promo 123.536 Ft



1-0, 2-0, 3-0, 4-0



NH5074

Mathieu-Kocher Perma Sharp™ Needle Holder, 14 cm

Retail 109.865 Ft

Promo 60.426 Ft

For more information on HuFriedyGroup products or to place an order,
call at 00800 48 37 43 39 or visit **HuFriedyGroup.eu**

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/1123