

# SAFESCRAPER

## TWIST *curve*

# VOLUMIZER



### Extra volume version

Developed from the long-standing experience of the SAFESCRAPER family.

The VOLUMIZER offers the same performance while collecting a larger volume.


Enables rapid and accurate collection of autologous cortical bone.

Quality of harvested bone means presence of many viable osteocytes.

Promotes graft acceptance, neovascularization and bone regeneration.

Collected bone already contains blood and is ready to use.

Double grooves on the blade allow for three collecting surfaces.



# SAFESCRAPER

## TWIST *curve*

# VOLUMIZER



- Based on state-of-the-art technology and current medical practices
- Ergonomic curved body easier for surgeons to use
- Cutting performance reduces patient discomfort
- Blade has a cutting range of 160°
- Works on plane, concave and convex surfaces
- Effective for both small and large bone defects
- Devices are disposable and come ready to use in single sterile packages

**Ref. 5501** SAFESCRAPER TWIST *curve* VOLUMIZER (chamber capacity 2.5 cc) 3 disposable units, sterile packaging



META

Art and Innovation in Medical Technology

**META TECHNOLOGIES S.R.L.**

VIA E. VILLA 7 42124 REGGIO EMILIA ITALY  
TEL. +39 0522 502311

► [info@metahosp.com](mailto:info@metahosp.com)

► [www.metahosp.com](http://www.metahosp.com)

