


Validity: June 10th - July 22nd 2026 | Hungary

promotion



UP TO
-45%

*Surgical
Instruments
Special*

**Swiss Perio
Microsurgical
Scissors
18 cm
SPSPV**

Page 3

A selection
of Scissors,
Tissue Pliers,
Needle Holders
and Atraumair
Extraction Forceps.



Scissors HU-FRIEDY

-35%



#16 Goldman-Fox Scissors, Curved, Super-Cut 12,5 cm, 1 serrated blade, Black Line
S16SCX

58.577,48 Ft
90.119,20 Ft



#17 Iris Scissors, Straight, Super-Cut 11,5 cm, 1 serrated blade, Black Line
S17SCX

44.577,00 Ft
68.580,00 Ft



#18 Iris Scissors, Curved, Super-Cut 11,5 cm, 1 serrated blade, Black Line
S18SCX

44.577,00 Ft
68.580,00 Ft



#14 LaGrange, Curved, Super-Cut 11,5 cm, 1 serrated blade, Black Line
S14SCX

44.577,00 Ft
68.580,00 Ft



Metzenbaum Scissors, Super-Cut 14,25 cm, 1 serrated blade, Black Line
S19SCX

77.597,00 Ft
119.380,00 Ft



Goldman-Fox Scissors, Straight, Super-Cut, 13,5 cm, 1 serrated blade, Black Line
S20SCX

77.597,00 Ft
119.380,00 Ft



Metzenbaum Perma Sharp™ Scissors, Curved, Delicate 18 cm
S5069

77.398,88 Ft
119.075,20 Ft

Black Line

- More contrast between the instrument and the tissue.
- Avoid light reflection.
- More sharpened than standard.

Scissors HU-FRIEDY

-35%





Castroviejo
Microsurgical
Scissors, 18 cm,
Standard
SPV

193.332,10 Ft
297.434,00 Ft

Castroviejo
Microsurgical
Scissors
18 cm, 1 serrated
blade, Black Line
SPVX

213.144,10 Ft
327.914,00 Ft






Swiss Perio
Microsurgical
Scissors, 45°
18 cm
SPSPV45

257.721,10 Ft
396.494,00 Ft

Instruments developed
in collaboration with
Swiss Perio Group and
Dr. Rino Burkhardt


“
My surgical
instruments are
the sensory-
motor extension
of my hands.
”

Dr. Rino Burkhardt 



Swiss Perio
Microsurgical
Scissors
18 cm
SPSPV

213.144,10 Ft
327.914,00 Ft


**Discover
more**



Tissue Pliers HU-FRIEDY

-35%



**#42 Adson 1x2
Tissue Pliers**
12 cm
TP42

14.859,00 Ft
22.860,00 Ft



**Adson 1x2 Perma
Sharp™ Tissue Plier**
12,07 cm
TP5042

69.342,00 Ft
106.680,00 Ft



**Cushing Perma
Sharp™ Tissue Plier**
17,78 cm
TP5060

60.426,60 Ft
92.964,00 Ft



**General Perma
Sharp™ Tissue Plier**
15,24 cm
TP5070

56.134,00 Ft
86.360,00 Ft




**Perma Sharp™
Tissue Plier**
15,30 cm
TP5072

60.426,60 Ft
92.964,00 Ft

Tissue Pliers HU-FRIEDY

-35%



Corn Suture Pliers are ideal for use with guided tissue regeneration procedures.

Swiss Perio
Microsurgical Plier,
Anatomical
18,10 cm
TPDAPV

138.684,00 Ft
213.360,00 Ft

1x2 Curved
Diamond Dusted
Stainless Steel,
17,78 cm
TPDSCPV

179.959,00 Ft
276.860,00 Ft

Corn
Suture Plier
15.24 cm
SP20SLEM

208.191,10 Ft
320.294,00 Ft

#20 Corn
Suture Plier
15,24 cm
SP20

15.849,60 Ft
24.384,00 Ft



Needle Holders


HU-FRIEDY

-40%




Mathieu-Kocher Perma Sharp™ Needle Holder
1-0, 2-0, 3-0, 4-0
14 cm
NH5074E

71.932,80 Ft
119.888,00 Ft





Diamond Dusted



Diamond Dusted Castroviejo Microsurgical Needle Holder, Straight, 6-0, 7-0, 8-0, 18 cm
NHDPV

251.460,00 Ft
419.100,00 Ft


Perma Sharp™ Cortellini Castroviejo Microsurgical Needle Holder, Straight Round Handle 5-0, 6-0, 7-0, 8-0, 18 cm
NH5024CORT

152.247,60 Ft 253.746,00 Ft

[Discover more](#)

Fine Swedish Pattern Perma Sharp™ Needle Holder
5-0, 6-0, 7-0, 8-0
15 cm
NH5052

59.436,00 Ft
99.060,00 Ft



Olsen-Hegar Needle Holder Scissors
3-0, 4-0, 5-0, 6-0,
17 cm
NHOH

48.768,00 Ft
81.280,00 Ft



Needle Holders HU-FRIEDY

-40%



**Perma Sharp™
Castroviejo Needle
Holder, Straight**
5-0, 6-0, 7-0, 8-0
14 cm
NH5020

122.834,40 Ft
204.724,00 Ft



**Diamond
Dusted**

**Diamond Dusted
Castroviejo
Microsurgical Needle
Holder, Curved** 5-0,
6-0, 18 cm
NHDCPVN

256.032,00 Ft
426.720,00 Ft



**Perma Sharp™
Diamond Dusted
Castroviejo
Microsurgical Needle
Holder, Straight
Round Handle**
6-0, 7-0, 8-0, 9-0
18 cm
NHM-5026R

138.684,00 Ft
231.140,00 Ft




**Diamond
Dusted**

**Simion
Perma Sharp™
Castroviejo
Needle Holder**
3-0, 4-0
18 cm
NH5024SIM

152.247,60 Ft
253.746,00 Ft

**Perma Sharp™
Castroviejo
Needle Holder,
Curved, Round
Handle**, 5-0, 6-0,
7-0, 8-0, 9-0, 10-0,
18 cm
NH5024RC

152.247,60 Ft
253.746,00 Ft



**Mayo-Hegar Perma
Sharp™ Needle
Holder** 1-0, 2-0,
3-0, 4-0, 18 cm
NH5044

59.436,00 Ft
99.060,00 Ft



Atraumair Extraction Forceps HU-FRIEDY

-45%



#51 Atraumair Forceps, Serrated Upper Root FX51XS



#35 Atraumair Apical Forceps Upper Can/Premolar FAF35XS



#35 Atraumair Apical Forceps Upper Can/Premolar F35XS



#23 Atraumair Forceps Cowhorn, Lower Molars F23XS



#217 Atraumair Forceps 1st and 2nd Molars F21XS



#36 Atraumair Forceps Lower Incisors F36XS



#73 Atraumair Forceps Lower Molars FX73XS

120-853,20 Ft **66.469,26 Ft**

Atraumair Extraction Forceps HU-FRIEDY

-45%



#1 Atraumair Standard Forceps Upper Incisors
F1XS




#67 Diamond Dusted Apical Forceps Upper 3rd Molar
F67XS



#45 Atraumair Forceps, Serrated Lower Root
FX45XS



#13 Atraumair Forceps Lower Premolars
F13XS



#44 Atraumair Forceps Universal, Upper Root
FX44XS



#36 Atraumair Apical Forceps Lower Premolars/ Incisors
FAF36XS



79 Diamond Dusted Forceps Lower Molars
FX79XS



#1 Atraumair Apical Forceps Upper Incisors
FAF1XS

120:853,20 Ft **66.469,26 Ft**



Instrument Lubricant HU-FRIEDY

Instrument care and maintenance products ensure the proper functioning of instruments, while also minimize wear.



DID YOU KNOW?
 All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**

Instrument Lubricant Spray, Keeps hinged instruments moving freely, 237 ml ILS


18.034,00 Ft



Are you Creating or Renovating a Practice?



Would you like to be organized and do you have many questions about storage, the sterilization chain and ergonomics?



Take advantage of our support and let us guide you!


**Free call
00800 48 37 43 39**

Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.



Our Contacts



Dealer informations

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. Available while stocks last. All prices are VAT included.

©2026 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/0626